


# STRAIN SPEC SHEET

# AMB-1

*Magnetospirillum Magneticum*




## STRAIN OVERVIEW

*Magnetospirillum magneticum* AMB-1 naturally takes in iron to form magnetosomes, a membranous organelle that allows the organism to align with a magnetic field. This organelle helps the strain navigate towards optimal, microaerophilic environmental conditions. By varying culturing conditions, researchers are able to control magnetosome formation. Researchers have used magnetic bacteria to determine magnetic poles in meteorites and rocks, to synthesize gold nanoparticles, and to serve as nanorobot. Potential applications include removal of heavy metal and radionuclides, enzyme immobilization, and oxygen depletion measurement.<sup>1</sup>


**MORPHOLOGY:** spiral

**DOUBLING TIME ON PLATES:** 8 hours

**COLONY FORMATION:** 48 hours



## CULTURING



### CHEMICALS FOR ENRICHED MAGNETIC SPIRILLUM GROWTH MEDIA (E-MSGM)


- 90% total volume milliQ water
  - L-ascorbic acid 0.035 g/L
  - Sodium acetate 0.07 g/L
  - Sodium thiosulfate 0.1 g/L
  - Sodium nitrate 0.12 g/L
  - Succinic acid 0.37 g/L
  - L-tartaric acid 0.37 g/L
  - Potassium phosphate monobasic 0.68 g/L
  - Wolfe's Mineral Solution 5 mL/L [ATCC #MD-TMS]
  - Yeast extract 0.1 g/L
  - Polypeptone 0.2 g/L
  - L-cysteine 0.05 g/L
- pH adjusted to 6.9 with 10M NaOH

### FERRIC MALEATE SUPPLEMENT NEEDED FOR MAGNETOSOME FORMATION.

Iron (III) chloride 0.486 g/L + DL-Malic acid 1.207 g/L  
Prepare with milliQ water to make 100X stock solution.

**30°C**  
**OPTIMAL TEMPERATURE**

# QUANTIFYING CELLS



# PLATING TECHNIQUES

For agar plates, add 10.0 g of agar per liter of E-MSGM. Add 1% of total volume E-MSGM in Wolfe's Vitamin Solution [ATCC #MD-VS] and 1% of total volume E-MSGM in iron malleate solution.

## METHOD 1

Bacteria can be plated on 1% agar directly, as done for E. coli. Plates must be wrapped using an excess of Parafilm M® to reduce oxygen content in environment. Colony formation in 20-24 hours.

## METHOD 2

Cells can be plated on 1% agar dish and grown in an anaerobic chamber. Colony formation in 10-12 hours.

# TRANSFORMATION

**NOTE: AMB-1 MUST BE GROWN AEROBICALLY PRIOR TO TRANSFORMATION TO PREVENT MAGNETOSOME FORMATION. MAGNETOSOMES WILL CAUSE CELLS TO DIE DURING ELECTROPORATION. SIMPLY INOCULATE 500 UL CELLS IN 5 ML OF E-MSGM SUPPLEMENTED WITH 50 UL OF IRON MALLEATE AND 50 UL OF VITAMIN SOLUTION IN CULTURE TUBE ON SHAKER AT 30°C. IN 48 HOURS, CELLS WILL BE CONCENTRATED ENOUGH FOR TRANSFORMATION.**

## PARTS FOR VECTORS

### ORI + REP GENE

Size of essential region: 2,973 bp  
Extracted from pMGT.<sup>2</sup>  
Available through Penn iGEM.

### PMSP3

Size: 475 bp  
Promoter available as BBa\_K624016.

### PMSP3 RBS

Size: 31 bp  
Available as BBa\_K624012.

*Promoters of different strengths also available through BioBrick Registry.*

## PREPARING CELLS

Combine 10 tubes (3.5 mL each) of aerobically grown bacteria into a 50 mL falcon tube. Count cells and determine concentration. Must be at magnitude of at least  $10^7$  cells/mL to proceed. Centrifuge at 3700 RPM for 15 minutes. Pour off supernatant. If supernatant is murky, spin for longer. Resuspend pellet with 25 mL of 10 mM TES buffer containing 272 mM sucrose (pH 7.5). Centrifuge at 3700 RPM for 10 minutes. Resuspend pellet in total volume of TES+sucrose buffer needed to concentrate cells to  $10^9$  cells/mL. Split into 50  $\mu$ L aliquots. Store at -80°C.<sup>2,3</sup>

## ELECTROPORATION

Place 0.1-cm electroporation cuvette on ice. Add 1.5  $\mu$ L of DNA to 50  $\mu$ L aliquot of competent cells. Stir gently with pipet. Transfer solution to electroporation cuvette. Add 500  $\mu$ L TES + sucrose buffer to an Eppendorf tube as recovery media. Electroporate with a Gene Pulser® (Bio-Rad Laboratories, Richmond, Calif.) at following settings: Capacitance 25  $\mu$ F // resistance 200  $\Omega$  // 1mm cuvette // 10 kV/cm (1000V). Transfer 55  $\mu$ L of electroporated solution into Eppendorf tube with recovery media. Recover at 30°C in a shaker at 100 RPM for 16 hours before plating.<sup>2,3</sup>

# SAVING CELLS

Inoculate 500  $\mu$ L cells into a sterile 50 mL falcon tube (1:100 dilution) with E-MSGM supplemented with 1% Wolfe's Vitamin Solution and 1% iron malleate. Grow 50 mL culture for 48 hours at 30°C (oxygen depletion in full falcon tube is sufficient). Centrifuge at 3500 RPM for 10 minutes. Make glycerol stock with final glycerol concentration of 30%.

<sup>1</sup> Yan, Lei, Shuang Zhang, Peng Chen, Hetao Liu, Huanhuan Yin, and Hongyu Li. "Magnetotactic Bacteria, Magnetosomes and Their Application." Microbiological Research 167.9 (2012): 507-19.

<sup>2</sup> Okamura, Y., H. Takeyama, T. Sekine, T. Sakaguchi, A. T. Wahyudi, R. Sato, S. Kamiya, and T. Matsunaga. "Design and Application of a New Cryptic-Plasmid-Based Shuttle Vector for Magnetospirillum Magneticum." Applied and Environmental Microbiology 69.7 (2003): 4274-277. Web.

<sup>3</sup> Schultheiss, D., Schüler, D. "Development of a genetic system for Magnetospirillum gryphiswaldense." Archives of Microbiology 179.2 (2003): 89-94. Web. doi: 10.1007/s00203-002-0498-z

# SPECIAL FEATURES

