

Colony PCR

- First alternative:
 - Pick one colony with a sterile tip and elute it in 100 μl ddH₂O or medium
 - Store the colony at 4 °C while colony PCR is running
- Second alternative:
 - Pick one colony with a sterile tip and streak cells at a marked position on a new plate
 - Put tip in PCR tube already containing the reaction mixture
- For both alternatives continue as follows:
 - One reaction mix contains:
 - 5 μl 5x buffer
 - 1 μl MgCl₂ (25 mM stock)
 - 0.5 μl 10mM dNTPs
 - 0.25 μl primer mix (prefix/suffix primers or sequencing primers)
 - 17.625 μl ddH₂O
 - 0.125 μl GoTaq polymerase (Promega)
 - 0.5 μl template
- PCR program:
- Cell lysis and initial denaturation: 5 min, 95 °C
- 30 cycles of:
 - 10 s, 95 °C
 - 30 s, annealing temperature
 - 1 min / 1 kb of expected product, 72 °C
 - Final elongation: 5 min, 72 °C
- Gel electrophoresis: check the fragment size
- First (as above) alternative:
 - Plate the correct colony
- Second (as above) alternative:
 - Use cells from the right positions to start liquid cultures or streak them on a new plate

