

Fast cell lysis for SDS-Page

- Measure OD₆₀₀
 - Transfer two samples a 1,5 ml in a microcentrifuge tube
 - Centrifuge at top speed for 5 minutes at room temperature
 - Discard supernatant quantitative
 - Store pellet at -20 °C
 - Thaw pellet and resuspend in **Fast Cell Lysis sample buffer**. (30 µl **Fast Cell Lysis sample buffer**. per OD₆₀₀)
 - Heat for 5 - 10 minutes at 95 - 100 °C
 - Centrifuge at top speed for 15 minutes at room temperature
 - Transfer supernatant to a new microcentrifuge tube
 - Analyze samples by **SDS-PAGE**. Use 3 - 10 µl per sample

