

Mini Prep

1. inoculation
 - a. touch culture and grow on sterile media
 - b. *if transformation colonies*
 - i. 20 uL of **STERILE** water, touch colony and mix into sterile water
 - ii. put 3.5 mL of LB medium into culture tube
 - iii. put 15 uL of the 20 uL colony solution into the culture tube
 - iv. put on shaker, keep tube at angle for more aeration
 - v. for remaining 5uL, go to colony PCR (optional)
2. centrifuge culture (make pellet)
 - . use 15 mL conical centrifuge tube to not lose cap in centrifuge
 - a. spin down at (3,000 RPM for 5 mins for a good centrifuge) 5,000 RPM for 7 mins on old BE centrifuge at Room temperature
 - b. pour off supernatant
3. remove RNA
 - . take P1 out of fridge, check the cap sticker to make sure RNase is added
 - a. add 250 uL P1 buffer and resuspend (pipet mix)
 - b. after using P1, keep it at 4C in fridge
 - c. transfer to new centrifuge tube 1.5 mL (LABEL!)
4. lyse cells with detergent
 - . add 250 uL P2 buffer
 - a. invert 4-6 times to mix
 - b. CHECK:
 - . DO NOT leave lysis buffer in pellet for more than 5 mins
 - i. viscous
5. ;precipitate unwanted content with neutralization buffer (heavy genome sinks to bottom, plasmid stays on top)
 - . add 350 uL N3 buffer
 - a. invert 4-6 times **vigorously**
 - b. centrifuge 13,000 RPM for 10 mins
 - c. remove tube from centrifuge gently!!, carry rack over
6. filter supernatant to purify plasmid
 - . remove 900 uL supernatant onto column
 - a. Be careful not to get ANY of the white pellet
 - b. centrifuge 13,000 RPM for 1 min
 - c. pour off flow-through in liquid waste
7. wash with alcohol PE buffer to keep plasmid stuck on filter
 - . check the cap sticker to make sure ethanol is added
 - a. add 750 uL PE buffer to column filter
 - b. centrifuge 13,000 RPM for 1 min, no liquid should be left on filter
 - c. pour off Ethanol flow-through in collection tube into liquid waste
8. evaporate alcohol by spin
 - . centrifuge 13,000 RPM for 1 min
9. add Elution Buffer (EB: H2O+ salts) to filter to elute DNA
 - . **move filter over brand new Eppendorf 1.5mL tube**
 - a. add 50 uL EB to filter
 - b. Let SIT ON BENCH on filter for 1 min
 - c. centrifuge 13,000 RPM for 1 min

10. Nanodrop measures DNA concentration at 260nm and 280nm (gives ratio close to 1.8-2.0)
- . 1.5 uL
 - a. wipe and dry cover and sensor with H₂O
 - b. it will ask you what kind of sample? click nucleic acid
 - c. click no for the next question and make sure the arm is down
 - d. Load a blank: EB buffer (1.5 microliters) and small tips
 - e. load it right onto the dot and gently close the lid, click play
 - f. Hit blank
 - g. how will you know it blanked correctly? clean it off and add another drop, treat is as your first sample, it should be a flat line within accuracy of 2
 - h. add actual sample now and record concentration