



Why to do this :

1. To determin the quantity of curli

What you need :

1. Culture media : LB

- 10 g bactotrypton
- 5 g yeast extract
- 5 g NaCl
- 0,5 mL NaOH 10N
- Qsp 1 L

2. Antibiotics concentrations

Chloramphenicol (Cm) : 2 mg/mL

Tetracycline (Tet) : 1 mg/mL

Kanamycin (Kann) : 5 mg/mL

Ampicillin (Amp) : 10 mg/mL

→ 50 µL antibiotic / 5mL medium

3. Apparatus : Agilent 7500cx

How to do :

1. Bacteria culture

- a) Prepare 50mL of liquid culture
- b) Take 20mL of this culture (if you are testing several strains at once, adjust the quantity in order to have the same OD for each strain)

2. Preparing the bacteria

- a) Centrifuge 10 min at 10000g
- b) Add 5mL of ultrapure water
- c) Remove supernatant
- d) Centrifuge 10 min at 10000g
- e) Remove the water (be careful not to dissolve the pellet, do this step one more time if needed)
- f) Add 20mL of Nickel solution at 100uM
- g) Leave them in contact for 30min, on the rotating.
- h) Centrifuge 10 min at 10000g
- i) Remove the non-colored supernatant, centrifuge again if needed until the supernatant is totally removed
- j) Keep the pellet
- k) Measure the OD at 600nm of the supernatant for each strain.

3. Step 2bis

- a) The pullet (from 5 to 15 mg of body weight material) is then mineralized with (1,5ml) nitric acid 69% at 105°C during 3 hours
- b) The mineralization tubes are filled with deionized water to 50 ml
- c) Prior to ICP analysis, sample are diluted from 100 to 10000 times with a mixture HNO₃ 0.5% / In 2% (Indium is used as an internal standard in order to have any signal variation)
- d) ICP-MS analysis were carried out with a calibration curve going from 1 ppb to 100 ppb (1 to 100 µg/L)