



# GUARD

## Protocols



### Assembly of new promoters from oligo's:

Oligo's Klenow Polymerase protocol:

- Make 10 uM stock from 100 uM stock
- Mix the following in a 0.6 mL sterile tube
  - 10  $\mu$ L 10X NEBuffer 2.1 (containing BSA)
  - 4  $\mu$ L oligo 1 10 uM (typically 1  $\mu$ g or more)
  - 4  $\mu$ L oligo 2 10 uM (typically 1  $\mu$ g or more)
  - 75  $\mu$ L deionized sterile H<sub>2</sub>O
- Anneal the two oligos together by either placing the mixture in a thermal cycler ([MJ Research](#), PTC-200) at 94°C for 5 mins, a cool down for 0.1°C/sec to 5°C below the melting temperature of the primers, hold that temperature for 5 mins, then cool down at 0.1°C/sec to 37°C. Alternatively, the tube can be placed in a beaker of boiling water and let cool to room temperature.
- Add 1  $\mu$ L Klenow 3'  $\rightarrow$  5' exo<sup>-</sup> polymerase to mixture. Vortex polymerase before pipetting to ensure it is well-mixed.
- Add 2.5  $\mu$ L 10mM dNTPS (equal to 0.25 mM final concentration of each dNTP). *Recommend using a thermal cycler for the following incubation steps.*
- Incubate 1 hr at 37°C.
- Heat inactivate polymerase by incubating at 75°C for 20 minutes. This inactivation temperature might be higher than the melting temperature of your annealed and extended primers. It may be prudent to ramp the temperature down from 75°C. *See [Restriction Digest](#) for more information on the following steps.*
- Add 1  $\mu$ L EcoRI and 1  $\mu$ L SpeI
- Incubate for a minimum of 2 hrs.
- Heat inactivate restriction enzyme by incubating at 80°C for 20 mins.
- [Purify DNA](#) as necessary and measure concentration with nanodrop

