


Preparation of *B. subtilis* competent cell

Materials

- Sterile 1.5 mL microtubes
- LB broth
- T base – Autoclave
 - $(\text{NH}_4)\text{SO}_4$ 2g
 - $\text{K}_2\text{HPO}_4 \cdot 3\text{H}_2\text{O}$ 18.3g
 - KH_2PO_4 6g
 - Trisodium citrate . $2\text{H}_2\text{O}$ 1g
- SpC – Made fresh on day
 - T base 20ml
 - 50%(w/w) glucose 0.2ml
 - 1.2%(w/w) $\text{MgSO}_4 \cdot 3\text{H}_2\text{O}$ 0.3ml
 - 10%(w/w) Bacto yeast extract 0.4ml
 - 1%(w/w) tryptone 0.5ml
- SpII – Made fresh on day
 - T base 200ml
 - 50%(w /w) glucose 2ml
 - 1.2%(w/w) $\text{MgSO}_4 \cdot 3\text{H}_2\text{O}$ 14ml
 - 10% Bacto yeast extract 2ml
 - 1%(w/v) tryptone 2ml
 - 0.1 M CaCl_2 1ml

Apparatus

- Centrifuge
- Spectrofotometer
- Shaker
- Flow chamber
- Autoclave



Method

1. Streak out the strain to be made competent on an LB agar plate as a large patch and incubate overnight at 30°C;
2. The following morning scrape the cell growth off the plate and use to inoculate fresh, pre-warmed, SpC medium (20 mL) to give an OD₆₀₀ reading of about 0.5;
3. Incubate the culture at 37 °C with vigorous aeration and take periodic OD readings (OD₆₀₀) to assess cell growth;
4. When the rate of cell growth is seen to depart from exponential (i.e. no significant change in cell density over 20-30 min) inoculate 200 mL of pre-warmed, SpII medium with 2 mL of stationary-phase culture and continue incubation at 37 °C with slower aeration;
5. After 90 min incubation, pellet the cells by centrifugation (4,000 g, 25min) at room temperature;
6. Carefully decant the supernatant into a sterile container and save;
7. Gently resuspended the cell pellet in 2,25 mL of the saved supernatant, mix gently.