



Why to do this :

1. To cut plasmids on a specific site (depending on the restriction enzyme used) in order, later, to verify the plasmid or to insert a part

What you need

1. For 2 μL DNA in a total volume of 10 μL
 - 0.5 μL of restriction enzyme
 - 2 μL of DNA
 - 1 μL of buffer 10X
 - 7.5 μL sterile water
 - 2 μL of buffer STOP 6X (30% glycerol, 0,25% BBT, 0,25% cyanol xylen)

How to do :

- a) Add 2 μL of DNA, 1 μL of buffer 10X, 7.5 μL sterile water (to reach a volume of 10 μL) and at the end 0.5 μL of each enzyme.
 - b) Incubate for 1h30 at 37°C.
 - c) Add 2 μL of buffer STOP 6X
-