

E. 9 Entrapping of *E. coli* in alginate beads

16.09.2014

Inoculation of over night cultures of JM 109 C2 Lac-end.

17.09.2014

Measurement of cell number: $5,8 \cdot 10^{10}$ using a Hemocytometer.

Washed in 0,9% saline solution

1 mL of that washed cultures was added to 200 ml of Saline solution. → plating efficiency test (1:10 – 1:1000000) to check how many cells survived the preparation so fare.

2 g alginate into 200 mL water, mixed for a long span of time, add 20 mL of diluted saline culture.

0,05 M CaCl_2 solution was prepared and the alginate culture solution was added → alginate beads formed.

Inoculation of:

Beads in 3 mL 2x YT

Beads in 3 mL rumen liquid

Over night (24h)

18.09.2014

Centrifugation of beads in rumen liquid and 2x YT for 5 min at 800 rpm. Supernatants were transferred into new falcons and centrifuged with *E. coli* in 2x YT and rumen liquid for 10 min at 4000 rpm.

Beads were washed in 0,9% NaCl and incubated in 3 mL 1x PBS for 1h at 37 °C and 250 rpm.

Supernatants were decanted and pellets were washed two times with 3 mL 1x PBS. Pellets were resuspended in 3 mL 1x PBS and then plated on Agar-Plates (Ca & Amp, Dilution 1:1000 & 1:10.000).

Beads were resuspended in 3 mL 1x PBS and plated on Agar-Plates (Ca & Amp, Dilution 1:100 & 1:10.000).

Incubation of plates over night at 37°C.

Results:

19.09.2014

Niels Rüdiger

The experiment from 18.09.2014 is repeated.

20.09.2014

Niels Rüdiger

The experiment from 18.09.2014 is repeated again.

→ Result: