

Use of Escherichia coli genetically modified as a cooper biosensor for educational and research purposes

Risk analysis of Zamorano iGEM's project

Valeria Soledad Guerrero Kesselman y Oliver David Chamorro Ojeda

Presented to Servicio Nacional De Sanidad Agropecuaria, SENASA,
Tegucigalpa, Honduras. September 9, 2014.

Table of Contents

| | |
|--|----|
| Personnel Involved..... | 3 |
| Evaluation Purpose | 3 |
| Genetic material description..... | 4 |
| Transformation methods | 5 |
| Evaluation location..... | 5 |
| Biosecurity measurements..... | 6 |
| Access:..... | 6 |
| Personnel protection:..... | 6 |
| Procedures: | 6 |
| Programed destination | 7 |
| Containment measurements | 7 |
| Final disposition method..... | 8 |
| Annex 1..... | 9 |
| ¿What are the contents of a plasmid? | 9 |
| How to make competent cells?..... | 10 |
| Transformation Protocol..... | 10 |
| Growing conditions of transformed bacteria..... | 11 |
| Plasmid: Copper biosensor..... | 11 |
| Annex 2..... | 12 |
| How does a biosensor based on a microorganism works? | 12 |
| Annex 3..... | 13 |
| Good aseptic practices. | 13 |
| Good pipetting techniques..... | 13 |
| Annex 4..... | 14 |
| Biological waste management. | 14 |
| Bibliography | 16 |

SERVICIO NACIONAL DE SANIDAD AGROPECUARIA (SENASA)

Biotechnology and Biosafety National Committee

Personnel Involved

Name : Zamorano iGEM Team

Students : Oliver David Chamorro Ojeda, Valeria Soledad Guerrero Kesselman.

Supervisor : Dinie Espinal, M.Sc. Maria Mercedes Roca, PhD.

Address : Escuela Agrícola Panamericana Zamorano.

Residencial la Hacienda. Casa #5 Bloque E Calle Pastizales.
Esquina Excel Automotriz. Tegucigalpa Honduras.

Phone Number : 9777-6349/9818-8064

Evaluation Purpose

To introduce the new organism has 2 goals:

A) Educational: As a laboratory module for students from Zamorano University on the learning by doing program. The learning objectives include to build student capabilities on basic bacteriology, aseptic techniques, general principles of biosafety inside the lab, general principles of genetics and molecular biology and to introduce the students on topics of genetic engineering and into synthetic biology as a cutting edge applied science.

B) Research: Biosensor detector of heavy metals in water, as a biological technique to indicate the presence of copper traces in water due to the abuse of agrochemicals. A mechanism of detection it's the first step for the managing and responsible use of such compounds and the preservation of water resources in a specific area.

Both objectives were established due to the feasibility of its usage. The bacteria that changes color in presence of copper becomes a biosensor to detect contamination caused by this metal.

For the developing of the project, there were used bacteria *Escherichia coli* from the non-pathogenic strain DH5 Δ Z1. It was added a plasmid with protein producing genes that changes the bacteria color in presence of copper.

To achieve this objective it was necessary to insert specific promoters for the expression pigmentation genes. As a molecular marker used to detect transformed bacteria, it was used standard antibiotic resistance genes specifically for ampicillin and tetracycline.

Genetic material description

As detailed before, the organism used for the project "Construction of bioscience knowledge" was a non-pathogenic *E coli* of the DH5ΔZ1 strain whose genome was intact.

To achieve the pigmentation change on the bacteria, it was used the plasmid 5006 pd, derived from the pJ8004 plasmid (Figure 1), synthesized by the company DNA 2.0. This plasmid provides resistance to two antibiotics ampicillin (AmpR+) and tetracycline (TerR+) on lab conditions (See annex 1).




Figure 1. Description of the genetic construction showing a pink pigmentation gene.

The plasmid contains: Replication origin site (ColE1 origin), ampicillin (AmpR+) and tetracycline (TetR+) resistance gene, cooper detection sequence (CopA promoter), pigmentation gene (pink, green or fluorescent orange), Rhamnose promotor (F- endA1 glnV44 thi-1 recA1 relA1 gyrA96 deoR nupG Φ 80dlacZΔM15 Δ(lacZYA-argF, U169, hsdR17(rK- mK+), λ-, LacIq+)).

The replication origin site ColE1 present in the plasmid allows approximately 15 to 20 copies of the plasmid for each cell. (copy average).

The bacteria developed with naturally in the absence of cooper, with no phenotypical changes detected in the organism. Once the modified bacteria was in contact with cooper diluted in water, the color changing reaction manifested depending on the pigmentation gene inserted (Pink, green or fluorescent orange).

SYSTEM USED TO PRODUCE THE REGULATED ARTICLE

The pink color derives from the registered gene with trademark CupidPink <https://www.dna20.com/products/protein-paintbox?exp=2> #collapse10 from the company DNA 2.0.

The green color derives from the registered gene with trademark <https://www.dna20.com/products/protein-paintbox?hl=gfp#collapse27> from the company DNA 2.0.

The orange color derives from the registered gene with trademark YukonOFP <https://www.dna20.com/products/protein-paintbox?hl=gfp#collapse44> from the company DNA 2.0.

The gene reactions are under the control of two promoters: the CupidPink gene uses rhaPBad (ramnose active promoter) and the COPA promoter (inducted by cooper), both are necessary for the reaction to take place. The gene is under a constitutive promoter so it is always expressed as long as there is metal present.

Transformation methods

The genetic transformation of the DH5ΔZ1 strain was achieved in collaboration with the Endy lab and Stanford University, Palo Alto, California, EEUU. The genetic material was imported to Zamorano University, Honduras, under an import permit emitted by SENASA from the Honduran Secretary of Agriculture on May 2013.

A standard protocol of bacterial transformation was used with viable cells for a lab exercise with students. The organism practices were conducted in the Zamorano University facilities.

In the annex 2, its detailed information about the transformation method and how the biosensor works.

Evaluation location

Zamorano University, Honduras, Environmental and Development Research Lab.

Biosecurity measurements

The bacteria used it's a biosecurity level I (This category demands less restrictive biosecurity measurements). To guarantee the students and the personnel safety who worked on the modification and test of the bacteria, the following biosecurity were taken:

Access:

Only the authorized personnel was allowed entrance to the laboratory.

The laboratory doors remained closed during the experiment.

It was not allowed the entrance to kids nor animals on the laboratory facilities.

Personnel protection:

1. At all times there were used lab coats while working on the lab.
2. There were used protection gloves for all the procedures that may involve direct or accidental contact with blood, body fluids and any other potentially infectious materials. Once used the gloves were removed aseptically followed by a hand washing.
3. All personnel had to wash their hands after handling lab material, as well as when leaving the working area.
4. There were used security goggles, face protectors and other devices when needed to protect the eyes from spilling, impacts or UV radiation.
5. The laboratory jackets and the rest of the protection equipment were exclusively used for authorized areas inside the lab in order to maintain the aseptic integrity of the practice.
6. There were used toe protected shoes.
7. In the working areas there was forbidden to eat, drink, smoke, use of make up or handling contact glasses.
8. It was forbidden to keep food or drinks for human consumption in the working areas of the lab.
9. The protection clothing was stored on a different storage than the regular use clothing.

Procedures:

1. It was strictly forbidden to create vacuum during pipetting using the mouth.
2. There was no material in contact with the mouth.

3. All technical procedures were conducted in order to reduce at minimum the formation of aerosols and fine drops.
4. It was mandatory to let the lab supervisor know in case of a spill, accident or real exposition or potential to infectious materials. During the development of the project, none of the situations cited above occurred.
5. All contaminated fluids (chemically or physically) were decontaminated before sent to the sewer. Zamorano University has its own treatment mechanism for waste water flow inside the campus as detailed in the bio-contention normative.
6. The written documents planned to go out from the lab, were protected from contamination while in the lab.

For further information regarding the procedure, please check annex 3 (Good aseptic/pipetting practices).

Programed destination

The bacterium was handled completely isolated inside the lab. It was developed on solid medium using petri dishes with nutritious agarose agar and liquid medium using controlled growth on solutions.

During storage the bacteria was keep on exclusively and special refrigerators, hermetically closed in temperatures from 5-8 C° with 45.50% of humidity so the bacteria can enter in a controlled latency status. at the moment of use, there were removed and exposed to the lab conditions (25C°, 50-70% humidity) during a few moments, just to extract the samples and then put back on cold storage in the same medium.

The lab is innocuous, there was not detected the presence of animals or insects that may consist of contamination vectors and spreading of the bacteria. For its manipulation it was required the use of gloves and tweezers to avoid direct contact between the investigator or the student with the bacteria.

Containment measurements

1. The laboratory was clean and free of non-related to the experiment materials at all times.
2. In the eventuality of accidental spill of potentially dangerous materials the procedure to follow was the decontamination of the area using alcohol and chloride. There were no incidents as mentioned during the experiment.
3. All the instruments, samples and growth medium containing biologic material must have been decontaminated before being

discarded. In the need of re use the material, it was clean and previously sterilized.

4. The windows that could be opened were equipped with filters making impossible the access of arthropods.

Final disposition method

Once the experimental procedures with the bacteria in study were concluded, this ideal environment for the developed organism was destroyed adding isopropyl alcohol to the probing tubes and petri dishes that were used. The addition of the compound killed all the microorganisms present on the growth medium due to cell membrane disintegration and consequent denaturalization of proteins inside them.

After this all waste was put together in special containers separating the liquid waste from the solid waste. Next there were eliminated along with the University's medical clinic waste.

Finally all the instruments used were sanitized with chlorine water and brushes to maintain the asepsis inside the lab. Then the instruments were exposed to high temperatures and pressure inside the autoclave.


The correct handling of bio-waste through applied bio containment procedures was one of the key factors to ensure the fulfilling of the biosecurity norms established and the integrity of all participants. In annex 4 there is more information about the handling of biological waste.

Annex 1

¿What are the contents of a plasmid?

The picture below illustrates the components and generalities of a plasmid.

Que es un plásmido?




- ADN circular
 - No es parte del cromosoma
 - Se replica autónomamente
 - Ocurre de forma natural
- Ej. *Agrobacterium tumefaciens* y transformación de cultivos Bt

Cellular preparation (bacteria) components containing plasmid (DNA)


Below there is a graphic explanation of a competent bacteria and the methodology of its transformation with lab procedures.

“Bacterial Competence” y transformación

“Competence”: habilidad para absorber ADN foráneo



How to make competent cells?



Transformation Protocol

To achieve a correct genetic transformation, a series of steps are follow as described below:

Protocolo de transformación genética

- Rotule los tubos y platos (3X) con el nombre del grupo .
- Mezcle 100 μ L 2xTSS con 100 μ L células bacterianas (E.coli)
- Ponga el tubo en hielo
- Adicione 10 μ L del plásmido (ADN) a las células
- Mezcle breve y suavemente invirtiendo el tubo
- Incube en hielo por 5 min
- Adicione 400 μ L medio LB

Condiciones de incubación :

Esparza 100 μ L de las células con las esferas de vidrio (suavemente de arriba a abajo)


Rotule los 3 platos petri • 15 min

- Plato con antibiótico
- Plato con cobre
- Plato con antibiótico + cobre

Growing conditions of transformed bacteria.

The bacterial cells once transformed, need to habit a special culture medium for its correct and controlled development. On the following image it's possible to see the best conditions for its

Condiciones de crecimiento de bacterias transformadas




Alternativa : platos petris con kanamicina en lugar de ampicilina

growth.


Plasmid: Copper biosensor

On the following image are shown the characteristics of the plasmid used in the project, this parts guarantee the color change of the bacteria under the presence of copper dissolved in water.

Plásmido: biosensor de cobre



Genetic transformation of plants in classic biotechnology



The bacteria *Agrobacterium tumefaciens* is a genetic vector widely used in molecular biology. *A. tumefaciens* is a common pathogenic bacteria that attacks solanaceae plants as tomato, causing scars. It has a natural plasmid able to introduce the bacterial DNA. The plasmid Ti, contains a gene that causes a tumor (Ti = Tumor inducing) to the infected plant cells. In the lab it is used a strain of *A. tumefaciens* with a Ti plasmid without the tumor gene.

Annex 2

How does a biosensor based on a microorganism works?

The plasmids are widely used on molecular biology for transformations on classic biotechnology and also in synthetic biology. In classic biotechnology genes from a living being are inserted with a sequence already existing in nature. In synthetic biology, the plasmids can also be used as genetic vectors to introduce fragments of artificial DNA (artificial genes) synthesized from chemical compounds (nucleotides A, C, G and T) in a lab. The information from this synthetic genes is obtained from bioinformatics data and also from the creativity of a gene designer. In synthetic biology "standard parts" are used following the basic principles of engineering. With this artificial or synthetic genes, microorganism genomes can be programmed to have a metabolism and function different from natural microorganisms. The transformation technique with a plasmid is the same in classic biotechnology as in synthetic biology.

Annex 3

Good aseptic practices.

Bacteria are everywhere. Microbiology it's a field of biological science with applications in many fields such as human medicine and veterinary, agriculture, environmental science, industry and bioenergy production. Learning how to properly work applying asepsis it's the basis of microbiology.

To avoid the contamination of petri dishes and probing tubes:

1. Team learned how to autoclave equipment and reactivities as well as growing media.
2. Minimize the time lab equipment its exposed to open air.
3. It was used sterilized equipment and surfaces (by fire, autoclaving or chloride).
4. If a sterile equipment was set in contact with a non-sterilized surface it was discarded immediately.
5. It was used a new pipette tip every time.

Good pipetting techniques

Push softly trough the pipet's end.

Push the plunger until the stop position #1.

Insert the tip on the liquid.

Softly let go the plunger to the initial position.

Take the tip out form the liquid and make sure there is no liquid outside the tip nor bubbles inside the tip.


Insert the tip in the probing tube or the deposit where the liquid is going to be.

Push the plunger to stop position # 1. (Most of the content may come out).

Continue pushing the plunger to stop position # 2 (blow out).

Take of the tip from the liquid and softly let the plunger go until its original position.

Micropipetas



Annex 4

Biological waste management.

Even if the *E. coli* strain used for this project requires the biosafety level number 1, the biological waste will be handled as if it was the highest biosafety level (4) to illustrate and learn the biological compounds waste practices.

Note: The biosafety scale goes from 2 to 4.

As an example: The Ebola virus has a biosafety level of 4 and requires laboratories, equipment and highly specialized personnel for its management. The tomato mosaic virus has a biosafety level of 1 for humans (since it is not pathogenic) but the skin contact must be avoided if the person has tomato related activities.

It is not necessary to use lab coat nor gloves in this lab.

It is not allowed the consumption of food or drinks in the lab.

Avoid to take objects to the mouth, such as pencils or fingers.

Do not take biological material outside the lab.

Immediately inform the personnel in charge of the lab about any accident or chemical reactive spill in the lab.

Handle the bacterial culture with caution, avoiding placing bacteriological loops or pipette tips on the lab seat. Sterilize bacteriological loops using fire and discard the pipette tips or any other disposable object that may have been in contact with bacterial culture.

Use autoclave ready bags for biological waste such as petri dishes and other materials like pipette tips, and Eppendorf tubes.

Autoclave biological waste in specialized autoclaving bags.

Clean all surfaces on the work area using disinfectant.

Add chloride in the liquid bacterial cultures when no longer in use, if not possible, the procedure is to autoclave the growing media and the glassware immediately.

Carefully wash the hands with water and soap when the practice it's done.

Bibliography

ArgenBio. (9 de Diciembre de 2010). *Consejo Argentino para la información y el Desarrollo de la Biotecnología*. Obtenido de <http://www.argenbio.org/index.php?action=notas¬e=5467>

Collado, J., & Moreno, G. (nd). *Escherichia coli en Estudios Post-Genómicos Microbianos*. Centro de Investigación sobre fijación de Nitrógeno. Recuperado el 2014, de <http://www.biblioweb.tic.unam.mx/libros/microbios/Cap19/>

Colwell, R. (2009). *Understanding Biosecurity protecting against the misuse of science in today's world*. Washington: National Research Council of the National Academies.

Galli, L. (2012). *Estudio de los factores de adherencia de cepas de Escherichia Coli productoras de toxina Shiga aisladas de bovinos*. Universidad Nacional de la Plata, Facultad de ciencias veterinarias. Buenos Aires: SEDICI.

Kido et al. (2012). *Guía para la Evaluación de Riesgo Ambiental de Organismos Genéticamente Modificados*. Sao Paulo: International Life Sciences Institute do Brazil.

Lara, A. (2011). *Producción de proteínas recombinantes en Escherichia coli*. Universidad Autónoma Metropolitana - Cuajimalpa, Departamento de Procesos y Tecnología . México DF: Scielo. Obtenido de http://www.scielo.org.mx/scielo.php?pid=S1665-27382011000200006&script=sci_arttext

Ledermann, W. (2007). *Una historia personal de las bacterias*. Santiago de Chile: Sociedad chilena de infectología. Obtenido de http://books.google.hn/books?id=pydkwNBz9dEC&pg=PA247&lpg=PA247&dq=theodore+von+escherich&source=bl&ots=TMhMpVY-Pb&sig=jYB-amqknupEWe_6LcHVsHu85AU&hl=es&sa=X&ei=7xi2U72yA-qgsQT-zID4Aw&sqi=2&ved=0CHEQ6AEwCQ#v=onepage&q&f=false

LUTZ R, Bujard H: Independent and tight regulation of transcriptional units in *Escherichia coli* via the LacR/O, the TetR/O and AraC/I1-I2 regulatory elements. *Nucleic Acids Res* 1997, 25:1203-1210

Martinez, A., & Gosset, G. (2007). *Ingeniería metabólica de bacterias*. México D.F. Obtenido de <http://es.scribd.com/doc/232538278/Ingenieria-Metabolica-de-Bacterias>

MCLAB. (2014). *Molecular Cloning Laboratories*. Obtenido de <http://www.mclab.com/Dh5-Alpha-Competent-E.-Coli.html>

Microbe Wiki. (2014). Obtenido de [https://microbewiki.kenyon.edu/index.php/DH5-Alpha E.coli](https://microbewiki.kenyon.edu/index.php/DH5-Alpha_E.coli)

Ministerio de Salud de Chile. (2013). Obtenido de División de Planificación Sanitaria.

Moredo, F. (2012). *Prevalencia de Escherichia coli enterotoxigénico y Escherichia coli productor de toxina Shiga en cerdos sin manifestación clínica de diarrea de la provincia de Buenos Aires*. Universidad Nacional de la Plata, Facultad de Ciencias Veterinarias. Buenos Aires: SEDICI.

OMS. (2005). *Manual de bioseguridad en el laboratorio* (Tercera edición ed.). Ginebra: Organización Mundial de la Salud.

Secretaría del Convenio sobre la Diversidad Biológica. (2000). *Protocolo de Cartagena sobre la Seguridad de la Biotecnología del Convenio sobre la Diversidad biológica: textos y anexos*. Montrea: Secretaría del Convenio sobre la Diversidad Biológica.

SENASA. (2014). *Gobierno de la República de Honduras - Secretaría de agricultura y ganadería*. Obtenido de <http://www.senasa->

sag.gob.hn/index.php?option=com_content&task=view&id=101&Itemid=56
1

ORGANIZACIÓN MUNDIAL DE LA SALUD (OMS). 2005. Manual de bioseguridad en el laboratorio. Edit. Organización Mundial de la Salud. Ginebra. 210p