

## Protocol for double digestion (50μl system)

Pipette the following into a 1.5ml microfuge tube:

Enzyme A	2μl
Enzyme B	2μl
10× buffer	5μl
DNA	0.5-1μg
dd H <sub>2</sub> O	rest of the volume

incubate at recommended temperature (37°C) for at least 2 hour;

Purify the digestion product;

### Notes:

The enzymes used here are TaKaRa enzymes (*EcoRI*/ *XbaI*/ *SpeI*/ *PstI*), and buffer M or H is suitable for most of double digestion;

For 20μl reaction system, the suggested amount of each restriction enzyme is 1μl;