


Protocol for small segments synthesis by PCR


Description:

There are some parts we cannot find in the Kit Plates and whose length is too short to get them from purification. So, we designed a method to synthesize this kind of small segments: design 4 primers which have almost 20 base pairs overlapping sequence and conducted twice PCR, and finally put the segment into a vector by double digestion and ligation.

The whole gene:



4 designed primers:



First PCR:

Mixture: primer 1, primer 2, primer 3, primer 4, dNTP, 10× Buffer;
Pfu DNA polymerase;
T_m : 35°C~50°C;
Circulation: 22;

Gel Purification with AxyPrep DNA Gel Extraction Kit;

Second PCR :

Template: product of first PCR;
Primer 1, primer 4;
Pfu DNA polymerase;
dNTP;
10× Buffer;
T_m : 60°C~65°C;
Circulation: 34;

Gel Purification with AxyPrep DNA Gel Extraction Kit;

Double enzyme digestion and ligation

Note:

Length of primers: 32~36bp

Length of overlapping sequence: 17~22bp