

QuikChange Site Directed Mutagenesis

Primer Design Guidelines

- Both of the mutagenic primers must contain the desired mutation and anneal to the same sequence on opposite strands of the plasmid.
- Primers should be between 25 and 45 bases in length, with a melting temperature (T_m) of $\geq 78^\circ\text{C}$. Primers longer than 45 bases may be used, but using longer primers increases the likelihood of secondary structure formation, which may affect the efficiency of the mutagenesis reaction.
- The following formula is commonly used for estimating the T_m of primers:
$$T_m = 81.5 + 0.41(\%GC) - (675/N) - \% \text{ mismatch}$$
 - N is the primer length in bases
 - values for **%GC** and **% mismatch** are whole numbers
- For calculating T_m for primers intended to introduce insertions or deletions, use this modified version of the above formula:
$$T_m = 81.5 + 0.41(\%GC) - (675/N)$$
where N does not include the bases which are being inserted or deleted.
- The desired mutation (deletion or insertion) should be in the middle of the primer with ~ 10 – 15 bases of correct sequence on both sides.
- The primers optimally should have a minimum GC content of 40% and should terminate in one or more C or G bases.

PCR Reaction

- Use 125 ng of each primer. To convert nanograms to picomoles of oligo, use the following equation:

$$X \text{ pmoles of oligo} = (\text{ng of oligo}) / (330 \times \# \text{ of bases in oligo}) \times 1000$$

For example, for 125 ng of a 25-mer:

$$(125 \text{ ng of oligo}) / (330 \times 25 \text{ bases}) \times 1000 = 15 \text{ pmole}$$

- Use standard Phusion PCR protocol with following modifications:
 - (i) elongation time ~ 1 minute for 1 kb
 - (ii) 12 cycles
 - (iii) Annealing temperature 60°C

It usually works well to try different template DNA concentrations (e.g. 5, 10, 20 and 50 ng).

As a control, prepare a reaction without Phusion (should give no colonies)

DpnI digestion

1 µl DpnI/PCR reaction
Incubate 60 min at 37°C

***E. coli* transformation**

According to a standard protocol, with 10 µl PCR reaction

Protocol generously provided by the lab
Prof. Thorsten Mascher
Großhadernerstr. 2-4
82152 Planegg-Martinsried
www.syntheticmicrobe.bio.lmu.de