

## Date: 7/30/13 People in lab: Emily Puleo, Kelsey Crossen

---

**Title:** Preparation of *P. aeruginosa*

**Start Time:** 3:00 PM

**Purpose:** To make a stock of *P. aeruginosa* for future lab work.

**Protocol:** Inoculated LB broth with *P. aeruginosa* from plated bacteria, no antibiotic. Inoculated another tube of LB broth with liquid *P. aeruginosa* culture, 10ml into 10ml.

Left both cultures in 37C incubator overnight.

Spun down 1ml cultures of *P. aeruginosa*, discarded supernatant, and re-suspended in 1ml Milli-Q water.

Placed cell solutions in 1.5ml Epi tubes and placed in -20C freezer.

**Products:**

Sample Label	Description	Source Label	Quantity
<i>P. aeruginosa</i> 7/30	<i>P. aeruginosa</i> culture	<i>P. aeruginosa</i>	2
<i>P. aeruginosa</i> 7/30	<i>P. aeruginosa</i> culture, frozen	<i>P. aeruginosa</i>	10

**Stop Time:** 5:00 PM

**Next:** Amplify DNA segments from prepared *P. aeruginosa* samples.

## Date: 7/31/13 People in lab: Emily Puleo, Kelsey Crossen

---

**Title:** PCR of norCB, dilution of norCB primers

**Start Time:** 4:57 PM

**Purpose:** PCR norCB from *P. aeruginosa*, adding iGEM prefix and suffix; to prepare primers for future lab work

**Protocol:** LTM ed. 2 pg. 48 **Exceptions:** 1. Less than 0.1 uM final primer concentration.

**Products:**

Sample Label	Description	Source Label	Quantity
7/31/13 A norCB 1	Amplified norCB, with prefix/suffix, 7.5uL template	<i>P. aeruginosa</i> 7/30	1
7/31/13 A norCB 2	Amplified norCB, with prefix/suffix, 5.63uL template	<i>P. aeruginosa</i> 7/30	1
7/31/13 A norCB 3	Amplified norCB, with prefix/suffix, 3.75uL template	<i>P. aeruginosa</i> 7/30	1
7/31/13 norCB F	Forward primer for norCB, 0.73 uM	norCB Forward Primer	1
7/31/13 norCB R	Reverse primer for norCB, 0.87 uM	norCB Reverse Primer	1

**Results:** Sample 3 had too little product and was discarded.

**Notes:** Different amounts of template used to find optimum amount needed.

Forward primer:  $(24.2 \cdot 10^{-9}) / (0.33 \cdot 10^{-3}) \text{ mol/L} = 73 \text{ uM}$ ,  $(73 \text{ uM})(10 \text{ uL}) = (0.73 \text{ uM})V_2$ ,  $V_2 = 1000 \text{ uL}$

Reverse primer:  $(28.7 \cdot 10^{-9}) / (0.33 \cdot 10^{-3}) \text{ mol/L} = 87 \text{ uM}$ ,  $(87 \text{ uM})(10 \text{ uL}) = (0.87 \text{ uM})V_2$ ,  $V_2 = 1000 \text{ uL}$

Used Dr. Shannon's thermocycler, copied "igem\_pcr" program

**Stop Time:** 9:00 PM

**Next:** Run samples on a gel, extract the stronger band from between them, sequence the result.