

October 7, 2014

**Continue testing sensor**

- Transform plasmids into Top10 Cells
- Overnight cultures - dilute 100  $\mu$ l in 2 ml
  - Grow at 37 °C for 1.5 hours
  - Add 1  $\mu$ l of beta-estradiol in ethanol (10 mg/ml)
    - cultures treated with  $\sim$ 20  $\mu$ M of estrogen
- Incubate at 30 °C
- Grow overnight

October 8, 2014

**Continue testing sensor**

- Set up cultures of WT sensor, DEAD Sensor, 115 YFP, 115 YFP and T7RFP term, T7RNAP YFP and T7RFP term, and Top10 cells only at 37 °C, 30 °C and 20 °C to grow overnight

**Transformation of Killer Red Codon Optimized and Supernova Codon Optimized into Top10 Cells**

- 100 µl of cells for remaining DNA (about 1 - 3 µl)
- Plate 400 µl transformants on LB + Cam plates

iGEM Transformation Protocol: <http://parts.igem.org/Help:Protocols/Transformation>

October 9, 2014

**Continue testing sensor**

- Cells at 20 °C grew too slowly to be read
- Cells at 30 °C and 37 °C read in TECAN at yellow and red wavelengths

**Table 1. TECAN Parameters used to Read Samples**

<b>Fluorescent Protein</b>	<b>Excitation Wavelength (nm)</b>	<b>Emission Wavelength (nm)</b>	<b>Gain (if not optimal)</b>
YFP	514	527	90
RFP	584	607	175

Set up cultures of 115 YFP, 115 YFP and T7RFP term, T7RNAP YFP and T7RFP term in triplicate at 37 °C, 30 °C and 20 °C to grow overnight