



Why to do this :

1. Insert plasmids notably into NM522 strains, or any other E.Coli strain.
2. It has been found more efficient than the CaCl_2 transformation protocol, thus we recommend using it.

What you need :

1. Culture media : LB

- 10 g bactotrypton
- 5 g yeast extract
- 5 g NaCl
- 0,5 mL NaOH 10N
- Qsp 1 L

2. Antibiotics concentrations

Chloramphenicol (Cm) : 2 mg/mL

Tetracycline (Tet) : 1 mg/mL

Kanamycin (Kann) : 5 mg/mL

Ampicillin (Amp) : 10 mg/mL

→ 50 μL antibiotic / 5mL medium

3. Material :

- 100 μL of TSS
- 1 μL of DNA
- 900 μL LB medium
- Ice
- Heating bath (42°C)
- 5mL of NM522 cell culture

How to do :

1. Cell culture

- a) Monitor the OD600 of a 5mL culture of NM522 cells (or the strain you want to transform) in LB medium at 37°C. Proceed to the next step when it reaches 0.3-0.4 (not higher than 0.6), which takes approximately 3 hours.
- b) Take 1 mL from the previous culture and centrifuge at 8000 rpm for 2mn

2. Transformation

- a) Resuspend into 100 μL of TSS. Starting from this step, keep the bacteria on ice.
- b) Incubate on ice for 5-10mn.
- c) Add 1 μL of the DNA to be transformed
- d) Incubate on ice for 10mn at least
- e) Heat shock the bacteria by heating them at 42°C for 50s
- f) Incubate on ice for 2mn
- g) Add 900 μL LB and mix by inverting the tube
- h) Incubate for 1h at 37°C
- i) Plate 100 μL on LB medium with the appropriate antibiotic