

Lab Notebook

Beta System

1C3+1K3 preparation

05/08/2014

Primer name	forward/ reversed	Tm	Add DMSO?
Suffix f	forward	61	yes
Prefix r	reversed	62	yes

In the distribution the concentration 25 ng/μl

For final concentration of 2 ng/μl add 4 μl to 46 μl UPW

Reaction mix:

component	Volume[μl]	Mix x2
phusion reaction buffer(x5)	10	20
dNTPs(10 mM)	1	2
Suffix f -forward primer	2.5	5
Prefix r -reverse primer	2.5	5
plasmid (2 ng/μl)	5	-----
Phusion hot start II	0.5	1
DMSO	1.5	3
UPW	27	54
tot	50	-----


PCR program

stage	Temp C ⁰	time
Initial denaturation	98	30 sec
35 cycles	98	10 sec
	61	30 sec
	72	1.20 min
Final extension	72	10 min
hold	4	

Part name	1C3 (lin)	1K3 (lin)
date	05/08/2014	05/08/2014
concentration	49 [ng/μl]	57 [ng/μl]
-20 °C storage	iGEM B sys lin.S	iGEM B sys lin.S

Run gel agarose 1% for product validation

Expected product: 1C3: 2070 bp, 1K3: 2204 bp



Blunt ligation

08/08/2014

Reaction content	1C3 49[μl]	1K3 57[μl]
T4 ligase buffer(x10)	2.5	2.5
T4 kinase (PNK)	1	1
Plasmid [ng/μl]	2	2
PEG 4000 (50%)	2.5	2.5
DDW	17	17
total	25	25

Incubate at 37 °C for 1 hr starting at 14:20

Inactivation 65 °C for 20 min

Ligation started at 15:50 for 2 hr at 25 °C in dry both.

Storage in -20 °C

Part name	1C3 (cir)	1K3 (cir)
date	08/08/2014	08/08/2014
-20 °C storage		

Transformation to top 10

09/08/2014

20 min frost on ice 100 mL
 5 mL of Gibson product (cir)
 30 min on ice
 1 min heat shock 42 °C
 2 min on ice


incubation overnight.

T7 poly preparation

In order to accommodate the T7 polymerase with type 2 Gblocks , we will design 2 PCR reaction that will build Gib3-T7 RNA polymerase-Gib4-Gib2 by 3 PCR reaction
 Then the linear construct will be verified by sequencing

24/07/2014

Incubation of plate with Top 10, pGSV T7 egfp o\n (from Inbal GS)




25/07/2014

starter incubation in shaker 37 oC 250 rpm o\n at 16:40

26/07/2014

09: 46 GS has been made
10:18 MP has finished



Part name	pGSV T7 egfp (cir)
date	26/07/2014
concentration	703 [ng/ μ l]
-20 °C storage	iGEM B system \C1
Part name	pGSV T7 egfp (cir)
date	26/07/2014
Cell type	Top 10
resistance	AMP
-80 °C storage box name	iGEM 2014 B
Location	B1

Oligo PCR

16/09/2014

name	seq	length	Tm1	Tm2	neb	%GC
Oligo T7 f	CCA ATA TCC TTA GCT GAT CAC CAT TAA AGA GGA GAA ACC ACC AAT GAA CAC GAT TAA CAT CGC TAA GAA CG	71	56	68	St.1:66 St.2:72	40
Oligo T7 r	CCG TGT GAT CAA TTA AGG ATT CCC TTC AAA TTA ACA GGA ATC GGT TAC GCG AAC GCG AAG TCC GAC	66	62	70	St.1:66 St.2:72	47

PCR mix

component	Volume[ml]	In?
phusion reaction buffer(x5)	10	
dNTPs(10 mM)	1	
Oligo T7 f (forward primer)	2.5	
Oligo T7 r (reverse primer)	2.5	
T7 poly st.3 (2ng/ml)	5	
Phusion hot start II	0.5	
DMSO	1.5	
UPW	27	
tot	50	


PCR program

stage	Temp C ⁰	time
Initial denaturation	98	30 sec
5 cycles	98	10 sec
	66	30 sec
	72	1:30 min
30	98	10 sec
	72	30 sec
	72	1:30 min
Final extension	72	10 min
hold		

Part name	3,2 T7 poly (lin)
date	16/9/2014
concentration	228 [ng/ml]
-20 ⁰ C storage	iGEM B system

Run gel agarose 1% for product validation

Expected product: T7 RNA Poly: 2670 bp

**PCR**

06/08/2014


Primer name	forward/ reversed	Tm	Add DMSO?
Suffix r	forward	61	yes
Prefix f	reversed	62	yes

component	Volume[ml]	Mix x6
phusion reaction buffer(x5)	10	60
dNTPs(10 mM)	1	6
Suffix f -forward primer	2.5	15
Prefix r -reverse primer	2.5	15
template (2 ng/ml)	5	--\--
Phusion hot start II	0.5	--\--
DMSO	1.5	9
UPW	27	162
tot	50	

PCR program

stage	Temp C ⁰	time
Initial denaturation	98	30 sec
35 cycles	98	10 sec
	61	30 sec
	72	40 sec
Final extension	72	10 min
hold	4	

Run gel agarose 1%

**Gibson assembly**

08/08/2014

Insertion of plux-ter to pSB1C3 by Gibson assembly
Calculation of the concentrations:

$$\left[\frac{pmol}{\mu l} \right] = \frac{\left[\frac{ng}{\mu l} \right] \times 10^3}{size \times 650}$$

Reaction mix:

name	size (bp)	ng/ μ l	pmol/ μ l	Gibson assembly master mix(x2)	pSB1C3 (pmol)	pSB1C3 (ng/ μ l)	Gblock (0.15pmol)	pSB1C3 (0.05pmol)	UPW	total
Pconst-LuxR	1114	84	0.116	15	0.036	49	1.3	1.4	2.3	20
PT7-LuxI	928	84	0.139	15	0.036	49	1.1	1.4	2.5	20
PCI-LacI	1464	70	0.07	15	0.036	49	2.1	1.4	1.5	20
Ptet-tetR	1001	83	0.127	15	0.036	49	1.2	1.4	2.4	20
Plac-Immbda CI	1149	97	0.129	15	0.036	49	1.2	1.4	2.4	20
POmpC-amilCP	1097	34	0.047	15	0.036	49	3.2	1.4	0.4	20

Incubation at 50 $^{\circ}$ C for 1 hr.

Plux-G.Luc

21/08/2014

PCR

Primer name	forward/reversed	Tm	Add DMSO?
Suffix r	forward	61	yes
Prefix f	reversed	62	yes

-

component	Volume[μ l]
phusion reaction buffer(x5)	10
dNTPs(10 mM)	1
Suffix f -forward primer	2.5
Prefix r -reverse primer	2.5

Plux-green luciferase (2 ng/μl)	5
Phusion hot start II	0.5
DMSO	1.5
UPW	27
tot	50

-
PCR program

stage	Temp C ⁰	time
Initial denaturation	98	30 sec
35 cycles	98	10 sec
	61	30 sec
	72	40 sec
Final extension	72	10 min
hold	4	

-20 storage:

name	size (bp)	date	ng/μl	iGEM B sys lin.s
Plux-green luciferase	1314	21/08/2014	166	

Gibson assembly

Insertion of plux-ter to pSB1C3 by Gibson assembly

Calculation of the concentrations:

$$\left[\frac{pmol}{\mu l} \right] = \frac{\left[\frac{ng}{\mu l} \right] \times 10^3}{size \times 650}$$

Reaction mix:

name	size (bp)	ng/ml	pmol/ml	Gibson assembly master mix(x2)	pSB1C3 (pmol)	pSB1C3 (ng/ml)	Gblock (0.15pmol)	pSB1C3 (0.05pmol)	UPW	total
Plux-green luciferase	1314	166	0.19	15	0.036	49	0.8	1.4	2.8	20

Incubation at 50 °C for 1 hr.

Store in -20 °C

Part name	Plux-green luciferase
date	
-20 °C storage iGEM B system	

Colony PCR

Primer name	forward/reversed	First Tm	Tm	Add DMSO?
prefix f	forward	50	63	yes
suffix r	reversed	50	63	yes


Reaction mix: make total mix and then divide 20 [μl] to each PCR tube

component	Volume[μl] for 1 colony	For total mix 6 reaction	In?
Taq ready mix (x2)	10	60	V
DMSO	0.5	3	V
prefix f (10 ng/μl)	1	6	V
suffix r (10 ng/μl)	1	6	V
colony	colony	----	
UPW	7.5	45	V
tot	20		

PCR program

stage	Temp C ⁰	time
Initial denaturation	94	3 min
35 cycles	94	30 sec
	57	30 sec
	72	50 sec
Final extension	72	10 min
hold	4	

Run gel agarose 1%, expected product : 1314 bp



Colony PCR

28/08/2014


Primer name	forward/ reversed	First Tm	Tm	Add DMSO?
prefix f	forward	50	63	yes
suffix r	reversed	50	63	yes

Reaction mix: make total mix and then divide 20 [ml] to each PCR tube

component	Volume[ml] for 1 colony	For total mix 32 reaction	In?
Taq ready mix (x2)	10	320	V
DMSO	0.5	16	V
prefix f (10 ng/ml)	1	32	V
suffix r (10 ng/ml)	1	32	V
colony	colony	----	
UPW	7.5	240	V
tot	20		

code	name	length
A	Pcat-LuxR	1114
B	PT7-LuxI	928
C	PompC-AmilCP	1097
D	PCI-LacI	1464
E	Plac-CI	1149
F	Ptet-tetR	1001

A: Pcat-LuxR, expected product: 1114 bp




Although we have positive result, it seem that the colonies are unpure or the PCR T_m was Low. In order to support the positive result restriction reaction was made.

Reaction mix:


	Volume [ml]	in
1C3:Pcat-LuxR (114 [ng/ml])	5	V
EcoRI HF	0.5	V
Cutsmart buffer (x10)	2	V
MBW	12.5	V
total	20	+

Run agarose gel 1%, positive : 3000 bp negative: 2000 bp




Colony 1 has been selected to be sent to sequencing

B: PT7-LuxI, expected product: 928 bp




Colony 3 has been selected to be sent to sequencing

C: PompC-AmilCP expected product: 1097 bp




Colony 3 has been selected to be sent to sequencing

D: PCI-LacI expected product: 1464 bp




Colony 1 has been selected to be sent to sequencing

E: Plac-CI expected product: 1149 bp



Colony 4 has been selected to be sent to sequencing

F: Ptet-tetR expected product: 1001 bp



Colony 1 has been selected to be sent to sequencing

Colony PCR- 2

01/09/2014


Primer name	forward/ reversed	First Tm	Tm	Add DMSO?
prefix f	forward	50	63	yes
suffix r	reversed	50	63	yes

Reaction mix: make total mix and then divide 20 [ml] to each PCR tube

component	Volume[ml] for 1 colony	For total mix 25 reaction	In?
Taq ready mix (x2)	10	320	V
DMSO	0.5	16	V
prefix f (10 ng/ml)	1	32	V
suffix r (10 ng/ml)	1	32	V
colony	colony	----	
UPW	7.5	240	V
tot	20		


code	name	length
A	PompC-AmilCP	1097
B	Pcat-LuxR	1114
C	PT7-LuxI	928
D	PCI-LacI	1464
E	Plac-CI	1149

A: PompC-AmilCP expected product: 1097 bp




Unspecific band in all colony

B: Pcat-LuxR, expected product: 1114 bp




Col 4 seem to be positive and clean

C: PT7-LuxI, expected product: 928 bp




Unspecific band in all colony

D: PCI-LacI expected product: 1464 bp



Colony 1 and colony 2 contains Unspecific band with the positive band
Colonies 3,4 are negative


E: Plac-CI expected product: 1149 bp



Colony 2 seem to be positive and clean

Restriction validation

28/08/2014



Restriction validation


03/09/2014

code	name	colony	Status
A	PompC-amilCP	1	x
		2	x
B	Pcat-LuxR	2	x
		4	V
C	PT7-LuxI	2	x
		3	x
D	PCI-LacI	1	x
		2	x
E	Plac-CI	2	V
		4	V

Cut with 2 restriction enzyme:

NcoI and BglII

Positive: 2500 bp and 800 bp negative: 2070 bp



Gel map

Ladder	1C3	C9												
Ladder	A5	A6	A7	A8	A9	B5	B6	B7	B8	B9	C5	C6	C7	C8

Restriction validation

07/09/2014

code	name	colony	Status
A	PompC-amilCP	5	x
		6	x
		7	x
		8	x
		9	x
B	PT7-LuxI	5	V
		6	x
		7	x
		8	x
		9	x
C	PCI-LacI	5	x
		6	V
		7	x
		8	x
		9	V
cont	pSB1C3	----	ref

Cut with 2 restriction enzyme:

NcoI and BglII

Positive: 2500 bp and 800 bp negative: 2070 bp

PompC-amilCP - no positive

PT7-LuxI: colony 5 positive

PCI-LacI: colonies 6,9 Positive

Colony 9 will be taken

Gel map

Ladder	1C3	C9												
Ladder	A5	A6	A7	A8	A9	B5	B6	B7	B8	B9	C5	C6	C7	C8



Gibson assembly

14/09/2014

Calculation of the concentrations:

$$\left[\frac{pmol}{\mu l} \right] = \frac{\left[\frac{ng}{\mu l} \right] \times 10^3}{size \times 650}$$

cod e	name		size (bp)	ng/ ml	pmol/m l	Gibson assembly master mix(x2)		Plasmid (0.1pmol)	T7 poly (0.1pmol) 2739 bp	UPW	total	
X1	1C3: PCI-LacI (lin)		3491	137	0.06	15		1.66	0.78	2.56	20	
X2	1C3: Plac-Immbda CI (lin)		3176	122	0.06	15		1.69	0.78	2.53	20	
X3	1C3: Ptet-tetR (lin)		3028	175	0.09	15		1.12	0.78	3.09	20	
cod e	Reactio n name	Inser t size (bp)	plasmidsi ze (bp)	Inser t ng/ml	Plasmi d ng/ml	Insert pmol/ml	Plasmi d pmol/ml	Gibson assemb ly master mix(x2)	Plasmid (0.05pm ol)	Insert (0.15pm ol)	UP W	tot al
Y1	1C3: PCI-tetR	777	2163	141	128	0.28	0.09	15	0.55	0.54	3.91	20
Y3	1C3:Pla c-tetR	777	2226	141	167	0.28	0.12	15	0.43	0.54	4.03	20
Y5	1C3:Pte t-LacI	1245	2168	76	136	0.09	0.10	15	0.52	1.60	2.88	20
Y4	1C3:Pte t-Cl	867	2168	149	136	0.26	0.10	15	0.52	0.57	3.91	20
Y6	1C3:PT 7-amilCP	819	2137	63	162	0.12	0.12	15	0.43	1.27	3.30	20
Y2	1C3:Plu x-amilCP	819	2296	63	181	0.12	0.12	15	0.41	1.27	3.32	20

Incubation at 50 °C for 1 hr.
Store in -20 °C

Transformation to top 10

20 min frost on ice 100 mL
5 mL of 1AK3:Plux-ter (cir)
30 min on ice
1 min heat shock 42 °C
2 min on ice

incubation overnight.

Restriction validation



16/09/2014



Gibson assembly

23/09/2014

Calculation of the concentrations:

$$\left[\frac{pmol}{\mu l} \right] = \frac{\left[\frac{ng}{\mu l} \right] \times 10^3}{size \times 650}$$

code	Reaction name	Insert size (bp)	plasmid size (bp)	Insert ng/ml	Plasmid ng/ml	Insert pmol/ml	Plasmid pmol/ml	Gibson assembly master mix(x2)	Plasmid (0.05pmol)	Insert (0.15pmol)	UPW	total
Z1	Plux-g.luc-tetR	711	3341	15	82	0.03	0.04	15	1.32	4.62	0.00	20.95
Z2	Plux-g.luc-Cl	801	3341	31	82	0.06	0.04	15	1.32	2.52	1.16	20.00
Z3	Plux-g.luc-LacI	1179	3341	27	82	0.04	0.04	15	1.32	4.26	0.00	20.58
Z4	Plux-amilCP-Cl	801	3071	31	71	0.06	0.04	15	1.41	2.52	1.07	20.00
Z5	Plux-amilCP-lacI	1179	3071	27	71	0.04	0.04	15	1.41	4.26	0.00	20.66
Z6	Plux-amilCP-tetR	711	3071	15	71	0.03	0.04	15	1.41	4.62	0.00	21.03
M1	Plux-mcherry	727.00	2434	137.00	133.00	0.29	0.08	15.00	1.78	0.17	3.04	20.00

Incubation at 50 °C for 1 hr.
Store in -20 °C

Plux-mCherry

Reaction: Che A

PCR mix

component	Volume x1[ml]
phusion reaction buffer(x5)	10
dNTPs(10 mM)	1
che Af-forward primer	2.5
Che A r-reverse primer	2.5
Template (2ng/ml)	5
Phusion hot start II	0.5
DMSO	1.5
UPW	27
tot	50

PCR program

stage	Temp C ⁰	time
Initial denaturation	98	30 sec
35 cycles	98	10 sec
		sec
	72	1:40 min
Final extension	72	10 min
hold	4	--\--

Part name			
date			
concentration	[ng/ml]	[ng/ml]	[ng/ml]
-20 °C storage	iGEM B system \	iGEM B system \	iGEM B system \

Primers

name	seq	length	TM 1	Tm 2	Tn neb
Che G F	CCA CCA ATG GTG AGC AAG GGC G	22	64	64	1s: 55 2s: 72
Che G R	CAG GAA TCG GTT ACT TGT ACA GCT CGT CC	29	53	62	1s: 55 2s: 72
Che P F	GGA CGA GCT GTA CAA GTA ACC GAT TCC TGT TAA TTT GAA GG	41	54	64	1s: 66 2s: 72
CheP R	CGC CCT TGC TCA CCA TTG GTG GTT TCT CCT CTT TAA TGG	39	55	67	1s: 66 2s: 72


Transformation to top 10


20 min frost on ice 100 mL
 5 mL of product (cir)
 30 min on ice
 1 min heat shock 42 °C
 2 min on ice

incubation overnight.

Restriction validation

24/09/2014





Gibson assembly

06/10/2014

code	Reaction name	Insert size (bp)	plasmid size (bp)	Insert ng/ml	Plasmid ng/ml	Insert pmol/ml	Plasmid pmol/ml	Gibson assembly master mix(x2)	Plasmid (0.05pmol)	Insert (0.15pmol)	UPW	total
M2	Plux-mcherry-Cl	801	3116	31	43.00	0.06	0.02	15.00	7.07	0.84	0.00	22.91
M3	Plux-mcherry-LacI	1179	3116	27	43.00	0.04	0.02	15.00	7.07	1.42	0.00	23.48
M4	Plux-mcherry-tetR	711	3116	15	43.00	0.03	0.02	15.00	7.07	1.54	0.00	23.61
M6	Plux-mcherry-luxI	669	3116	70	43.00	0.16	0.02	15.00	7.07	0.31	0.00	22.38
M8	Plac-Cl-mcherry	798	3176	67	48.00	0.13	0.02	15.00	6.45	0.39	0.00	21.84

Colony PCR

Reaction mix: make total mix and then divide 20 [ml] to each PCR tube

component	Volume[ml] for 1 colony	For total mix 15 reaction	In?
Taq ready mix (x2)	10	150	V
DMSO	0.5	7.5	V
1C3 f (10 ng/ml)	1	15	V
1C3 r (10 ng/ml)	1	15	V
colony	colony	----	
UPW	7.5	112.5	V
tot	20		

