



## Why to do this :

1. Determine the best way to kill bacteria

## What you need :

1. Culture media : LB

- 10 g bactotrypton
- 5 g yeast extract
- 5 g NaCl
- 0,5 mL NaOH 10N
- Qsp 1 L

2. Antibiotics concentrations

Chloramphenicol (Cm) : 2 mg/mL

Tetracycline (Tet) : 1 mg/mL

Kanamycin (Kann) : 5 mg/mL

Ampicillin (Amp) : 10 mg/mL

→ 50 µL antibiotic / 5mL medium

3. Apparatus : UV exposure, PCR, microscope.

## How to do :

Be aware of being in sterile conditions!

1. Spots on LB plates to count survival bacteria

- a) Prepare 5 dilutions from the culture exposed to UV/temperature and the control:
  - 100 times : put 10µl culture into 990µl sterile water
  - 10<sup>4</sup> times : put 10µl previous culture into 990µl sterile water
  - 10<sup>5</sup> times : put 20µl previous culture into 180µl sterile water
  - 10<sup>6</sup> times : put 20µl previous culture into 180µl sterile water
- b) Drop off 10 µl of each concentration on LB plates (containing antibiotic appropriate with the strain) for the culture exposed to UV/temperature and the control.
- c) Incubate plates at 37°C, when drops start drying

2. Extraction of genomic DNA from concentrated cultures

- d) Preheat a bath at 99°C
- e) Centrifuge 800 µl liquid cultures
- f) Centrifuge 800 µl liquid cultures (in micro tubes) 5 minutes at 12000 rpm and discard half of the supernatant (to concentrate bacteria). Resuspend by pipetting up and down.
- g) Put it at 99°C for 5 min then vortex a little bit or pipette up and down. Repeat this operation.
- h) Centrifuge (in micro tubes) 5 minutes at 14000 g
- i) Take 20µl supernatant and add 4µl Buffer stop. Run an electrophoresis gel

### 3. PCR

a) Prepare in a micro tube a premix composed with

- Water	16,25 µl
- Primer Forward (FW)	1 µl
- Primer Reverse (RV)	1 µl
- Deoxy-nucleotide-triphosphate (dNTP)	0,5 µl
- Enzyme Q5 Buffer	5 µl
- Enzyme Q5	0,25 µl

b) The volume indicated must be multiply by  $n+1$  where  $n$  is the number of samples (to be sure to pipette the same volume for each).

c) In case of a culture being exposed to UV/temperature and its control,  $n$  would be equal to 3, leading to the following quantities:

- Water	48,75 µl
- Primer Forward (FW)	3 µl
- Primer Reverse (RV)	3 µl
- Deoxy-nucleotide-triphosphate (dNTP)	1,5 µl
- Enzyme Q5 Buffer	15 µl
- Enzyme Q5	0,75 µl

Be aware to put premix micro tube in ice after adding Q5 enzyme !

d) For each sample:

- Pipette 24 µl premix and put it in PCR tube. Keep them in ice.
- Add 1 µl extracted DNA.
- Centrifuge (by putting PCR tube in open microtube).
- Put it in PCR machine, already programmed (29 cycles + conservation at 4°C).

#### 4. Back Light coloration

- e) Centrifuge in micro tubes about 200 µl liquid culture (exposed to UV/temperature + control) 5 min at 15000g. Discard the supernatant.
- f) Resuspend in 300 µl NaCl 0,85 % (first washing)
- g) Centrifuge again 5 min at 15000g. Discard the supernatant
- h) Resuspend in 300 µl NaCl 0,85 %
- i) Add 0.2 µl premix AB from Back Light kit
- j) Centrifuge 5 min at 15000g. Discard the supernatant
- k) Resuspend in 300 µl NaCl 0,85 % (second washing)
- l) Let it rest in dark for about 15 minutes

#### 5. Preparing microscope slides (sterile conditions)

- a) Put slides in alcohol and hold them in fire with pliers by inclining slides down not to be burnt  
Resuspend in 300 µl NaCl 0,85 % (first washing)
- b) Pipette 25 µl coloured cultures on slides
- c) Put microscope slides at 60°C in Pasteur oven for a few minutes, until drops dry (to prevent bacteria from moving on slides)
- d) Observe slides with an epifluorescence microscope