

Preparing your samples:

- Use 0.5ml Eppies
- Add 250ng of DNA
- Add primer to a concentration of 5pmol
- Add dH₂O to a final volume of 8 microliters
- Label tubes (Same label as online order)

Ordering the sequencing:

- Go to Etonbio.com
- Log in and enter payment information
- Select plasmid sequencing
- Fill out forms (sample names matching tube labels)
- Print order form
 - (Samples should be in the mailbox by 3pm)

Shipping the samples:

- Place eppie tubes in a 50ml tube
- In the mailroom, place the 50ml tube in an envelope with the printed order form
- Label the envelope DNA sequencing iGEM Duffylab
- Put in box labeled Eton Biotech