

June 23, 2014

**Estrogen Receptor Intein Construction: PCR**

Purpose: To synthesize the estrogen receptor intein

**Table 1. PCR Protocol to Synthesize Estrogen Receptor**

Reagent	Volume (μl)	Volume (μl)
Water	12.4	5.4
5X HF PCR Buffer	4	4
dNTPs	0.4	0.4
Primer 1 (10 uM Nterm DNA)	1	5
Primer 2 (10 uM Cterm DNA)	1	5
Phusion Taq enzyme	0.4	0.4
	Total: 20.4	~20.4

**Table 2. PCR Conditions for First Round**

Temperature ( °C)	Time
98	2 minutes
98	10 seconds
60	15 seconds
72	1.5 minutes
72	1 minute
4	Hold

- Repeat for 35 cycles
- Take 2.5 μl DNA from this PCR reaction for the next PCR

**Table 3. PCR Protocol to Synthesize Estrogen Receptor**

Reagent	Volume (μl)
Water	31
5X HF PCR Buffer	10
dNTPs	1

Primer 1 (10 uM Nterm DNA)	2.5
Primer 2 (10 uM Cterm DNA)	2.5
DNA	2.5
Phusion Taq enzyme	0.5
	Total: 50

- Same Conditions
- Take 2 - 3  $\mu$ l DNA and run in 1 % agarose gel

**Results:** No bands present

**Overnight:**

- Cultures for Top10 Cells expressing fluorescent proteins

June 24, 2014

**PCR to Synthesize Estrogen Responsive Intein**

- Repeat protocol from previous day (tables 1, 2, and 3 from June 23, 2014)

**PCR to Prepare Fluorescent Protein DNA for Submission to iGEM Registry**

**Table 1. Protocol for PCR Cocktail**

Reagent	Volume (μl)
Water	155
5X HF PCR Buffer	50
dNTPs	5
Primer 1 (Prefix)	12.5
Primer 2 (Suffix)	12.5
Phusion Taq enzyme	2.5
	Total: 237.5

- Add 47.5 μl of PCR Cocktail to each tube (5 tubes)
- Add 2.5 μl of DNA to corresponding tube
- PCR Conditions in Table 2. (June 23, 2014)

**Overnight:**

- Streak out MACH cells with fluorescent proteins
- Culture of MACH cells containing WT Lac

June 25, 2014

### PCR for RFP with His tag

**Table 1. Protocol for Single PCR Reaction**

Reagent	Volume (µl)
Water	31
5X HF PCR Buffer	10
dNTPs	1
Primer 1 (10 uM Prefix)	2.5
Primer 2 (10 uM Suffix)	2.5
DNA	2.5
Phusion Taq enzyme	0.5
	Total: 50

- PCR Conditions in Table 2. (June 23, 2014)

### PCR for Killer Red and Super Nova with BioBricks Prefix and Suffix

**Table 2. Protocol for Single PCR Reaction**

Reagent	Volume (µl)
Water	31
5X HF PCR Buffer	10
dNTPs	1
Primer 1 (10 uM Prefix)	2.5
Primer 2 (10 uM Suffix)	2.5
DNA (Killer Red or Super Nova respectively)	2.5
Phusion Taq enzyme	0.5
	Total: 50

- PCR Conditions in Table 2. (June 23, 2014)

### PCR for Killer Red and Super Nova with Ribosome Binding Site

**Table 2. Protocol for Single PCR Reaction**

Reagent	Volume (µl)
Water	31
5X HF PCR Buffer	10
dNTPs	1
Primer 1 (10 uM XbaPreRBS)	2.5
Primer 2 (10 uM PstSuf6His)	2.5
DNA (Killer Red or Super Nova respectively)	2.5
Phusion Taq enzyme	0.5
	Total: 50

- PCR Conditions in Table 2. (June 23, 2014)

**Miniprep of WT Lac Plasmid**

Protocol for GeneJet Plasmid Miniprep:

<http://www.thermoscientificbio.com/uploadedfiles/resources/k0502-product-information.pdf>

Yield = 225.5 ng/µl

**Digestion of Wild Type Lac Vector and Fluorescent Proteins for Ligation****Table 3. Protocol for Digestion of WT Lac Plasmid**

Reagent	Amount (µl)
10X FastDigest Buffer	2
Plasmid DNA	8
Restriction Enzyme PstI	2
Restriction Enzyme SpeI	2
Sterile Water	6
<b>Total Volume</b>	20

- Digest at 37 °C for 30 - 45 minutes
- Gel Photo (Name: 062514b)

**Table 4. Protocol for Digestion of Fluorescent Proteins, Killer Red and Super Nova Plasmids**

Reagent	Amount (µl)
10X FastDigest Buffer	2
Plasmid DNA	14
Restriction Enzyme XbaI	2
Restriction Enzyme SpeI	2
<b>Total Volume</b>	20

- Digest at 37 °C for 30 - 45 minutes
- Gel Photo (Name: 062514c)
  - All fluorescent proteins and Killer Red (prefix/suffix) showed bands
  - None for Killer Red (RBS) or Super Nova
- Protocol: <http://www.thermoscientificbio.com/uploadedFiles/Resources/fast-digestion-dna.pdf>

### **Ligation of Fluorescent Proteins into WT Lac and BioBrick Vectors**

**Table 5. Ligation Reaction of Fluorescent Proteins**

Reagent	Amount (µl)
Vector DNA (WT Lac or BioBrick Vector)	2
Insert DNA (Fluorescent Proteins)	6
Ligation Buffer	1
T4 Ligase Enzyme	1
<b>Total Volume</b>	10

- Ligate for 10 minutes at room temperature

### **Transformation of pSB3K3 and WT Lac Plasmid into MACH Cells**

**Purpose:** To be used in cloning

- 50 µl cells for each ligation reaction
- Plate 400 µl transformants on LB + CAM plates
- Incubate at 37 °C overnight

iGEM Transformation Protocol: <http://parts.igem.org/Help:Protocols/Transformation>

June 26, 2014

**Overnight:**

- Cultures of the transformants (FP's in BBKR and WT Lac Backbones)