

Name: Jennifer Zhang

Date: 6/13/14

WEEKLY ONE SENTENCE: Test our RBS primers and promoters with machinery

CONFIRMATION?

gel ✓

GEL LANE ORDER:

DAILY ONE SENTENCE: ~~test~~ insert promoter M1

in front of gel purified ~~the~~ RBS M2 and L1 (w/ machinery) (and purify RBS primers w/ machinery (w/ 50°C PCR protocol))

METHOD(S) OF CHOICE:

PCR
gel extraction

GEL PICTURE / OTHER SPACE:

• ran a gradient PCR w/ R2.1 with ~~temp~~ 54, 50, 58 to see if different temps will optimize

gel: 1kb 54 56 58 60 62 64

EXPECTED RESULT:

promoter M1 → band of ~1kb
RBS primers → band of ~1kb



(~~check~~ if APE file exists on Drive)

TEMPLATE(S): RBS M2/L1 with machinery (gel purified)

OLIGO(S): promoter M1 → O:... new
O2U3

THERMOCYCLER SETTINGS: 94 - Z; {94-15", 58-15", 58-1.5", 58-4", 4 forever

SOME BS 6/12? w/ 1m

POSITIVE:

NEGATIVE: no polymerase ✓

Name:

Cohen Tam

Date:

6/13/14

WEEKLY ONE SENTENCER:

CONFIRMATION?

GEL LANE ORDER:

gel

DAILY ONE SENTENCER:

troubleshooting the
RBS primer PCR & inserting promoter primers

METHOD(S) OF CHOICE:

PCR, gel extract, gel electrophoresis

GEL PICTURE / OTHER SPACE:

EXPECTED RESULT:

band at ~900-1000 for RBS primer
& promoter primer

☒ check if APE file exists on Drive)

TEMPLATE(S):

OLIGO(S):

963, RBS primers (74-79), promoter primers (80-85)

THERMOCYCLER SETTINGS:

94-2', { 94'-15" } 68-1.5' } x3, 68-4', 4

CONTROLS:

POSITIVE:

NEGATIVE:

No polymerase

~~no template~~

Name: Katie Fisher

CONFIRMATION?

Date: 6/13/14

thin gel ok.

WEEKLY ONE SENTENCER:

GEL LANE ORDER:

clone diff. combinations of RBS primers + promoter
into RBS plasmid

DAILY ONE SENTENCER:

~~Run~~ Gel extraction for PCR run @ 65° for annealing
Run PCR w/ gel extract products + promoters + emp. insert promoters + then
gel extraction

GEL PICTURE / OTHER SPACE:

METHOD(S) OF CHOICE:

gel extraction

gel purification

PCR

EXPECTED RESULT:

band btw 900 + 1000 base pairs for first PCR product
w/o promoter

1100 base pairs for other 2 PCR products

(☒ check if APE file exists on Drive)

TEMPLATE(S): mcherry + mcherry M2 + L1 from yesterday

OLIGO(S): promoter + 303

THERMOCYCLER SETTINGS:

CONTROLS:

use day before #'s

POSITIVE:

NEGATIVE:

PCR ~~product~~ product run w/ no
polymerase ~~at all~~
no DNA template