

General Labjournal

25.06.2014

Grow on shakeflasks the 25 samples that we requested from the iGEM HQ. The biobrick numbers of the requested samples: K1172305, K1172306, K1172401, K1172403, K1172404, K917013, K917014, K1127006, K1172303, K1172304, K1179019, K917012, K917009, K917006, K917003, K1172405, K1172501, K1172502, K1172503, K1172504, K1172505, K1172507, K540000, K342003, K540001.

26.06.2014

Miniprep + Glycerol stock of the 25 samples that we requested from the iGEM HQ.

07.07.2014

Prepare solutions for competent cells preparation (MgCl₂ and CaCl₂)

08.07.2014

Prepare all the sequencing mixtures. The parts to be sequenced were: K1172306, K1172403, K1172404, K1172303, K917009, K917006, K917003, K1172502, K540000 and K342003.

15.07.2014

- Made competent cells of C43(DE3) and BL21(DE3) for transformation
- Cultivate the transformed constructs from 14.07.2014 in LB medium for miniprep
 - Also cultivate mKate (amp), eGFP (amp) and Rhamnose promoter - Bba_K914003 (cam)
 - And the strain AYCE189 from Lu lab is also cultivated (for a glycerolstock)

16.07.2014

- Miniprepped the samples cultivated on 15.07.2014, except for the strain AYCE189.
- Measured the concentration of the isolated DNA with nanodrop:
 - PAYC002 122.8 ng/ul
 - rr12y(rii)g 52.6 ng/ul
 - PAYC003 187.0 ng/ul
 - rrjt12(11)g 32.3 ng/ul
 - PAYC005 27.8 ng/ul
 - PAYC008 22.3 ng/ul
 - PAYC006 37.0 ng/ul
 - PAYC007 22.3 ng/ul
 - l5023 22.5 ng/ul
 - C640 14.3 ng/ul
 - mKate 39.2 ng/ul
 - eGFP 165.1 ng/ul
 - Rhamnose 47.7 ng/ul
- Cultivated the samples again in shakeflasks.
- Made glycerolstocks of all the cultivated samples from 15.07.2014.

Test the competency of C43(DE3)

- 30 ul competent C43 cells + 100 ng pUC19 (Amp^R)(190 ng/ul), plated in duplo on Cam (neg. control), Amp (pos. control) LB plates and on LB plate without antibiotics

- 30 ul competent C43 cells plated in duplo on Cam (neg. control), Amp (neg. control) LB plates and on LB plate without antibiotics (pos. control)

17.07.2014

Results of the plates for positive and negative control of C43

- | | |
|---|-----------------|
| • C43 + pUC19(AmpR) + Amp plates | Much growth |
| • C43 + pUC19(AmpR) + Cam plates | No growth |
| • C43 + pUC19(AmpR) + without antibiotic | Much growth |
| • 2x diluted C43 + pUC19(AmpR) + Amp plates | Single colonies |
| • 2x diluted C43 + pUC19(AmpR) + Cam plates | No growth |
| • 2x diluted C43 + pUC19(AmpR) + without antibiotic | Much growth |
| • C43 + no plasmid + Amp plates | No growth |
| • C43 + no plasmid + Kan plates | No growth |
| • C43 + no plasmid + Cam plates | No growth |

We decided to not dilute the competent cells for the future experiments.

CFU for pUC19(AmpR): 1120 colonies on the plate.

Test the competency of BL21(DE3)

- 30 ul competent BL21 cells + 100 ng pUC19 (AmpR)(190 ng/ul), plated in duplo on Cam (neg. control), Amp (pos. control) LB plates and on LB plate without antibiotics
- 30 ul competent BL21 cells plated in duplo on Cam (neg. control), Amp (neg. control) LB plates and on LB plate without antibiotics (pos. control)
- 2x diluted 30 ul competent BL21 cells + 100 ng pUC19 (AmpR)(190 ng/ul) plated in duplo on Amp (pos. control) LB plates

18.07.2014

Results of the plates for positive and negative control of BL21(DE3)

- | | |
|--|-----------------|
| • BL21 + pUC19(AmpR) + Amp plates | Much growth |
| • BL21 + pUC19(AmpR) + Cam plates | No growth |
| • BL21 + pUC19(AmpR) + without antibiotic | Much growth |
| • 2x diluted BL21 + pUC19(AmpR) + Amp plates | Single colonies |
| • BL21 + no plasmid + Amp plates | No growth |
| • BL21 + no plasmid + without antibiotic | Much growth |
| • BL21 + no plasmid + Cam plates | No growth |

22.07.2014

Prepared antibiotics, 500ul in each eppendorf cup. Saved in 'Antibiotics' box in the freezer with other iGEM materials. 1000x stock solutions.

- Chloramphenicol diluted in EtOH
34 mg/ml in EtOH -> 0,391 gram diluted in 11,5 ml EtOH
- Kanamycin diluted in H₂O
10 mg/ml in H₂O -> 0,095 gram diluted in 9,5 ml H₂O
- Ampicillin diluted in H₂O
100 mg/ml in H₂O -> 0,92 gram diluted in 9,2 ml H₂O

Pre-culture DH5-alpha has been prepared, 100 mL LB medium is cultivated and placed in the shaker. Overnight at 37 degrees and 180 rpm.

23.07.2014

29-7-2014								
Making the Trace Elements Solution								
Chemical in paper - 1	Chemical in our buffer -2	g/mol - 1	g/mol -2	Weight (gram) in article	Weight (gram) in our buffer	Company	CAS number	Cupboard number
Na2EDTA	Na2EDTA.2H2O	336,2	372,24	2,26	2,517	Sigma-Aldrich ¹	6381-92-6	E25
MgSO4.7H2O		246,47		24,89	24,89	JT Baker ²	10034-99-8	M03
MnSO4.4H2O	MnSO4.H2O	169,02	223,08	0,029	0,022	Sigma-Aldrich ¹	10034-96-5	M07
NaCl		58,44		0,058	0,058	JT Baker ²	7647-14-5	S17
FeCl2	FeCl2.4H2O	126,73	198,81	0,068	0,107	Tu Delft afd. FD-LMS	231-729-4	I04
CoCl2	CoCl2.6H2O	129,81	237,93	0,065	0,119	JT Baker ²		C08
ZnSO4.7H2O		287,54		0,029	0,029	Merck		Z03
CuSO4.5H2O		249,68		0,005	0,005	JT Baker ²	7758-99-8	C12
H3BO3		61,83		0,35	0,35	USB	10043-35-3	B11
Na2MoO4	Na2MoO4.2H2O	205,91	241,95	0,08	0,094	Sigma-Aldrich ¹	10102-40-6	S24
NiCl2.6H2O		237,7		0,119	0,119	Aldrich	7791-20-0	N41
Na2SeO4		188,9		0,028	0,028	Sigma		S55
- FeCl2.4H2O is diluted in 22,5 mL HCL								
- Adding order: EDTA + 800 mL H2O -> FeCl2.4H2O + 25 mL 6 M HCl -> MgSO4.7H2O -> etc								
- Precipitation after adding MnSO4.4H2O -> o/n on a stirrer -> precipitation still there (30-7-2014) -> failed..								

30.07.2014

Esra

Making the M4 minimal medium Trace Elements Solution (third attempt) + buffer (exact contents will be updated)

-Changed adding order

-Adding order 1- EDTA, 2- Mg, 3-Mn, 4-NaCl, 6-CoCl etc.

-Added number 5- (FeCl diluted in 22.5 mL HCl) at the end!!

-All added except FeCl (in HCl) -> no precipitations!! All is diluted well!

Inoculated the transformed iGEM registry constructs

Resistance	Transformation	Grown in (ml LB):
Cam	DH5a + k823017	5ml & 10 ml
Cam	DH5a + k808000	5ml & 10 ml
Kan	DH5a + l20260	5ml & 10 ml
Cam	DH5a + 80017	5ml & 10 ml
Amp	DH5a + J231100	5ml & 10 ml
Amp	DH5a + pUC19	5ml & 10 ml

20.08.2014

Janna

Made medium M4 with different carbon sources.

1. 400 ml 40 mM D/L-lactate M4

342 ml MilliQ

10 ml Buffer 40 x

4 ml 0.1 M CaCl_2

40 ml 0.4 M D/L-lactate

4 ml Trace elements 100 x

2. 400 ml 40 mM glycerol M4

342 ml MilliQ

10 ml Buffer 40 x

4 ml 0.1 M CaCl_2

40 ml 0.4 M Glycerol

4 ml Trace elements 100 x

3. 400 ml 40 mM glucose M4

342 ml MilliQ

10 ml Buffer 40 x

4 ml 0.1 M CaCl_2

40 ml 0.4 M Glucose

4 ml Trace elements 100 x

21.08.2014

Janna

Tested competent cells BL21 culture 1, BL21 culture 2 and C43. I used pUC19 as test DNA (ampR). Made the following combinations:

1. 30 µl BL21.1 with 1 µl MilliQ
2. 30 µl BL21.1 with 1 µl pUC19
3. 30 µl BL21.2 with 1 µl MilliQ
4. 30 µl BL21.2 with 1 µl pUC19
5. 30 µl C43 with 1 µl MilliQ
6. 30 µl C43 with 1 µl pUC19

There were no colonies on the plates, so the cells are not competent.

Cristy: Made competent cells of CsgB (approximately 30 eppendorfs containing 100ul competent cells). OD600: 0,567 and 0,589

21.08.2014

Joan

Prepare samples for Melbourne iGEM team:

- pET23b - Ulp1-His6 (AmpR)
- BBa_K1022107:pcI-Ulp in pSB1C3 col2
- BBa_K1022113:pBAD-Ulp-TT in pSB1C3 col2

26.08.2014

Chloramphenicol diluted in EtOH

34 mg/ml in EtOH -> 0,3912 gram diluted in 11,5ml EtOH. Divided in aliquots of 500 ul.

Refilled stock with 22 new cups.

Janna

Making M4 with D/L-lactate. Same protocol as 20/8, only filter sterilized the whole bottle after making it.

28.08.2014

Janna

Transformation of pUC19 in C43+ET20 as test

Made the following transformations:

0.5 ul pUC19 plasmid (concentration of 100.2 ng/ul) added to 30 ul of C43+ET20 competent cells

1 ul MilliQ added to 30 ul of C43+ET20 competent cells as a negative control

Followed the transformation in home-made competent cells protocol and made six plates with each 100 ul:

pUC19 in C43+ET20 (ampR and camR):

Ampicillin - some growth

Chloramphenicol - some growth

Amp and Cam - lots of growth

MQ in C43+ET20 (camR):

Ampicillin - no growth

Chloramphenicol - some growth

Amp and Cam - no growth

This means the cells are competent, but the efficiency is low.