

E. Coli - CUPID 


Shely Agustina 21113038
Kartika Fandika 21113043
Ramadhani Safitri 21112014
Risha A Pratiwi 20613027

WHY

E. Coli - CUPID ? 

Cuprum Indicator

FOOD DEMANDS



INTENSIVE FARMING




CHEMICALS PESTICIDES USAGE



Chemicals Pesticides :

A chemicals compound used for destroying insects or other organisms harmful to cultivated plants or to animals



- **Advantages:** Quick result and Easy to apply
- **Disadvantages:** High intensity usage is toxic for human and environment

Copper based pesticides is one of the most used pesticides due to broad spectrum function

Copper-based pesticide	Use pattern(s)
Copper (metallic)	Algaecide, antifouling paint
Copper (metallic in the form of chelates of copper citrate and copper gluconate)	Algaecide, bactericide, fungicide
Copper carbonate	Algaecide, herbicide, wood preservative
Copper ethanolamine complex	Algaecide, wood preservative
Copper ethylenediamine complex	Herbicide
Copper hydroxide	Antifouling paint, bactericide, fungicide, plant growth regulator, wood preservative
Copper naphthenate	Wood preservative
Copper oxychloride	Algaecide, bactericide, fungicide
Copper salts of fatty and rosin acids	Bactericide, fungicide
Copper sulfates	Algaecide, bactericide, desiccant, fungicide, herbicide
Copper triethanolamine complex	Algaecide
Copper oxides	Algaecide, antifouling paint, wood preservative

- **Copper Hydroxide**
Agrocide 77WP,
Champion77WP, Kocide54WG
- **Copper Oxide**
Nordox 56WP & 86WP
- **Copper Oxine**
Hylite 150EC
- **Copper Hexa Fluorocylcate**
Cefka 97SP
- **Copper Oxichloride**
Kasumin 5/75WP, Probox 50WP
- **Copper Oxisulphate**
Sultricob 93 WP
- **Copper Sulphate**
Betawood 97WG

COPPER EXCESSIVE SYMPTOMS in PLANTS

Chlorosis



Silene vulgaris shows that soil contains high concentrations of cuprum



Necrosis



COPPER TOXICITIES in HUMAN

Main symptoms of Copper poisoning

Systemic

- Chills
- Fever
- Pain

Eyes and skin

- Yellowing (jaundice)

Circulatory

- Anemia
- Shock

Mouth

- Metallic taste

Muscular

- Convulsion
- Ashes
- Weakness

Gastric

- Vomiting
- Nausea
- Abdominal pain
- Burning sensation

Kidneys

- No urine production

Intestinal

- Diarrhea (may be bloody or bluish)




Table 4.1 Critical tissue concentrations for Cu toxicity in a variety of species

Species	Plant part	Critical toxicity (mg.kg ⁻¹)	Method of cultivation	Reference
Twelve spp. of tropical grass	Shoots	18 ¹	Soil pot experiment	(Plenderleith and Bell, 1990)
<i>A. hypogaea</i>	Shoots	230 ²	Soil pot experiment	(Borkert <i>et al.</i> , 1998)
<i>Carthamus tinctorius</i> (safflower)	Leaves	10 ¹	Sand culture	(Pandey and Sharma, 1999)
<i>G. max</i>	Shoots	140 ²	Soil pot experiment	(Borkert <i>et al.</i> , 1998)
<i>O. sativa</i>	Whole plant	35 ³	Solution culture	(Lidon and Henriques, 1992)
<i>O. sativa</i>	Shoots	<20 ²	Soil pot experiment	(Borkert <i>et al.</i> , 1998)
<i>T. aestivum</i>	Shoots	75 ⁴	Solution culture	(Wheeler and Power, 1995)
	Roots	300		
<i>Vigna mungo</i> (blackgram)	Leaves	67 ¹	Soil pot experiment	(Kalyanaraman and Sivagurunathan, 1993)
	Stems	50		
	Roots	41		

¹Tissue concentration at 90% of maximum yield

²Intersection of regression lines: "shoot growth independent of leaf metal concentration" and "growth inhibition at higher metal levels"

³Nutrient calibration curve of for total tissue concentration and root growth plotted on a log-log graph

⁴Estimated as a narrow range of plant concentrations where relative yield decreases rapidly

Bio Brick Design:

E. Coli - CUPID




OBJECTIVES:

COPPER INDICATOR IN WATER

USING ENGINEERED *Escherichia coli*

PARTS USED:

Design *E. coli* CUPID using assembly method RCF 10



Forward terminators

-?-		Name	Description	Direction	Efficiency Fwd. Rev.		Chassis	Length
1★		BBa_B0010	T1 from <i>E. coli</i> <i>rrnB</i>	Forward				80
A	^	BBa_B0012	TE from coliphage T7	Forward	0.309[CC]	-0.308[CC]		41
A	X	BBa_B0013	TE from coliphage T7 (+/-)	Forward	0.6[CC]	-1.06[CC]		47
1★	W	BBa_B0015	<u>double</u> terminator (B0010-B0012)	Forward	0.984[CC] 0.97[JK]	0.295[CC] 0.62[JK]	<i>E. coli</i>	129
1★		BBa_B0017	double terminator (B0010-B0010)	Forward				160

PARTS USED:

Protein coding sequences/Selection markers

[< Back to Protein coding sequences](#)

Selection [markers](#) are [protein coding sequences](#) that confer a selective advantage or disadvantage to host chassis. For example, a common type of prokaryotic selection marker is one that confers resistance to a particular antibiotic. Thus, cells that carry the selection marker can grow in media despite the presence of antibiotic. Most plasmids contain antibiotic selection markers so that researchers can ensure that the plasmid is maintained during [cell replication](#) and division. (Cells that lose a copy of the plasmid will soon either die or fail to grow in media supplemented with antibiotic.) A second common type of selection marker, often termed a positive selection marker, are those that are toxic to the cell. Positive selection markers are frequently used during cloning to select against cells transformed with the cloning vector and ensure that only cells transformed with a plasmid containing the insert.

[More...](#)

-?-	Name	Protein	Description	Tag	Direction	UniProt	KEGG	Length
1★	BBa_T9150	PyrF	orotidine 5	None	Forward	P08244	eco:b1281;	741
1★ W	BBa_J31002	AadA-bkw	kanamycin resistance backwards (KanB) [cf. BBa_J23012 & BBa_J31003]			P0AG05	none	816
A	BBa_J31003	AadA2	kanamycin resistance forward (KanF) [cf. BBa_J23012 & BBa_J31002]			P0AG05	none	816
A	BBa_J31004	CAT-bkw	chloramphenicol acetyltransferase (backwards, CmB) [cf. BBa_J31005]			P62577	none	660
1★ W	BBa_J31006	TetA(C)-bkw	tetracycline resistance protein TetA(C) (backwards) [cf. BBa_J31007]			P02981		1191
1★ W	BBa_J31005	CAT	chloramphenicol acetyltransferase (forwards, CmF) [cf. BBa_J31004]			P62577	none	660
1★ W	BBa_J31007	TetA(C)	tetracycline resistance protein TetA(C) (forward), [cf. BBa_J31006]			P02981		1191
1★	BBa_K145151		ccdB coding region					306
A W	BBa_K143031		Aad9 Spectinomycin Resistance Gene					771
1★	BBa_K156011		aadA (streptomycin 3'-adenyltransferase)					789
A W	BBa_K389005		Kanamycin resistance					819

PARTS USED:

Promoters/Catalog/Metal sensitive

This set includes [promoters](#) that are sensitive to various metals. The promoters are typically regulated by a [receptor protein](#) that binds to the metal ion or complex. There are promoters that are sensitive to iron, lead, and copper.

-?-	Name	Description	Promoter Sequence	Positive Regulators	Negative Regulators
1★	BbA_I721001	Lead Promoter	... gaaaacctgtcaatgaagagcgatctatg		
	BbA_I731001	FeaA promoter	... ttatglttgagactatagctgaacacacac		
	BbA_I760005	Cu-sensitive promoter	atgacaaaattgtcat		
	BbA_I765000	Fe promoter	... accaattcctgggaacggtcagggcactaa		
	BbA_I765007	Fe and UV promoters	... ctgaaagcgcataccgctatggaggggggt		
1★	BbA_J3902	PrFe (PI + PII rus operon)	... tagatatgcctgaaagcgcataccgctatg		

Anderson RBS family members

Identifier	Sequence ^a	Predicted Strength ^b	
		Mean	CV
Master Sequence	TCTAGAGAAAGANNNGANNNACTAGATG		
BBa_J61100	TCTAGAGAAAGAGGGGGACAAACTAGATG		
BBa_J61101	TCTAGAGAAAGACAGGACCCACTAGATG		

Part:BBa_K592010

Designed by: Lei Sun Group: iGEM11_Uppsala-Sweden (2011-09-18)



amilGFP

Not Released

Sample Pending

Experience: Works

8 Uses

Get This Part

amilGFP, yellow chromoprotein



This chromoprotein from the coral *Acropora millepora*, amilGFP, [naturally exhibits](#) strong [yellow color](#) when expressed. The color is readily visible to naked eye both in LB-culture and on agar plates. Color development can be seen in less than 24 hours of incubation.

Usage and Biology



This part is useful as a [reporter](#).




iGEM11_Uppsala-Sweden: Expression of chromoproteins. The images above show *E. coli* constitutively expressing amilGFP BBa_K592010 (yellow), amilCP BBa_K592009 (blue), and RFP BBa_E1010 (red).




Desain Bio Brick CUPID




Cut with EcoRI and SpeI




Cut with EcoRI and XbaI



Bio Brick Design:



Overview CUPID System



FUTURE DEVELOPMENT

FUTURE RESEARCH:

- Quantitative Measurement of copper content
- To enhance *E. coli* resistance to high excessive copper



FUTURE PRODUCT

EASY AND CHEAP KIT

Blotter paper color indicator

REFERENCES

Fishel, F.M. 2011. Pesticide Toxicity Profile: Copper-based Pesticides. University of Florida IFAS extension.

Franke. S., et al, 2003. Molecular Analysis of the Copper-Transporting Efflux System CusCFBA of *Escherichia coli*. Journal of Bacteriology: 3804- 3812

Rademacher. R., et al. 2012. Review: Copper-responsive gene regulation in bacteria. Microbiology 158:2451-2464

Reichman. S. M., 2002. The responses of plants to Metal Toxicity: A review focusing on Copper, manganese and Zinc. Australian Minerals & Energy Environment Foundation

World Health Organization. 2004. Copper in drinking water. WHO publication

Stern. B. R.,et al. 2007. Copper and Human Health: Biochemistry, Genetics, and Strategies for Modelling Dose response Relationships. Journal of Toxicology B.10:157-222

[Http://parts.igem.org/main_page](http://parts.igem.org/main_page)



T

H A

N K

Y

O

U

2005

MALL BUGGY

I'M HERE TO CONQUER

FORE 2005

We're on Track For a Cure

**TAMBAHAN KALAU ADA
PERTANYAAN YANG BERHUBUNGAN**





Table 1

Essential elements of copper homeostasis in *E. coli*

Homeostatic mechanism	Regulated by	Function
CopA	CueR (sensing cytoplasmic Cu[I]) and CpxR (sensing cell envelope stress)	Detoxification of cytoplasmic Cu[I]
CusCFBA	CusRS (two-component regulation system, sensing periplasmic Cu[I])	Detoxification of periplasmic (and possibly cytoplasmic) Cu[I]
CueO	CueR (sensing cytoplasmic Cu[I])	Protection of periplasmic proteins from copper-induced toxicity
PcoABCD	PcoRS, (two-component regulation system, sensing periplasmic Cu[I]) (and CusRS)	Protection from extreme periplasmic copper stress
PcoE	CusRS (PcoRS)	Periplasmic copper chaperone, copper binding



COLOUR
FIGURE

Fig. 1. Copper homeostasis in *E. coli*. Cu^+ and Cu^{2+} ions enter the periplasm (shown in grey), probably via porins spanning the outer membrane. ComC lowers the permeability of the outer membrane for copper ions. It is unknown whether ComC interacts with porins. At low ambient copper concentrations, *comC* transcription is repressed by ComR (not shown). Cu^+ ions proceed from the periplasm into the cytoplasm by an unknown mechanism, while Cu^{2+} ions do not cross the cytoplasmic membrane. Copper efflux involves the ATPase CopA and the multicomponent system CusCFBA. The multicopper oxidase CueO oxidizes periplasmic Cu^+ to Cu^{2+} . For further details, see text.




Fig. 2. Functional model of the Cus efflux complex. The four-part Cus complex consists of the inner membrane pump CusA, the periplasmic protein CusB and the outer membrane protein CusC forming a channel bridging the periplasmic space. Entry of copper may occur from the periplasm (A), from the cytoplasm (B) or via the copper-binding chaperone CusF from the periplasm (C).




FIGURE 1. Schematic drawing of the key elements of copper homeostasis and defects in Wilson and Menkes diseases. Dietary copper in the intestine may be taken up by intestinal cells via DMT1 and also by an as yet unidentified system. Metallothioneins (MT) can buffer excess copper inside these and other cells. The Menkes copper ATPase (MNK) exports copper from intestinal cells into the circulation. In the blood, copper can form complexes with small molecules like histidine and is bound to proteins such as albumin or transcurein. The liver takes up circulating copper by a first-pass process, probably involving CTR1. In hepatocytes (liver cells), copper is transported into the trans-Golgi network by the Wilson ATPase (WND) for incorporation into ceruloplasmin (CP). Ceruloplasmin, a multicopper ferroxidase with a role in iron homeostasis, is secreted into the circulation via the secretory pathway. Excess copper is excreted into the bile by hepatocytes. This process also requires the activity of the Wilson ATPase, which undergoes copper-induced trafficking from a trans-Golgi to a periplasmic location. The transport of copper across the blood–brain barrier (BBB) appears to be catalyzed by the Menkes ATPase expressed in cerebrovascular endothelial cells. The lack of function of the Menkes copper ATPase in Menkes disease (—MD) leads to copper accumulation in the intestine and concomitant copper deficiency in most tissues. Systemically administered copper cannot be transported to the brain. A defect in the Wilson copper ATPase in Wilson disease (—WD) results in copper accumulation in hepatocytes.

TABLE 1. A Postulated Spectrum of Copper Metabolism (from Aggett, 1999)

Dose range	Approximate daily intakes	Health outcomes
Toxic	>5.0 mg/kg body weight	Death Gross dysfunction and disturbance of metabolism of other nutrients; hepatic “detoxification” and homeostasis overwhelmed
	100 µg/kg body weight	Gastrointestinal metallothionein induced (possible differing effects of acute and chronic exposure);
Adequate	34 µg/kg body weight	Plateau of absorption maintained; homeostatic mechanisms regulate absorption of copper
	11 µg/kg body weight	Hepatic uptake, sequestration and excretion effect homeostasis; glutathione-dependent uptake of copper; binding to metallothionein; and lysosomal excretion of copper
	9 µg/kg body weight	Biliary excretion and gastrointestinal uptake normal
Deficient	8.5 µg/kg body weight	Hepatic deposit(s) reduced; conservation of endogenous copper; gastrointestinal absorption increased
	5.2 µg/kg body weight	Negative copper balance
	2 µg/kg body weight	Functional defects, such as lysyl oxidase and superoxide dismutase activities reduced; impaired substrate metabolism; Peripheral pools disrupted; gross dysfunction and disturbance of metabolism of other nutrients; death

METHODS TO DETECT PESTICIDES

Chromatographic Methods → Gas Chromatography (GC) & High Performance Liquid Chromatography (HPLC) , coupled with mass spectrometry (MS)

- Sensitive and realiable
- Complex
- A lot f time to realize the analysis
- Required highly trained technicians
- not in-field / on-site detection


METHODS TO DETECT COPPER IN WATER

The *Copper reagents and applications* table lists proprietary reagents and applications.

Table 1 Copper reagents and applications

Reagent	Form measured		Application
	Without pretreatment	With digestion	
CuVer 1™ ¹	Free	Total recoverable	water, wastewater
CuVer 2™	Total dissolved copper	Total recoverable	
Free copper reagent	Free	Total recoverable	hard water, wastewater, seawater
Porphyrin reagents	Free	Total recoverable	extremely low levels in water, wastewater and seawater

¹ CuVer is a trademark of Hach Company.




**Organophosphate
Pesticides (OP)**

**Carbamate
Pesticides (CP)**

**Types of
Chemical
Pesticides**

**Organochlorine
Insecticides (OI)**

**Pyrethroid
Pesticides (PP)**



COPPER FOR CROP PRODUCTION

- **Classified as a micronutrient**

Role in plant growth :

- promotes seed production and formation
- plays an essential role in chlorophyll formation
- essential for proper enzyme activity

COPPER IN SOILS

- Copper is not mobile in soils

Available Cu can vary from 1 to 200 ppm in both mineral and organic soils

- Related to soil pH



availability Cu



Crop	Plant Part Sampled	Time of sampling	Deficient	Low	Sufficient	High	Excessive
					-----ppm Cu-----		
corn	ear leaf	silking	<2.0	2.0-5.0	5.1-20.0	20.1-50	>50
soybeans	top trifoliolate	flowering	<5.0	5.0-9.0	9.1-30.0	30.1-50	>50
wheat	top leaves	boot	<3.0	3.0-5.0	5.1-20.0	20.1-50.0	>50


FUTURE DEVELOPMENT



Interference
usually means
a need for sample
pre-treatment

Requires simple read out
and data interpretation

FUTURE DEVELOPMENT



FUTURE DEVELOPMENT

Areas of Technological Opportunity in Biosensors

Materials for Supports and Electrodes

Nanohybrids, Nanocomposites, Graphene Carbon Nanotubes, Quantum Dots, Magnetic Nanoparticles, Metallic Nanoparticles, Nanowires, Nanorods, SAMs, etc.

Study and Design of Enzymes


Kinetics of Multiple Enzymes, Design of Genetically Modified Enzymes (GMEs), Characterization of GMEs, Identification of Mechanisms of Kinetics of GMEs, etc.

Novel Techniques of Detection

Automated Flow-Injection Analysis directly in Biological Processes, Exploitation of Mediators, Standardization for Guaranteeing No-Toxicity, Encapsulation of Pesticides, etc.



CUPID BIOSENSOR

Design principle





Biosensor example

a Transcriptional



b Translational




CUPID BIOSENSOR

Design principle



Biosensor example



Desain biologi sintetik CUPID *E. coli* ini menggunakan Standar RFC 10, yaitu desain standar perangkatian *BioBrick Parts* dengan menggunakan enzim restriksi dan ligasi yang membentuk situs *scar* atau M

Prefix: **GAATTC**GCGGCCGCT**TCTAGAG**
CTTCCGCGCCGGCGA**AGATCTC**

EcoRI

XbaI

Suffix: T**ACTAGT**AGCGGCCG**CTGCAG**
A**TGATCA**TCGCCGGC**GACGTC**

SpeI

PstI

NotI

Prefix cut with **XbaI**:

...GCT**T[^]CTAG** AG
...CGAA **GATC[^]TC**

Suffix cut with **SpeI**:

T**A[^]CTAG** TAGC...
AT **GATC[^]A**TCG...

SpeI + XbaI

...T**A[^]CTAG** AG...
...AT **GATC[^]TC**...

PENGEMBANGAN RANCANGAN

- Bakteri *E.coli* dibuat dalam bentuk *blotter*, sehingga mudah untuk diaplikasikan di lapangan
- Sistem CUPID dikembangkan menjadi metode deteksi kuantitatif

Level induksi (konsentrasi Cu) → level ekspresi

- Sistem CUPID dikembangkan untuk menjadi detektor senyawa logam lainnya, misalnya merkuri

