

# Improved Lab Notation

## Colour Code

The lab notation is colour coded to facilitate the understanding of the information. It is not strictly necessary that the Notation is colour coded for complete understanding. However, it has been found to be more visually friendly and promote faster understanding of the protocol.

The colours have been chosen to be colour blind friendly. Specially, it has been designed for colour-blind people that have protanopia, which is the most common type of colour blindness, having an incidence of up to 8% in males and 0.5% in females.

Action	Pink
Reagents	brown
Materials	black
Sterility of Protocol	blue
Comments	Yellow

## Upper and Lower Cases

Actions and/or equipment are represented using the notation in upper case. The reagents can be represented in both upper or lower cases depending on how they are known to the science community and how they are likely to be most easily identified.

## Indentations

The lab notation has several indentations to help the understanding to the scientist of what step comes under each action. They can be regarded just like in a programming language or as subheading. Indentations in this lab notation normally follow an action

## Brackets

Brackets are incorporated to help characterize/define the reagent, material or ACTION being performed. Examples of how they are used are shown in the table above.

## Units

For improved standardization, all units used to define quantities must be SI units.

## Signs

The plus sign (+) is to be use only when one material/reagent is being added to another, and be used specifically and only in the order the action takes place.

The equal sign (=) identifies the result of an action and is important to name the mixture formed so that the resultant mixture can then be used more easily when it is used later in the protocol.

e.g: m\_cyl (ddH2O, 100ml) + c\_flask (LB, 2g)  
 = c\_flask(LBmix) - This can only mean that the contents of the measuring cylinder are being added to the contents of the conical flask. The result is a conical flask containing an LB mix

Material	Code	Explanation	Example
Conical Flask	c_flask (cont(s), w, mv)	c_flask(content(solid), weight)  c_flask (content(liquid), v))	
Scalpel	scalpel		
Autoclave Tape	aclv_t		
Water Bath	H2O_Bath		
Distilled water	ddH2O		
Duran Bottle	durant(cont, v)	durant (contents, volume)	durant (LB, 500ml, 100ml)
Eppendorf Tube	epp(l, v)	epp (label, volume using)	epp (lipase, 1ml)
Agar Plate	plate (r, d)	agar plate (resistance (if any), date)	plate (cm, 31/7/2014)
Measuring Cylinder	m_cyl (cont, v)	cylinder (contents, volume using)	m_cyl (H2O, 1000ml, 10ml)
Well Plate	w_well_plate (n)	number of wells_well_plate(number of wells actually used)	24_well_plate (10)
Slides	slide		
Electrophoresis Gel Tray	electro_gel_tray		
Electrophoresis Gel Tank	electro_gel_tank		

Reagents	Code	Explanation	Example
Ethanol	eth		
Alpha cells	$\alpha$ -C		
Agar	agar		
Alpha select chemically competent cells	$\alpha$ -CComp		
Alpha select electrocompetent cells	$\alpha$ -EComp		
DH5 $\alpha$ Cells	DH5 $\alpha$ -C		
DH5 $\alpha$ chemically competent cells (E.Coli)	DH5 $\alpha$ -CComp		
Chloramphenicol	cm (date)		
Lysogeny broth	LB		
Ampicillin	am (date)	ap (date it was made)	
Kanamycin	kn (date)	kn (date it was made)	
Tetracyclin	tc (date)	tc(date it was made)	
Phosphate Buffer Saline	PBS		
DAPI Stain	DAPI		
Glycerol	Glycerol		
Triton	Triton		
Agarose	Agarose		
Tris Acetic acid EDTA	TAE		TAE(1x)
Ethidium Bromide	EtBr		

Command	Code	Explanation	Example
Autoclave	ACLV		
Centrifuge	CENT(cont, speed, t, T)	CENT (contents, speed (g), time, temperature)	CENT (overnight, 11000, 10mins, 4)
Cover	COV (mat)	COV(material)	COV(cotton, foil)
Scale	SCALE (s, w)	weigh out (sample, weight)	SCALE (agarose, 0.5g)
Freeze	FREEZE(s, T)	freeze (sample, temperature)	FREEZE (DH5α_CComp, -80)
Gel electrophoresis	GE (E, L, t, volt)	run a gel electrophoresis (type of gel, ladder, time, voltage)	GE (agarose, 1kb HyperLadder, 1hour, 100V)
Bioincubate	BIOINC (s, t, T, rpm)	incubate (sample, time, temperature, rpm)	BIOINC (overnight, 12hours, 37, 150)
Measure Optical Density	SPEC (s, OD)	spectrometer (sample, optical density used to measure (nm))	SPEC (overnight, 600nm)
Vortex	VORTEX(s, t)	vortex (sample, time)	VORTEX(Xbal construct, 5S)
Discard	DISCARD		
Microwave	MICRWV (p, t)	(power, time)	
Streak	STREAK		
Heat Shock	HEATSHOCK (temp, time, medium)		
Label	LABEL (material, name)	LABEL (what we are labelling, what the label is)	
Dry Block Heating Machine	DBH(cont, temp, t)	DBH(contents(label), temperature, time)	DBH(α-CComp, 42°C, 30s)
Swirl	SWIRL		
Pipette	PIP(cont, v)	PIP (content/label, volume)	PIP(α-C, 55μl)

Set	SET(t)	Allow to set for the minimum duration of t	SET(time)
Remove	REMOVE(i, loc)	REMOVE(item, location) ie: Remove( what?, from where?)	REMOVE(comb, electro_gel_tray)

Conditions	Code	Explanation	Example
Sterile	Ster		
Non-sterile	Non-Ster		