

Date: 01/30/14 People in lab: Levi Palmer

Title: PCR of norCB adding prefix/suffix then electrophoresis

Start Time: 10:15 AM

Purpose: To PCR amplify norCB and add prefix and suffix and confirm on a gel

Protocol: LTM ed. 2 **Exceptions:** 1. 7.5 uL of p. aeru suspension, 2. Final primer conc. <.1uM

Products:

Sample Label	Description	Source Label	Quantity
norCB 1	norCB gene with prefix and suffix	P aeru 7/31	2

Well	1	2	3	4	5	6	7	8
Sample		Ladder	Ladder	norCB1	norCB2			

Results: No product visible

Stop Time: 4:00 PM

Next: Retry with different temperatures