

Results shows that:

1. XynB osmotic shock supernatant has bleaching effect to kraft pulp.
 2. XynB treatment of pulp cannot boost bleaching effect of Sodium hypochlorite, possibly due to high ACTIVATOR concentration.
-
1. Make pulp from kraft, dry to make semi-finished product for further use
 2. XynB solution attained from osmotic shock supernatant
 3. Prepare Sodium hypochlorite solution with an active chloride content of 34.0 g/L - 46.0 g/L as chloride containing chemical decolorizer, designated "ACTIVATOR"
 4. Disperse 0.01g semi-finished product in each tube, add H₂O and XynB solution of each tube according to Table 1
 5. Mix contents of each tube, incubate at 60°C, 200 rpm overnight
 6. For each tube of Group A, centrifuge, discard supernatant. Add H₂O to 1ml, centrifuge, discard supernatant to wash away XynB. Wash for another time, add ACTIVATOR to each tube according to Table 1
 7. Incubate at 60°C, 200 rpm for 30 min
 8. Add 360 µL of each sample from Group A and Group B to a 96 well plate. Each sample is paralleled 3 times.
 9. Use Tanon 3500 gel image system, maximum diaphragm, exposure time = 0.2s to collect image of 96 well plate under white field
 10. Use accessory painting program of windows 8 (i.e. mspaint.exe) to collect central point of each well, use *color picker* tool to measure light intensity profile of each well, data summarized in Table 2

Group A	A1	A2	A3	A4	A5
XynB (mL)	0	0.2	0.4	0.6	0.8
H ₂ O (mL)	Add to 1				
ACTIVATOR(µL)	5				
Group B	B1	B2	B3	B4	B5
XynB (mL)	0	0.2	0.4	0.6	0.8
H ₂ O (mL)	Add to 1				

Table 1

Samples	Average light intensity
A1	28.0
A2	27.7
A3	29.0
A4	29.0
A5	25.3
BLANK	7.3
B1	22.0
B2	22.3
B3	25.3
B4	24.0
B5	27.3

Table 2

