

25 June 2014

To do list

1. Making LB agar

- Tools and Material :
 - ➔ Measure tube
 - ➔ 500 mL bottle
 - ➔ Autoclave band
 - Ps = will change to be black stripe if it's sterile
 - ➔ Stirrer
 - ➔ Scales
 - ➔ Aquades (dH₂O) , take from cabinet (above)
 - ➔ LB powder (@ chemical lab)
 - 3,5 gr / 100 ml
 - ➔ Petri dish
 - ➔ Autoclave
 - ➔ Micropipette
 - ➔ Alumunium foil
- Background :
LB agar is a general medium to grow E. coli on the petri dish.
- Purpose :
Prepare the LB agar for MIC
- Procedure :
 1. Prepare all of the tools and material
 2. Scale the LB powder (3,5 gr for 100 ml agar ; 20 ml/ petri dish)
 3. Mix LB agar powder with dH₂O in 500 mL bottle
 4. Cover it with alumunium foil and sign it with autoclave band
 5. Autoklave it (PS ! Only on 11:00 and 14:00 o'clock)
 - a. Put in pass through box DNA room
 - b. Put dawn beside the autoclave machine.
 6. Let the LB agar cooling (add antibiotic), (PS: prevent the out of antibiotic order)
 7. Pouring at BSL 2 in r. DNA
(PS: password "0001"; clean up the bench with alcohol; sign the petri dish on the beside site)
- Extraordinary occurence
 1. Millard reaction occured on the agar (brownny agar, because of too high temperature or repeating autoclave (protocol.org))
 2. Coagulating agar (need faster moves, solution: autoclave again)

Testified by

Supervisor

(.....)

(.....)

2. MIC

Tabel = variable → control 25, 50, 75, 100, 125, 150, 175

- Background : unknown chloramphenicol MIC
- Purpose : know the MIC of chloramphenicol succinate (equivalent chloramphenicol)
- Tools and material :
 - LB agar
 - petri dish
 - Micropipette
 - Chloramphenicol 1 gr
 - Aquades nuclease free 10ml

- Procedure (calculation)

Tabel

	Control OA OB	25	50	75	100	125	150	175
Medium (ml)	20	20	20	20	20	20	20	
Antibiotic (μL)	0	5	10	15	20	25	30	35

PS : dilute 1gr/10ml → 100 mg/ml → 0,1 mg/ml → 100 μg/ μL

25 mg/ml 20 ml = 500mg → 500 μg/(100 μg/ μL) = 5 μL

PS : single treatment to 25 μg/ml because of agar excess

- Procedure
 1. Mix 1 gr chloram with aquades nuclease free 10 ml in antibiotic tube
 2. Make adequate, take 1,7 mL tube, and save the excess in the refrigerator
 3. Mix with appropriate concentration the we want to the falcon when make the LB agar

Testified by

Supervisor

(.....)

(.....)

1. Mix LB agar powder with dH₂O in 500 ML bottle, with ratio 3,5 gr agar LB for 100 ml dH₂O
2. Covered 500 ML bottle with aluminium foil, pinned autoclave band → autoclave band will be black after the sterilization
3. Autoclave Sterilization Dilution Agar LB in dH₂O
 - Turn on the autoclave machine (autoclave schedule is everyday on 11.00 and 14.00)
 - Import the 500 mL bottle that contain the LB agar to the pass through box.
4. Let the LB agar cooling down (before it rigid)
 - If it was rigid, put LB agar into oven (5 Minutes)
5. Pour the LB agar solution into petri dish
 - Use falcon to amount the volume
 - Add antibiotic using micropipette, then shake agar+antibiotic slowly
 - Do the pouring in BSC

Testified by

Supervisor

(.....)

(.....)

26 June 2014

1. Make LB Broth (stock) 900ml

- Background :
LB broth is used to grow E. coli
- Purpose :
Make the stock if we want to culture the E. coli
- Tools and material :
Tools
 - ➔ 500 mL bottle
 - ➔ Magnetic stirrer
 - ➔ Stirrer
 - ➔ Autoclave band
 - ➔ ScalesMaterial
 - ➔ Trypton 9 gr, yeast 4,5 gr
 - ➔ NaCl
 - ➔ dH₂O
- Procedure :
 1. Mix tripton, yeast, and NaCl with ratio dH₂O (1:0,5:1:0,1) (Trypton:yeast:NaCl:dH₂O)
 2. Pour the dry powder into 500 mL bottle
 - ➔ Powder solution measured by scales, pour it to 500 mL bottle, if there's excess use dH₂O
 3. Add dH₂O till 900 ml
 4. Close the bottle and pinned the autoclave band
 - ➔ Autoclave band will be black
 5. Autoclave broth, before import into DNA room
 6. Let the LB broth cooling down ➔ put into fridge (kulkas -40°C)
- Result :
Transparent yellow LB both

Testified by

Supervisor

(.....)

(.....)

2. RFP Transformation dan spreading

- Background
Resistance chloramphenicol plasmid that contain RFP is inserted to the competent cell, so the transformation will be success or the resistance can be differed from the wild type, with the resistance ability to the antibiotic dosage
- Purpose
E. coli mutant which is resistant of chloramphenicol, was tagged by using fluorescent gene, then ready to isolate and replicate
- Tools
 - ➔ Micropipette
 - ➔ Microtube 1,7 mL tube 1,5 ml
- Material
 - ➔ Dry plasmid iGEM kit 2013 (3,5 g)
 - ➔ dH₂O
 - ➔ ice
 - ➔ SOC
- Procedure
 - ➔ Add 10 µL dH₂O into dry DNA well, solute in 5 minute (up down)
 - ➔ Take 1 µL DNA insert into 50 µL competent cell
- Result

Testified by

Supervisor

(.....)

(.....)

27 June 2014

1. Check and Analyze RFP Transformation

Result :

A lot of single colony that grow on LB agar, but we can't see the difference of transformation result that we expected.

Analyze :

We need to check which is the success RFP transformation to the competent cell. In other hand, there's no competent cell that know wild type transformation, make us must remake the transformation cell with better procedure and modification of IHVCB protocol lab which is appropriated with iGEM protocol.

2. Making LB agar 140 ml (MIC)

Background :

LB agar is used for E. coli growing medium

Purpose :

Make the stock for MIC E. coli

Tools

1. 500 mL bottle
2. Stirrer
3. Autoclave band
4. Scales
5. Alumunium foil
6. Petri dish
7. Falcon
8. micropipete

Testified by

(.....)

Supervisor

(.....)

28 June 2014

AB Aliquots Status

- ➔ there's difference of antibiotic usage
- ➔ 2 no rigid
(save in freezer -40°C)

AB (antibiotic)]

- ➔ boiling point 140° – 150 °C
 - ➔ Freezing point -40° sd -50 °C
- $\Delta T = (kf/b) \times M \times i$

Look the MIC result

30 June 2014

1. Made LB agar plate (220 ml) ➔ 11 plate
2. Replication
 - a. Tools and material sterilization
 - b. LB broth ➔ 40 ml
(Antibiotic 150, made 10 replica, falcon 10)
 - c. Bacteria 150 take 10 colony
(streak 3 times pakai tusuk gigi, put in to falcon + antibiotic, close)
 - d. Streak product on agar replica is incubated in static incubator
 - e. Isolation in the broth ➔ is incubated in static inkubator

Testified by

(.....)

Supervisor

(.....)

1 July 2014

Plasmid isolation

1. Centrifuge 4 ml replica broth with 12.000 g 4°C for 10 minutes
2. Pellet is settled

Testified by

(.....)

Supervisor

(.....)

Testified by

(.....)

Supervisor

(.....)

6 July 2014

Innoculate the big scale transformation

1. Made LB agar 2 petri dish + chloramphenicol 175 (1,4 gr + 40 ml dH₂O)
2. Calculate bacteria colony
 - a. Petri dish 75ng/ml
 $196 \times 4 = 784$
 - b. Petri dish 125 ng/ml
 $281 \times 4 = 1124$
 - c. Petri dish 150 ng/ml
 $320 \times 2 = 640$
 - d. Petri dish 175 ng/ml
 $306 \times 2 = 612$
3. Innoculate big scale (overnight culture)
LB broth stock

Testified by

(.....)

Supervisor

(.....)

7 July 2014

Protocol

Plasmid big scale isolation

- ➔ Overnight culture is putted into falcon 50 ml (2 falcon)
 - ➔ Centrifuge 3500 rpm 10 menit
 - ➔ Discard the supernatant, excess 1 ml (Pelet contain cell + plasmid)
 - ➔ Resuspension pelet with excess s/n and move it into 1,7 mL tube
 - ➔ Centrifuge 12000 rpm, 1 minute ➔ s/n discard
 - ➔ + P1 250 (125) μ L, up down
 - ➔ + P2 250 (125) μ L, invert 4 s.d. 6 x
 - ➔ + solution N 3 350 μ L, invert 4-6 x
 - ➔ Centrifuge 12000 rpm, 10 min ➔ s/n is taken
 - ➔ 850 μ L ➔ move the supernatant to the spin column based nucleic acid purification ➔ centrifuge 12000 rpm, 1 minute
 - ➔ Discard lysate to the column (repeat it twice ➔ centrifuge)
 - ➔ + PB 500 (200) μ L ➔ centrifuge 12000 rpm, 1minute ➔ make the lysate
 - ➔ + PE + Etoh 700 (500) μ L ➔ (centrifuge 12000 rpm, 1min ➔ discard the lysate) ➔ repeat twice
 - ➔ Move column into new microtube 1,7 mL
 - ➔ + 50 μ L 1/3 EB, let it 1 minute
 - ➔ centrifuge 12000 rpm, 1 minute
 - ➔ + 20 μ L 1/3 EB, let it 1 minute
 - ➔ centrifuge 12000 rpm, 1 minute
 - ➔ (Plasmid in the new microtube 1,7 mL)
- PE + Etoh ➔ 1 Vol PE : 4 Vol Etoh
 - 1/3 EB ➔ 1 Vol EB in total 3 vol \sum H₂O
Ex : 30 μ L EB + 60 μ L \sum H₂O

Dilute TE 1x

10 ml ➔ 200 μ L 50 x TE + 9800 μ L

5 ml ➔ 100 μ L 50 x TE + 4900 μ L

Testified by

Supervisor

(.....)

(.....)

$$50 \cdot x = \frac{1}{2} \cdot 100 \text{ ml}$$

$$X = \frac{1}{2} \cdot 100000 \text{ } \mu\text{L} / 50$$

$$X = 10 \text{ } \mu\text{L}$$

Nanodrop

	ng/ μL	Purity	260/230	H ₂ O
DB 1	72,5	2,03	3,5 μL	12,5 μL
DB 2	80,3	1,94	3,2 μL	12,8 μL

E1	PSB 1	3,5 μL		E2	PSB 1	3,2 μL
	dH ₂ O	12,5 μL			dH ₂ O	12,8 μL
	NEBBuf2	2,5 μL			NEBBuf2	2,5 μL
	BSA	0,5 μL			BSA	0,5 μL
	EcoR1	0,5 μL			EcoR1	0,5 μL
	TOTAL	19,5 μL			TOTAL	19,5 μL

Innoculate B. Substilis on petri dish + LB agar

Tools :

1. Petri dish
2. Tips + mcropipet 200 μL

Material :

1. 1,4 gr LB agar
2. 40 ml aquades
3. Bacillus substilis

Testified by

(.....)

Supervisor

(.....)

8 July 2014

Repeat PLASMID DIGESTION

Preparation running electrophoresis

Tools :

1. Tube cabinet
2. Micropipet
3. microtube 1,7 mL 6

Material :

1. Sample 2 PKS 1 dan PKS 2
2. digestion Enzyme @

RFP 720 bp

pSB1C3 2070 bp

	A	B	C
PKS 1	+EcoRI	+PstI	+EcoRI+PstI
O	O	O	O
PKS 2	+EcoRI	+PstI	+EcoRI+PstI
O	O	O	O

	PKS 1 + digestion Enzyme			PKS 2 + digestion Enzyme		
	A (μL)	B (μL)	C (μL)	A (μL)	B (μL)	C (μL)
DNA (250 ng/ μL)	3,5	3,5	3,5	3,2	3,2	3,2
Σ H ₂ O	12,5	12,5	12,5	12,8	12,8	12,8
NEB Buffer 3	2,5	2,5	2,5	2,5	2,5	2,5
BSA	0,5	0,5	0,5	0,5	0,5	0,5
EcoRI	0,5	-	0,5	0,5	-	0,5
PstI	-	0,5	0,5	-	0,5	0,5

Testified by

Supervisor

(.....)

(.....)

8 July 2014

Electrophoresis

	marker	PKS 1	1A	1B	1C	PKS2	2A	2B	2C
9 well									

Tools :

1. Micropipet
2. Glass bottle for agarose + aquades
3. Tips
4. Agarose molding
5. Electrophoresis tool
6. Gel Doc + computer

Material :

1. Agarose 0,8 gr
2. Aquades dH₂O 100ml
3. TAE 0,5
4. Loading dye
5. Ladder (marker)
6. DNA

Testified by

(.....)

Supervisor

(.....)

10 July 2014

- ➔ Nuclease from iGEM is coming ➔ save in the room temperature
- ➔ Plate agar ➔ RFP ➔ moved from incubator to the fridge

11 July 2014

1. Streak and spread E. coli in stab agar from iGEM to LB agar on petri dish
2. Regrow bacteria under petri RFP from incubator
 - Tools
 1. petri dish 4
 2. Micropipet
 3. Falcon 50ml
 4. Tube cabinet
 5. Microtube 1,7 mL tube
 6. Spatel
 7. Tips
 - MATERIAL
 1. Bacteria from iGEM contain plasmid nuclease
 2. Alcohol 70%
 3. LB agar 2,8 gram ➔ 80 ml aquades
 4. Antibiotic ➔ 4 x 35 μ L = 140 μ L
 - Procedure
 1. Scale 2,8 gr LB in chemical lab
 2. Take erlenmeyer, pour the agar into erlenmeyer
 3. Pour 80 ml sterile aquades ➔ cover with aluminium foil
 4. Pinned autoclave band on the aluminium foil
 5. Put erlenmeyer in pass through box
 6. Come in to DNA room ➔ take erlenmeyer LB agar and put into autoclave
 7. Autoclave schedule on 11.00 or 13.00
 8. After autoclave, let the agar cooling down
 9. Do the pouring in BSC
 10. Pour the agar into sterile falcon (20mL)
 11. Add chloramphenicol 30 μ L
 12. Agar LB + antibiotic ➔ pour into petri dish
 13. Let it until rigid
 14. Streak the bacteria contain chloramphenicol
 - a. Divide 10 quadrant (1 quadrant streak 3 times)
 - b. Petri dish packed by parafilm to prevent contamination
 - c. Put into static incubator for 1 day
 - d. Move it to fridge

Testified by

Supervisor

(.....)

(.....)

12 July 2014

Take out bacteria from static incubator → move into fridge.

15 July 2014

Sample		Conc	260/280	260/230
K (-)		-0,9	1,42	0,61
S100	3,4	72,7	1,82	1,33
TolC	6,8	26,6	1,77	1,1
HlyA	3,6	69	1,78	1,34
HlyBD	2,6	96,3	1,59	0,81
RBS	5,2	48,4	1,74	1,34
S119	6,1	40,7	1,8	1,15
D.Term	6	41,6	1,79	1,24
TolC A	4,7	53,2	1,74	1,18
RFP 1		3,2		
RFP 2				
RFP 3				
GFP 4	2,31	13,7		

Testified by

(.....)

Supervisor

(.....)

18 July 2014

Primer Dilution

- Purpose :
Dilute the primer
- Background:
Primer from IDT was in lyophilized primer, so we must dilute it
- Tools and material :
 1. Tube
 2. Micropipette
 3. TE buffer pH 8
 4. ddH₂O
- Procedure :
 - a. Making master stock (100μM)
Make the concentration to be 100 μM, for example, if on the tube is written 20,8 nmol, we must add $20,8 \times 10 = 208 \mu\text{L}$ TE buffer, so we get concentration of primer is 100 μM
 - b. Making working stock (10μM)
Dilution formula
 $V1.M1=V2.M2$
 $100.10=V2.10$
 $V2=100$

Testified by

(.....)

Supervisor

(.....)

19 July 2014

Nuclease

DNA Isolation : Nuclease
Nanodrop

	Σ H ₂ O		DNA	A260	A280	260/280	
0		1,2		0,024	0,002	11,9	0,38
Kontrol		4,5		0,089	0,031	2,84	0,81
6	9,3	37,3	6,7	0,746	0,833	1,95	1,80
4	10,3	43,2	5,7	0,865	0,438	1,97	1,97
3	12,1	68,5	3,9	1,369	0,714	1,92	2,10
2	9,3	37,3	6,7	0,746	0,365	2,04	1,77

DNA Isolation

26 microtube 1,7 mL

6 spin column based nucleic acid purification

Digestion

DNA

Σ H₂O

BSA 0,5

NEB Buffer 3 2,5

EcoRI 0,5

PstI 0,5

Electrophoresis

2 agar

Marker	O	Control	Nuc6	Nuc4	Nuc3	Nuc2	Marker
		-	EcoRI	EcoRI	EcoRI	EcoRI	

Marker	Nuc6	Nuc4	Nuc3	Nuc2	Nuc6	Nuc4	Nuc3	Nuc2	marker
	PstI	PstI	PstI	PstI	EcoRI PstI	EcoRI PstI	EcoRI PstI	EcoRI PstI	

Testified by

Supervisor

(.....)

(.....)

20 July 2014

- PCR Colony running on gell (metode sama dengan elektroforesis biasa)
- Innoculate for DNA isolation in volume 100 ml

Antibiotic concentration that we use: 175 ng/ μ L

Final volume 100ml = 10^5

$M1 = 100 \text{ mg/ml} = 100 \cdot 10^6 / 10^5 \cdot \text{ng}/\mu\text{L} = 100 \cdot 10^3 = 10^5 \text{ ng}/\mu\text{L}$

$V1 \cdot M1 = V2 \cdot M2$

$V1 \cdot 10^5 = 10^5 \cdot 175$

$V1 = 175 \mu\text{L}$

Making mix marker

DNA ladder 10 μ L

TE 1x 70 μ L

6x loading dye 20 μ L

➔ Total 100 μ L

Running list

M – Marker – U – S – U – S – Marker 1 – M2 – M3

U= universal (RFP)

S= substilisin

PCR result analysis

RFP colony + universal primer

If it's fail, the running :

- ➔ Primer was not pinned
- ➔ Wrong primer dilution (TE)
- ➔ PCR is hold on 15 hours in 0°C
- ➔ Minimum amount of colony

Testified by

(.....)

Supervisor

(.....)

21 July

DNA Isolation

K-

Nuc4

Nuc6

Nanodrop	ng/ μ L	A260	A280	260/280	260 620	DNA	Σ H ₂ O
K-	6,0	0,121	0,062	1,94	0,41	41 ~ 16	
Nuc4	148,9	2,979	1,490	2,0	1,83	1,67 ~ 1,7	14,3
Nuc6	101,3	2,026	1,058	1,92	1,52	2,46 ~ 2,5	13,5

BSA 0,5

NEB Buffer 3 2,5

EcoRI 0,5

PstI 0,5

Testified by

Supervisor

(.....)

(.....)

24 July 2014

PCR

Extraordinary occurrence :

1. Failed to make PCR mix

Make culture in LB Broth → falcon 50 ml

Grow B. Subtilis → for DNA isolation

Incubate → shaker incubator

Testified by

(.....)

Supervisor

(.....)

25 July

1. Make TEN reagent for B. Subtilis DNA isolation
2. DNA Plasmid isolation from B. Subtilis by heating
 - a. Take away overnight culture B. Subtilis from shaker incubator
 - b. Take it
 - c. → heating in 100°C for 10 minute PCR machine
 - d. Take it

Testified by

(.....)

Supervisor

(.....)

26 July 2014

Running PCR Result

Extraordinary occurrence :

- Yesterday agar was putted on electrophoresis chamber
- Lefted overnight → too soft, fricted agar, running was too fast band reach the border

Making 2 more agar

1. Substilis
2. Universal

Testified by

(.....)

Supervisor

(.....)

28 July 2014

1. MalS Transformation to the Top10
2. Stocked in -80°C
3. HlyA Transformation in LB broth

Testified by

(.....)

Supervisor

(.....)

30 July 2014

Learn new protocol → PCR : Hot Star !

Protocol → volume of mix PCR different with dream taq method

We use hot star PCR because of failure past PCR method (dimer)

Dimer

- forward Primer with reverse primer
- forward Primer with forward primer

The characteristic of dimer : After PCR running, there will be no underline on the top, but there's a thick sign at the under.

Touch down

- Not only use hot star method to prevent the dimer, but we also do the touch down / temperature optimisation → use 10 column / Gradient of temperature only 50-70°C

Column	Temperature (°C)
1	50
2	50,5
3	51,7
4	53,2
5	55,5
6	58,4
7	61,8
8	64,6
9	66,8
10	68,4
11	69,6
12	70,0

Testified by

Supervisor

(.....)

(.....)

30 July 2014

Protocol : Hot Star (12 reaction)

Volume total	125	μL (PCR Mix)
10 x taq buffer	12,5	μL (PCR Mix)
5 x solution MgCl ₂ 25 Mm	25,00	μL (PCR Mix)
dNTP 10 mM	2,5	μL (PCR Mix)
Primer mix (F t r)	6,25	μL
Taq Hot Star	0,63	μL
Template	6,25	μL
Σ H ₂ O	71,88	μL (PCR Mix)

- Running PCR result → HlyA + Nuclease → PCR banding dreamtaq + hotstar
- Transform HlyA (-PCR), MalS (+PCR), RFP (+)
- Subtilisin → protocol of lysis = make the pellet! not supernatant → PCR dreaming hotstar

50 μL competent cell + genofinterest

+RFP, control (+)

control (-) Mo → LB +antibio, +top 10 Wt

Testified by

(.....)

Supervisor

(.....)

31 July 2014

100 x TAE = 0

→ 500 → 2,5 mL TAE 100X → 500 mL – 2.5 = 497,5

50 x V = 0,5 X 500

V = (0,5 x 500)/50

v = 5 mL → 500 mL

Making LB Agar

300 mL aquades + 10,5 LB agar.

→ make 7 plate

2 RBS + @20 µL ampicillin

2 MalS + @40 µL chloramphenicol

2 RFP + @40 µL chloramphenicol

1 control (-)

Transformation of MalS, RFP, RBS.

MalS = 2013 plate 1 well 4P

RBS = 2014 plate 4 well 1N

RFP = Plasmid from DNA isolation

Testified by

(.....)

Supervisor

(.....)

31 July 2014

Transformation Protocol (Modified by Yuda)

1. Prepare the competent cell 50 μ L
2. Add 1 – 2 μ L DNA from kit
3. Incubate for 60 minutes in ice
4. Heat shock on 42⁰ C for 90 second
5. Incubate in ice 4⁰ C 5 minute
6. Add 200 μ L 50 C media
7. Incubate in shaker 2-3 jam
8. Plating → 200 ng / μ L add 40 μ L (kloramfenikol), 20 μ L (amphicilin)

Amphicilin stock concentration 100 mg/mL

MIC = 1000 mL = 1 mL

$$m1 \times v1 = m2 \times v2$$

$$100 \text{ mg/mL} \times 1 \text{ mL} = 1000 \text{ mL} \times M2$$

$$M2 = 0,1 \text{ mg/mL}$$

$$0,02 \text{ mL} = 20 \text{ } \mu\text{L}$$

$$100 \text{ v1} = 0,1 \times 20$$

$$v1 = 0,02 \text{ mL}$$

$$v1 = 20 \text{ } \mu\text{L}$$

Testified by

(.....)

Supervisor

(.....)

31 July 2014

Running Check Marker ✓
Running HlyA positive ✓
Running Universal works ✓

1 August 2014

- DNA isolation B subtilis Rao Method
- Check transformation
- PCR subtilis
 - malS was failed to be tranformed because we use kit 2013

1 August 2014

- HlyA, malS, HlyDB Transformation
- Isolation B. Subtilis. Rao Method

Testified by

(.....)

Supervisor

(.....)

4 August 2014

1. Boiling dan Sonication

Purpose : Get DNA genome

Boiling :

1. Bacteria is grown in LB and then shake 200 rpm overnight
2. Move it into microtube 125 ml, centrifuge 12.000 rpm 5 minute
3. Take the pellet and add 1/3 TE rinse 3x
4. Heat 100⁰ C 10 minute
5. Sonication
 - 20 seconds do the sonication
 - 10 seconds incubate in ice
 - repeat 6 times for the 5th step
6. Centrifuge 12.000 rpm 4⁰ C , take supernatant

Transformation	Well	Plate
Chez	18 G	1
GFP	13 L	4
R rose	15 A	4
R for LI		

Testified by

(.....)

Supervisor

(.....)

Loading dye 6x

$$m1.v1 = m2.v2$$

$$1.v1 = 6.mL$$

$$= 6 mL$$

$$TE = 5$$

$$1 X = 1/6 Ml$$

$$V1.M1 = V2.M2$$

5 ml LB broth Ampicilin 5 μ L

5 ml LB broth Chloramphenicol 10 μ L

Transformation	Well	Plate
CheZ		
GFP		
RNA rose		
RNA forU		
Cph8	15 P	2
PcyA	13 H	3
Ho7	17 P	3
OmPR promoter	3 I	3
Tol C	40	1

Testified by

(.....)

Supervisor

(.....)

5 August 2014

Nanodrop

No.	Sample	Concentration	260/280	Volume DNA yang dipakai
1	HlyA	49 ng/ μ L	1,97	5,2 μ L
2	HlyD	76,3 ng/ μ L	1,84	3,3 μ L
3	Strong Promoter	84,7 ng/ μ L	1,85	3,0 μ L
4	RBS	43,9 ng/ μ L	1,94	5,7 μ L
5	MalS	100,8 ng/ μ L	1,89	2,5 μ L

Resep Digesti

No.	Sample	DNA	dH ₂ O	NEB Buffer 2	BSA	EcoRI	Total
1	HlyA	5,2	10,8	2,5	0,5	0,5	19,5
2	HlyD	3,3	12,7	2,5	0,5	0,5	19,5
3	Strong Promoter	3,0	13,0	2,5	0,5	0,5	19,5
4	RBS	5,7	10,3	2,5	0,5	0,5	19,5
5	MalS	2,5	13,5	2,5	0,5	0,5	19,5

1. Marker 2,3 RBS 4,5 HlyA 6,7 MalS 8,4 HlyD 10,11 Strong Promoter

Nuclease : M-1,2,3,4,5-M-10,9,8,7,6-M

Isolasi Besar.

1. RBS
2. MalS
3. RFP
4. HlyA
5. HlyD,B
6. Nuclease
7. Strong Promoter

Testified by

(.....)

Supervisor

(.....)

Testified by

(.....)

Supervisor

(.....)

6 August 2014

Nuclease DNA Isolation

Nanodrop Nuclease DNA Isolation

Number	Nanodrop	mg/ μ L	A260	A280	260/280	260/320
1	R-	-	-	-	-	-
2	Nuclease 1	84,5	1,689	0,914	1,85	1,34
3	Nuclease 2	109	2,179	1,162	1,88	1,73

Mixing PCR

Hotstar 12 cycles (substilis)

1. Volume total	125	μ L (12 siklus)
2. 10 x taq buffer	12,5	μ L (PCR Mix appendorf)
3. 5 x solution $MgCl_2$ 25 mM	25,00	μ L (PCR Mix appendorf)
4. dNTP 10 mM	2,5	μ L (PCR Mix appendorf)
5. Primer mix (F + R)	6,25	μ L
6. Taq Hot Star	0,63	μ L + M13 + PKS
7. Template	6,25	μ L
8. Σ H ₂ O	71,88	μ L (appendorf PCR Mix)

Protocol PCR Hotstar taq

Temperature	Time	Cycles
95 ⁰ C (preheat)	15 minute	1 x
Denaturation 94 ⁰ C	0,5 s.d. 1 minute	35-40 x
*Annealing 50-68 ⁰ C	0,5 s.d. 1 minute	
**Elongation 72 ⁰ C	1 minute	
Last heating	10 minute	1 x

*T_m – 5⁰ C

**Elongation DNA waktu : 1 kb/ minute

Storage : 20⁰ C

Testified by

Supervisor

(.....)

(.....)

6 August 2014

Boiling dan Sonikasi Method V.2

1. Spin 500 μL decontaminated sample - for 15 minute
2. Centrifuge 10000 x gmg
3. Discard the supernatant
4. Take the pellet and add 100 μL DNAase free water ($\Sigma \text{H}_2\text{O}$)
5. Incubate for 20 minutes 95°C and boil it in waterbath
6. Sonication 15 minute
7. Centrifuge 5 minutes full speed 13000 \rightarrow take the supernatant
8. Use 5 μL supernatant for PCR template

MIXING PCR

Dream taq 5 cycles

Testified by

(.....)

Supervisor

(.....)

6 August 2014

Repeat the inoculation and streak replica of transformation 1st

Repeat the inoculation and streak replica of transformation 2nd

List colony transformation 1st :

1. MalS ✓
2. RBS ✓
3. RFP ✓
4. Nuklease → have been isolated today
5. Strong promoter ✓
6. HlyA ✓
7. HlyB,D ✓

List colony transformation 2nd:

1. CheZ ✓
2. GFP ✓ *Amphicilin
3. RNA rose ✓
4. RNA forU ✓
5. Cph8 ✓
6. PcyA ✓
7. HO7 ✓
8. OmpR promoter ✓
9. TolC ✓

Testified by

(.....)

Supervisor

(.....)

Streak replica → LB agar

Plate	Plasmid	Antibiotic
1	RRS + GFP	chloramphenicol
2		
3		
4		
5		
6		
7		
8		

Testified by

(.....)

Supervisor

(.....)

7 August 2014

Isolation result

Nanodrop : HlyA 1-5, Mol 5 1-5, RFP 1-5, HlyB,D 1-5

Number	Sample	DNA	Concentration (ng/ μ L)	Σ H ₂ O	A260	A280	260/280	260/230
1	K (-)		0,2		0,0004	0,0004	1,07	1,88
2	HlyA 1	26,5	9,4	6,6	0,188	0,094	1,89	2,38
3	HlyA 2	0,8	307,6		6,151	2,601	2,20	2,37
4	HlyA 3	12,5	20,0	3,5	0,400	0,207	1,93	1,86
5	HlyA 4	9,05	27,6	6,95	0,552	0,261	2,12	1,46
6	HlyA 5		0,1		0,002	0,004	0,46	0,1
7	MalS 1	6,5	38,4	9,5	0,768	0,383	2,00	1,94
8	MalS 2	3,25	76,8	12,75	1,536	0,815	1,88	2,22
9	MalS 3	3,23	77,2		1,544	0,817	1,89	1,79
10	MalS 4	8,38	29,8	7,62	0,596	0,310	1,92	1,79
11	MalS 5	5,19	48,1	10,81	0,963	0,499	1,93	1,61
12	RFP 1	4,1	59,8	11,9	1,195	0,592	2,02	2,17
13	RFP 2	1,9	177,5		3,550	1,846	1,92	2,20
14	RFP 3	2,36	105,8		2,115	1,092	1,94	2,04
15	RFP 4	1,90	131,3	14,1	2,626	1,402	1,87	2,22
16	RFP 5	3,99	62,6	12,01	1,251	0,645	1,94	2,17
17	HlyBD 1	5,89	42,4	10,11	0,874	0,427	1,98	1,38
18	HlyBD 2	8,25	30,3		0,605	0,304	1,46	1,77
19	HlyBD 3	9,46	26,4		0,527	0,276	1,41	1,94
20	HlyBD 4	18,3	13,6		0,272	0,130	2,09	1,65
21	HlyBD 5	6,52	38,3	9,48	0,767	0,395	1,94	1,52

Testified by

Supervisor

(.....)

(.....)

8 August 2014

Num	Sample	DNA	Concentration (ng/ μ L)	Σ H ₂ O	A260	A280	260/280	260/230
1	K (-)		0,1		0.002	0,010	0,20	1,22
2	St. Promoter 1	8,5	29,4		0,589	0,291	1,49	1,87
3	St. Promoter 2	9,15	37,3		0,746	0,381	1,95	1,84
4	St. Promoter 3	8,79	28,6		0,573	0,274	2,09	1,89
5	St. Promoter 4	2,63	94,7	13,37	1,894	0,995	1,90	1,99
6	St. Promoter 5	9,208	59,4	11,7	1,188	0,617	1,92	1,82
7	OMPR 1	7,12	35,1	8,8	0,702	0,352	1,99	1,85
8	OMPR 2	4,1	60,8	11,9	1,216	0,635	1,91	2,05
9	OMPR 3	6,8	36,5		0,730	0,364	2,01	1,94
10	OMPR 4	7,7	32,2		0,644	0,319	2,02	1,76
11	OMPR 5	8,3	30,1		0,602	0,302	2,00	1,80
12	HO7 1	4,17	59,9	11,83	1,199	0,626	1,91	2,14
13	HO7 2	9,6	25,8		0,516	0,255	2,02	1,95
14	HO7 3	11,6	21,5		0,430	0,192	2,23	1,92
15	HO7 4	5,28	47,3	10,72	0,445	0,491	1,93	2,04
16	HO7 5	15,7	15,9		0,317	0,141	2,25	1,87
17	RNA rose 1	5,9	42,1	10,1	0,842	0,411	2,05	1,61
18	RNA rose 2	4,72	52,9	11,28	1,059	0,586	1,98	1,91
19	RNA rose 3	8,74	28,6		0,572	0,260	2,20	1,79
20	RNA rose 4		27,9		0,557	0,262	2,13	1,64
21	RNA rose 5		13,2		0,265	0,109	2,42	1,38

Note : (Tol C 5, CheZ 4, CheZ 5, Rose 5) was missing, CPH8 only 2

Testified by

(.....)

Supervisor

(.....)

No.	Sample	DNA	Concentration (ng/ μ L)	Σ H ₂ O	A260	A280	260/280	260/230
1	K (-)		0,4		0,807	0,0002	-4,35	-2,93
2	RBS 2	14,97	16,7	1,03	0,334	0,175	1,90	1,61
3	RBS 3	3,931	63,6	12,069	1,272	0,681	1,87	1,81
4	RBS 5	9,58	26,1	6,42	0,522	0,270	1,93	1,81
5	RBS 4		102,2				1,83	1,09
6	TolC 1	9,857	27,3	6,843	0,545	0,286	1,91	1,52
7	TolC 2	3,444	72,6	12,556	1,452	0,770	1,88	2,00
8	TolC 4	4,386	57,0	11,6			1,93	1,68
9	GFP 4	2,31	108,2	13,7	1,141	0,589	1,86	2,15
10	GFP 3	4,55	55,0	11,45	1,101	0,568	1,94	1,81
11	GFP 5	5,28	47,3	10,72	0,946	0,490	1,93	1,85
12	RNA Rose 1	10,4	24,0	5,58	0,480	0,252	1,91	1,73
13	CheZ 1	3,99	62,7	12,0	1,254	0,664	1,89	2,07
14	CheZ 2	4,8	52,1	11,2	1,043	0,552	1,89	1,90
15	CheZ 3	4,31	58,0	11,69	1,160	0,626	1,85	1,78
16	RNA Rose 3	3,0	50,4	11			1,92	1,76
17	RNA Rose 5		37,1				2,01	1,67
18	Cph8 1	2,56	97,6	13,44			1,91	1,95
19	Cph8 2	9,3	26,9	6,7			1,08	1,11
20	Cph8 3	3,28	76,2	12,7			1,91	1,94
21	Cph8 4	3,7	67,3	12,3			1,90	1,68
22	Cph8 5		40,2				1,94	1,77
23	PcyA	4,2	59,5	11,8			1,91	1,89
24	PcyA	4,93	50,7	11,07			1,95	1,80
25	PcyA	7,08	35,3	8,92			2,05	1,04
26	PcyA	4,17	59,9	11,83			1,92	1,78
27	PcyA	4,1	61,1	11,9			1,92	1,98
28	RFP 4	1,81	130,1	14,2			1,80	1,93
29	RFP5	2,16	125,7	13,84			0,94	0,61

Testified by

(.....)

Supervisor

(.....)

8 August 2014

G Blocks

Ligase

1. Dilute TE 20 μ L
 - Concentration 200 ng/20 μ L = 10 ng/ μ L
2. Take 10 μ L (100 ng) \rightarrow double digest with Eco RI and pstI

Protocol	gBlocks	Vector (RFP)
DNA	10 μ L	250 ng \rightarrow 4 μ L
dH ₂ O	6 μ L	12 μ L
NEB2	2,5 μ L	2.5 μ L
BSA	0,5 μ L	0.5 μ L
EcoRI	0,5 μ L	0.5 μ L
PstI	0,5 μ L	0,5 μ L
Digested	20 μ L	20 μ L
	Inserted	Vector

3. PCR machine
37⁰ C 30 minute + 80⁰ C 20 minute
4. Take 4 μ L (20 ng) digested insert (untuk ligasi)
Take 4 μ L (50 ng) digested vector (untuk ligasi)

Protocol

digested vector	4 μ L
digested insert	4 μ L
T4 buffer	1 μ L
T4 ligase	0,5 μ L
H ₂ O	<u>0,5 μL</u> +
	10 mL

5. PCR machine \rightarrow 23⁰C 2 hours (16⁰C 16 hours) + 80⁰C 20 minutes
6. Take 4 μ L \rightarrow for transformation

Testified by

Supervisor

(.....)

(.....)

8 August 2014

Lab Plan 8 August 2014

1. Continue 8 genes isolation → 40 falcon
→ There's only 7 genes, except Nuclease, so it must be 35 falcons, but 3 had missed when washing, so there are 32 falcons left
2. Nanodrop in the next 2 pages
DNA isolation result → aliquot → digestion
3. Running electrophoresis → digestion result
→ 2 @1 part → take 1 (inoculate in LB Broth 100 ml)

Making LB Broth 1000 ml

Plan priority parts 2

- MalS
 - Strong Promoter
 - HlyA
 - HlyBD
 - RBS
 - TolC
 - RFP
4. PCR → touchdown MalS (DT) → running electrophoresis
 5. Restriction ligation of peptide 1018 + RFP vector

Testified by

(.....)

Supervisor

(.....)

Testified by

(.....)

Supervisor

(.....)

9 August 2014

1. DNA isolation MalS and Strong Promoter (after 12 o'clock → 2 hours)
2. Running digestion result (8 August 2014) (11 s.d. 12)
3. Digest DNA isolation result 2nd (11 s.d. 12)
4. Transform ligation result : peptide 1018 gBlocks (after 12 o'clock → 3 hours)
5. *B. subtilis* isolation (3 hours)

Running (8 August 2014)

Marker MalS MalS St. prom St. prom 1 2 3 4 m

Running (9 August 2014)

Marker(old). HO7₁ HO7₄ RNAFU₁ RNAFU₂ RFP₄ RFP₅ OMPR₁ OMPR₂ HlyB₁ HlyB₅ HlyA₃ HlyA₄

Marker(old) MalS₂ MalS₄ St. prom₄ St. prom₅ (PCR MalS₁) (PCR MalS₃) (PCR MalS₅) Marker(bar)

PCR Hifi PFX mix	(Nuclease) μ L	(MalS) μ L
10 x PFX	1,3	1,3
10 X PCR enhance	1,3	1,3
10 mM dNTP	0,39	0,39
50 mM MgSO ₄	0,26	0,26
Primer 10 μ L	0,39 (VF2 & VR)	0,39 (VF2 & MalS reverse)
Platinum PFX	0,1	0,1
Σ H ₂ O	8,26	8,26
Template	1,3	1,3

Program PCE machine

PCX Hifi

Testified by

(.....)

Supervisor

(.....)

9 August 2014

Digestion, Nanodrop
PCR Hifi Nuclease

Program protocol PFX Hifi (Nuclease, MalS) → Primer Universal Tm 58

	Temperature	Time	Cycles
Initiation Denaturation	94 ⁰ C	3 minute	1 x
Denaturation	94 ⁰ C	30 second	40 x
Annealing	opt(53 ⁰ C)	30 second	
Elongation	72 ⁰ C	1000 bp/min (2mm)	
Last Elongation	72 ⁰ C	10 minute	1 x

Testified by

Supervisor

(.....)

(.....)

11 August 2014

1: - MalS DNA isolation

- Subtilis PCR
- Digestion

2 : - Nanodrop

- make the Plate → spread transform
- Transform peptide 1018, added parts

Transformation kitplate 2014

Parts	Plate	Well	Antibiotic
Strong Promoter 100	4	17D	Amphicilin
Strong Promoter 119	4	17B	Amphicilin
Peptide 1018	-	-	Chloramphenicol
Double terminator	3	3F	Chloramphenicol
Lambda represor	2	4B	Chloramphenicol
Lambda Promoter	4	5L	Amphicilin

→ PCR DT → CDLony PCR

HlyA, RFD → primer = VF2, VR

MalS → VF2, MalS R

→Running PCR

Testified by

(.....)

Supervisor

(.....)

Testified by

(.....)

Supervisor

(.....)

11 August 2014

Electrophoresis

Marker; RNA Rose₂; RNA Rose₃; RBS B₂; RBS B₃; RBS B₅; GFP₃; GFP₄; PcyA₄; PcyA₅; CPH8₃; CPH8₄

Bacillus Subtilis

Column 3, 5, 7, 9, 11

PCR date 14 August 2014

List

- I. CheZ1, CheZ2, CheZ3 | GFP1, GFP2, GFP3 | RFP1, RFP2, RFP3 | RNAR1, RNAR2, RNAR3 |
720 706
- II. LR1, LR2, LR3 | OMPR1, OMPR2, OMPR3 | HB11, HB12, HB13 | CPH81, CPH82, CPH83 |
RNAF1, RNAF2, RNAF3 |
- III. PCY81, PCY82, PCY83

Testified by

Supervisor

(.....)

(.....)

12 August 2014

To do List Lab Target

1. making PQE 80L stock has done
2. running elektroforesis taq Hifi (11 August 2014)
3. streak bacteria target 6
 1. RB5 (Amphi)
 2. HlyB,D (Chloram)
 3. TolC (Chloram)
 4. HlyA (Chloram)
 5. Double Terminator (Chloram)
 6. St. Promoter 100 (Amphi) + 7 St. Promoter 119 Amphi 8 RFP
 7. MalS (chloram)
 8. Nuclease (chloram)
4. Isolate subtilis → PCR (Hotstar) ; Subtilisin
5. mix PCR for PCR colony DT → Hifi

Testified by

(.....)

Supervisor

(.....)

12 August 2014

Prepare agar in petri dish → make paper background

Make for 12 μL $\Sigma \text{H}_2\text{O}$ + bacteria colony

Innoculate 1 colony, and shake until the bacteria in $\Sigma \text{H}_2\text{O}$ 12 μL

RFP → is taken 3 colonies

1. RBS
2. TolC
3. HlyB
4. RFP
5. HlyA
6. D. Term
7. S. prom 100
8. S. prom 119
9. MalS
10. Nuclease

(Total : 30 PCR tubes contain template)

Testified by

(.....)

Supervisor

(.....)

13 August 2014

Running Electrophoresis (wian + hanifi)

Marker
RBS 1
RBS 2
RBS 3
TolC 1
TolC 2
TolC 3
HlyB 1
HlyB 2
HlyB 3
HlyA 1
HlyA 2
Marker

Marker
HlyA3
Double terminator 1
Double terminator 2
Double terminator 3
S100 1
S100 2
S100 3
S119 1
S119 2
S119 3
RFP
Marker

PCR DT
PCR Hifi
TolC 3
RBS 3
St. Prom 119 3
St. Prom 100 1
Double terminatot 1
HlyA 3
HlyB,D 1

Testified by

(.....)

Supervisor

(.....)

		B4		B4		B4		B4		B4	
			B2				B2				
		KoI		KoI		KoI		KoI		KoI	

Dream Taq

PCR mix dream taq → 50 µL = 5 reaction

10 x Dreamtaq Buffer	5 µL
2 . 10 mM dNTP	5 µL
Primer Forward	0,6
Primer Reverse	0,6
Dream taq polymerase	1,25
Σ H ₂ O	...

Program PCR dreamtaq	Temperature	Times	Cycle
Initial denaturation	95 °C	5 minute	1
Deanturation	95 °C	30 s	25-40
Annealing	(opt) 53 °C	30 s	
**Extension/elongation	72 °C	1000 bp/min	
Final extension	72 °C	10 minute	1

* T_m = X °C

** Annealing = X °C - 5 °C

Testified by

Supervisor

(.....)

(.....)

14 August 2014

Lab Target

1. DT
2. MalS
3. Nuclease
4. RFP
5. OMPR
6. GFP
7. CheZ
8. RNA rose
9. Cph 8
10. PcyA
11. H07
12. RNA for U

Testified by

(.....)

Supervisor

(.....)

15 August 2014

Electrophoresis PCR DT result
Take out the transformation

HlyBD DNA isolation
Nanodrop

Sample	Concentration (ng/ μL)	A60	A280	260/280	260/230
HlyBD 1	127,8	2,55	1,225	2,09	2,08
HlyBD 2	317,9	6,359	2,969	2,14	2,25

Testified by

(.....)

Supervisor

(.....)

16 August 2014

1. prepare Vibrio growth medium
2. do the DNA isolation Peptide 1018 & T5 promoter
3. nanodrop → digestion

1,7 mL tube used : 8

Spin column based nucleic acid purification used : 8

Testified by

(.....)

Supervisor

(.....)

	1018 1	1018 2	1018 3	1018 4	1018 5	T5 1	T5 2	T5 3	HlyB 1	HlyB 2
DNA										
$\sum \text{H}_2\text{O}$										
NEB 3	2,5	2,5	2,5	2,5	2,5	2,5	2,5	2,5	2,5	2,5
BSA	0,5	0,5	0,5	0,5	0,5	0,5	0,5	0,5	0,5	0,5
PstI	0,5	0,5	0,5	0,5	0,5	0,5	0,5	0,5	0,5	0,5

Sample	Concentration	A260/280	260/230
1018 1	97,6	1,93	1,81
1018 2	338,1	1,95	1,88
1018 3	96,0	1,97	1,49
1018 4	96,7	1,86	1,15
1018 5	149,5	2,01	1,84
T5 1	178,9	1,99	1,80
T5 2	107,7	1,93	1,57
T5 3	76,9	1,86	1,24

Peptide 1

- Take $\rightarrow 2,56\mu\text{L}$
- $\sum \text{H}_2\text{O} \rightarrow 13,45 \mu\text{L}$

Peptide 3

- Take $\rightarrow 2,6 \mu\text{L}$
- $\sum \text{H}_2\text{O} \rightarrow 13,4\mu\text{L}$

Peptide 4

- Take $\rightarrow 2,6 \mu\text{L}$
- $\sum \text{H}_2\text{O} \rightarrow 13,4 \mu\text{L}$

NEB Buffer 3 $\rightarrow 2,5 \mu\text{L} \times 6 = 15 \mu\text{L}$

BSA $\rightarrow 0,5 \mu\text{L} \times 6 = 3 \mu\text{L}$

PstI $\rightarrow 0,5 \mu\text{L} \times 6 = 3\mu\text{L}$

Testified by

(.....)

Supervisor

(.....)

DreamTaq

10x buffer	5 µL
dNTP	5 µL
Fwd	5 µL
Reverse	5 µL
H ₂ O	29,75 µL
DreamTaq	0,25 µL
Template	2 µL
TOTAL	50 µL

PCR Test

1	Subs + bac. Subs	10
2	Univ + plasmid kit	5
3	Univ + red codon	5
4	Subs + no DNA	10
5	Univ + no DNA	5
6	M13 + PKS	7

10x buffer	5 µL X 6
dNTP	5 µL X 6
Fwd	5 µL
Reverse	5 µL
H ₂ O	29,75 µL X 6
DreamTaq	0,25 µL
Template	

Testified by

(.....)

Supervisor

(.....)

	Peptide 1018	Linear pSB1C3	RFP
DNA	4 µL	10 µL	
dH ₂ O	12 µL	6 µL	
NEB 3	2,5 µL	2,5 µL	2,5 µL
BSA	0,5 µL	0,5 µL	0,5 µL
EcoRI	0,5 µL	0,5 µL	0,5 µL
PstI	0,5 µL	0,5 µL	0,5 µL
DpnI	-	0,5 µL	-
Total	20 µL	20,5 µL	20 µL

Testified by

(.....)

Supervisor

(.....)

19 August 2014

PCR Hifi MalS Nuclease

Testified by

(.....)

Supervisor

(.....)

20 August 2014

Ligase digestion result

Transform ligation result → try BL21 CP – RFP

Consultation and discussion with dr. Budi

Heat MaS, Nuclease

Assembly Step1 (and transform)

1. S100 + RBS (K)
2. S119 + RBS (K)
3. T5 + RBS (A)
4. 1018 + HlyA (A)
5. TolC + HlyBD (A)
6. S100 + Rose (A)
7. RBS + GFP (K)
8. RFP Amphi Backbone (A)
9. MaS + HlyA (no transform)
10. Nuclease + HlyA (no transform)
11. Subtilisin + HlyA (no transform)
12. RBS + 1018 (A)
13. RBS + CheZ (A)
14. RFP + GFP Backbone (A)

Testified by

(.....)

Supervisor

(.....)

21 August 2014

Isoamyl alcohol method → genome & subtilis isolation

Transformation result → incubator

10 → a bit growth of colony, small diameter, RFP was not expressed

13 → idem

BL21 CP → not turn red

TolC → HlyBD ✓

10 plate (3 chloram + 7 amphi)

Sample	Concentration	A260	A280	260/280	260/230
Blank	-0,1	-0,002	-0,011	0,15	-0,2
MalS	39,4	0,788	0,540	1,46	0,49
MalS	40,5	0,810	0,561	1,44	0,49
Nuclease	24,8	0,497	0,331	1,50	0,55
Nuclease	24,9	0,497	0,329	1,51	0,55

Testified by

Supervisor

(.....)

(.....)

22 August 2014

Sp100-RBS 3
S119-RBS 6
RBS-GFP 1
TolC-HlyBD 3
T5-RBS 3
1018-HlyA 3
Sp100-Rose 2
RFP-GFP 3
RFP-Amphi 1
RBS –cheZ 3

TolC

Testified by

(.....)

Supervisor

(.....)

23 August 2014

Digestion

- 1018 (3) DNA 2,6 ul + H₂O 13,4 ul
- RBS Hifi (48,4 ung/ul) DNA 5,2 ul + H₂O 10,8
- Linear Amp DNA 10 ul + H₂O 6 ul
- MalS Hifi new (130 ng/ul) DNA 2 ul + H₂O 14 ul
- Nuclease Hifi new (133 ng/ul) DNA 1,9 ul + H₂O 14,1 ul
- HlyA hifi DNA 3,6 ul + H₂O 12,4 ul
- RFP 2 (82,8 ng/ul) DNA 3 ul + H₂O 13ul

Nanodrop Purified PFX Hifi MalS Nuclease

Sample	Concentration (ng/ul)	A260	A280	260/280	260/230
Blank	0,5	0,01	0,013	0,72	0,61
MalS (1)	130,3	2,605	1,391	1,87	2,15
MalS (2)	133,6	2,672	1,430	1,87	2,07
Nuc (1)	133,8	2,675	1,441	1,86	2,1
Nuc (2)	138,9	2,777	1,492	1,86	2,09

PFX Hifi, success! (DNA concentrarion is high, Running electrophoresis: suitable length band)

LAB PLAN IGEM 2014

Making LB Broth stock 1L = devide 2 (500,500)

running PFX Hifi MalS, nuclease

Purify PFX Hifi MalS, nuclease

Running electrophoresis PCR colony

Confirmation → DNA isolation

DNA isolation

Nanodrop isolation

Digest, Ligase, Transform (3 Transformation)

1. P1018 – RBS
2. MalS – HlyA
3. Nuc – HlyA

DNA isolation list

- TolC HlyBD 2
- TolC HlyBD 3
- P1018 HlyA 2
- P1018 HlyA 3
- RFP Amphi 1
- Sp100 RBS 2
- Sp100 RBS 3
- T5 RBS 2
- T5 RBS 3
- Sp100 Rose 1

Testified by

Supervisor

(.....)

(.....)

- Sp100 Rose 2
- RFP GFP 1
- RFP GFP 3
- RBS CheZ 1

Isolation result

Nanodrop

Sampel	Conc (ng/ul)	A260	A280	260/280	260/230
Sp100-Rose 1	90,5	1,809	0,895	2,02	1,9
Sp100-Rose 2	63,7	1,275	0,644	1,98	1,83
T5-RBS 2	70,9	1,419	0,713	1,99	1,89
T5-RBS 3	93,3	1,867	0,917	2,04	1,96
RFP-GFP 1	111,8	2,237	1,110	2,01	1,89
RFP-GFP 3	61	1,221	0,622	1,96	1,77
	66,8	1,336	0,746	1,79	2,29
RBS-CheZ 1	88,9	1,779	0,885	2,01	1,92
TolC-HlyBD 2	125,8	2,517	1,253	2,01	2,08
TolC-HlyBD 3	135,8	2,716	1,34	2,03	2,06
	133,1	2,661	1,31	2,03	2,04
P1018-HlyA 2	165,8	3,316	1,66	2	1,95
	150,6	3,012	1,518	1,98	1,92
P1018-HlyA 3	125,2	2,504	1,28	1,96	2,01
	111,1	2,22	1,135	1,96	1,93
RFP Amphi 1	176,3	3,525	1,84	1,92	2,11
	176,4	3,529	1,849	1,91	209
RFP cheZ 3	88,6	1,772	0,901	1,97	1,99

Running Elektroforesis

Well 1

1. S100-Rose 1
2. S100-Rose 2
3. T5-RBS 2
4. T5-RBS 3
5. RFP-GFP 1
6. RFP-GFP 3
7. RBS-CheZ 1

Well 2

1. RBS-CheZ 3 → 2000
2. 1018-HlyA 2 → 2000-2500
3. 1018-HyA 3 → 2000-2500
4. RFP amphi 1 → 3000
5. TolC-HlyBD 2 → 2000-2500
6. TolC-HlyBD 3 → 2000-2500

SP Rose

S100-PCR

Rose-isolation

T5-isolasi

RBS-Hifi

CheZ-isolation

Testified by

Supervisor

(.....)

(.....)

1018-isolation
HlyA-PCR
TolC-PCR
HlyBD-isolation

Testified by

(.....)

Supervisor

(.....)

24 August 2014

DNA isolation

	Conc	260/280	260/230
S100- RBS 2	43,4	1,86	1,76
RBS-GFP 1	73,6	1,90	1,87
S100-rbs 3	70,5	1,97	1,89
Blank	0,7	0,6	0,45

DNA isolation

	Conc	260/280	260/230
RFP Chloram 1	46,7	1,69	0,79
RFP Chloram 2	33,4	1,68	0,61
RFP Amphi 1	32,8	1,49	0,44
RFP Amphi 2	46,6	1,56	0,47
GFP 1	27,3	1,76	0,7
GFP 2	24,3	1,87	0,79

Testified by

(.....)

Supervisor

(.....)

25 August 2014

Digestion

PART A (EcoR1+Spe1)

- RBS-CheZ 3
- S100-ROSE 1
- T5-RBS 3
- TolC Hifi

PART B (Xba1+Pst1)

- 1018
- GFP
- 1018-HlyA 2

PART B (Nhe1+Pst1)

- HlyBD 1

Mix PART A

- NEB Buffer 3 10ul
- BSA 10X 2ul
- EcoR1 2ul
- Spe1 2ul

Mix PART B

- NEB Buffer 3 7,5ul
- BSA 10X 1,5ul
- Xba1 1,5ul
- Pst1 1,5ul

Single Digest (Pst1) C

- S100-RBS 2
- RBS-GFP 1
- S100-RBS 3
- NEB 7,5
- BSA 1,5
- Pst1 1,5

	Template	DNA	H ₂ O
A	RBS-CheZ	2,9	13,1
B	S100-Rose 1	2,8	13,2
C	T5-RBS 3	2,7	13,3
D	TolC Hifi	11,3	4,7
E	1018 3	2,6	13,4
F	GFP 4	2,4	13,6
G	1018 HlyA 2	1,7	14,3
H	HlyBD 1 (isolasi)	3,9	10,1
I	RFP 2 isolasi kloram	3	13
J	RFP Amphi kecil 1	1,5	14,5
K	S100-RBS 2	5,8	10,2
L	RBS-GFP	3,4	12,6
M	S100-RBS	3,6	12,4

Testified by

Supervisor

(.....)

(.....)

27 August 2014

Assembly Step 2

Digest Step 2

A	B
S100-RBS 2	1018
S100-RBS 3	1018
S100-RBS	Nuc HlyA 2
S100-RBS	Nuc HlyA 3
S100-RBS	MalS HlyA 2

	$\sum H_2O$	DNA
S100-RBS 2	10,2	5,8
S100-RBS 3	12,4	3,6
Nuc HlyA 2	12,6	3,4
Nuc HlyA 3	14,3	1,7
MalS HlyA 2	12,7	3,3
RFP Amphi 1	14,6	1,4

Mix Digest

	A	B	C
NEB 4	7,5 ul	10 ul	2,5 ul
BSA 10X	1,5 ul	2 ul	0,5 ul
EcoR1	1,5 ul	-	0,5 ul
Spe1	1,5 ul	-	-
Xba1	-	2 ul	-
Pst1	-	2 ul	0,5 ul

Ligation Step 2

- S100-RBS – Nuc-HlyA I RFP Amph 5 1 ul Buffer BSA
- S100-RBS – MalS-HlyA I RFP amph 2,0 0,5 ul T4 Ligase
- S100-RBS – 1018 I RFP Amph 2,5 ul $\sum H_2O$

Digestion Step 1

	DNA	$\sum H_2O$
TolC	6,9	9,1
HlyBD	5,9	10,1
HlyBD	3	5

A

- NEB 1 2,5
- BSA 10X 0,5
- EcoR1 0,5
- Spe1 0,5

B

- Neb 1 2,5
- Bsa 10x 0,5

Testified by

Supervisor

(.....)

(.....)

- Nhe1 0,5
- Pst1 0,5

C

- NEB 1 1,3
- BSA 10X 0,3
- Pst1 0,25

D

- NEB 1 1,3
- BSA 10X 0,3
- Nhe1 0,25

Next Ligation Step 1
TolC-HlyBD

Running

M-TolC-HlyBD total-HlyBD Nhe1-HlyBD Pst1-Nuc HlyA 2-Nuc HlyA 3-Nuc hifi-MalS HlyA-M

Testified by

(.....)

Supervisor

(.....)

28 August 2014

Digestion check Nuc result

	DNA	Σ H ₂ O
1. Hifi HlyA	3,6	12,4
2. HlyA 4	9,1	6,9
3. Nuclease (stock 1)	3	13
4. Nuclease Hifi	1,9	14,1

PCR colony Electrophoresis

A1 S100-Rose II GFP 1

A2 S100-Rose II GFP 2

B1 RBS_CheZ II GFP 1

B2 RBS_CheZ II GFP 3

C1 S100-Rose II 1018 1

C2 S100-Rose II 1018 2

.....

C10 S100-Rose II 1018 10

D1 T5-RBS II 1018 1

D2 T5-RBS II 1018 2

.....

D10 T5-RBS II 1018 10

Testified by

(.....)

Supervisor

(.....)

29 August 2014

Running Elektroforesis

Marker – RFP Backbone 2 – HlyA –Nuclease –Nuc HlyA 3 – RFP Amphi 1

Nanodrop DNA isolation 28/8/14 (chloramphenicol)

	Concentration	260/280	260/230
T5-RBS II 1018 2	252,2	1,91	1,29
T5-RBS II 1018 3	60,2	2	1,19
T5-RBS II 1018 4	66,4	1,97	1,03
T5-RBS II 1018 7	82	2,06	1,46
S100-Rose II 1018 2	56,1	1,99	1,56
S100-Rose II 1018 5	47,1	2,07	1,13
S100-Rose II 1018 6	39,8	2,09	1,01
S100-Rose II 1018 8	72,7	2,02	1,35
RBS – CheZ II GFP 1	57,7	2,04	1,19
RBS – CheZ II GFP 2	17	2,22	0,73
RBS – CheZ II GFP 3	20,5	1,14	0,72
S100-Rose II GFP 1	151	2,19	1,83
S100-Rose II GFP 2	83	2	1,55

Testified by

Supervisor

(.....)

(.....)

30 August 2014

Running Isolation result

Control 1018-T5 RBS 1018 2-T5 RBS 1018 3-T5 RBS 1018 4-T5 RBS 1018 7-S100 Rose 1018 2-S100
Rose 1018 5

S100 Rose 1018 6-s100 Rose 1018 8-S100 Rose GFP 1-S100 Rose GFP 2-RBS CheZ GFP 1

Double digest 1 September 2014 ½ Reaction

			½ DNA	½ H ₂ O
1. T5-RBS - 1018 3	NEB 4	7,5 ul	2,1	5,9
2. T5-RBS – 1018 4	BSA 10X	1,5 ul	1,9	6,1
3. S100-Rose - 1018 8	EcoR1	1,5 ul	1,8	6,1
4. S100-RBS – 1018 2	PstI	1,5 ul	1,5	6,5
5. RBS – CheZ-GFP 1	@ 2ul		2,2	5,8

Running

1. T5-RBS (digestion)
2. T5-RBS 1018 3
3. T5-RBS 1018 4
4. HlyA (digestion/hifi)
5. S100-Rose-1018 8
6. S100-RBS-GFP 2
7. RBS-CheZ-GFP 1

Testified by

(.....)

Supervisor

(.....)

1 September 2014

Loading Dye 6X dan LD 10X

$M1.V1 = M2.V2$

$10XLD.X = 6XLD.60$

$10X = 360$

$X = 36$

$10XLD.X = 6XLD.100ul$

$10X = 600ul$

$X = 60ul$

$100-60ul = 40TE$

Making Ladder DNA (marker)

- | | |
|-----------------------|-----------------|
| 1. Goe Ruler Mix 50ug | = 10ul |
| 2. Loading Dye 6x | = 20ul |
| 3. TE 1X Buffer | <u>= 70ul</u> + |
| | = 100 ul |

Agenda Lab 1 September 2014

1. Nanodrop Subtilis → Running
2. Check Phenol in fumehood
3. (PCR Mix → Hotstar → Template B.Subtilis) 5x reaction → touchdown overnight
4. Inoculate P1018 → Column 3 → 20ml
5. Inoculate T5 RBS 1018 3, 4 → 20ml
6. Digestion → running electrophoresis
7. Check biofilm

Inoculate B. Subtilis, P. Aeruginosa, S. Aureus, E. WT, V. cholerae 20 ml

Check LB Broth

Make Plate Agar Stock

Running Electrophoresis 1 Sep 2014

M-T5 RBS-T5 RBS 1018 4-T5 RBS 1018 3-S100 Rose 1018 8-S100 Rose GFP 2-RBS CheZ GFP 1

Using 0,2-2 ul micropipet was not optimal, so we worry that the digestion was not perfect

HotStar PCR

Touchdown

Gradient Temp. Column

1. 1 50°C
2. 3 54,7 °C
3. 5 55,5 °C
4. 7 61,8 °C

Testified by

Supervisor

(.....)

(.....)

5. 10 68,4 °C

MIX

$\sum \text{H}_2\text{O} = 31,63 \text{ ul}$

10x taq buffer = 5,5 ul

5x Qsol 11ul

10mM dNTP mix 1,1 ul

Taq HotStar 0,28 ul

PRIMER

Substilis F 1,375 ul

Substilis R 1,375 ul

TEMPLATE

Genome 1/3 EB 5ul

MIX ALIQUOT 5 tube PCR → PCR Touch Down HotStar

Nanodrop Genome B. Substilis in 1/3 EB

	Concentration	A260	A280	260/280	260/230
Blank	0	1	-0,002	0,17	-0,05
1/3 EB B.subs 1	424	8,481	4,716	1,8	1,74
1/3 EB B.subs 2	306,5	6,129	3,442	1,78	1,76

Testified by

(.....)

Supervisor

(.....)

2 September 2014

Running Electrophoresis

Subtilis 1,2,3,4,5-Marker-HlyA-T5 RBS 1018 3- T5 RBS 1018 4-T5 RBS-S100 Rose 1018 8-S100 Rose
GFP 2-RBS CheZ GFP 1

Ligation

S100-RBS 3 II Nuc-HlyA 2 \rightarrow 16 °C overnight

S100-RBS 3 II Nuc-HlyA 2 \rightarrow 23 °C 1 hour

S100-RBS 3 II MalS-HlyA 2 \rightarrow 16 °C overnight

S100-RBS 3 II MalS-HlyA 2 \rightarrow 23 °C 1 hour

Transform

1. Kont (-) \rightarrow SOC <<
2. Kont (+) \rightarrow isi RFP Amphi 1
3. S100-RBS 3 II Nuc-HlyA 2 (23 °C)

To Do List

Electrophoresis \rightarrow staining

Ligase (AC, heatblock)

Transformation

Isolation

Testified by

(.....)

Supervisor

(.....)

3 September 2014

Take out the transformation

Making the replica

PCR colony, inoculate in to the falcon, 20ml LB Broth

Ligase 16 °C overnight → heatkill 80 °C for 20 minute transform

Take out the transformation, replicate S100 RBS Nuc HlyA → template

PCR mix, transform 16 °C start → put the mix in ice for 60 minute

Digest → ligase (30 °C 1 hour, 88 °C 20 minutes)

Take out the transformation from ice and start PCR

Ligase (23 °C 1 hour, 88 °C 20 minutes)

Transform 60 ° → 1 hour in ice

Digestion

	Eco + Spe= MalS	Xba+Pst= HlyA	Eco+Pst=RFP 2 chloramphenicol
NEB 4	2,5	2,5	2,5
BSA 10X	0,5	0,5	0,5
Enzim 1	0,5	0,5	0,5
Enzim 2	0,5	0,5	0,5
DNA	2	3,6	3
Σ H ₂ O	14	12,4	13

Transform : S100-RBS-Nuc-HlyA ligase overnight

Nanodrop

Sample	CONC	260/280	260/230
1018 3	30,1	2,07	1,44
T5 RBS 1018 3	33,1	2,17	1,69
T5 RBS 1018 4	37,4	2,06	1,39

Digestion

	½ DNA	½ H ₂ O	
1018 -3 (2 sept)	4,2	3,8	NEB 1,25 BSA 0,25 EcoR1 0,25
T5-RBS-1018 3 (2 sept)	3,8	4,2	NEB 3 5ul BSA 10X 1ul EcoR1 1ul Pst1 1ul @2 ul
T5-RBS-1018 4 (2 sept)	3,4	4,6	
T5-RBS-1018 4 (28 august)	2	6	

Running

M-HlyA-MalS-RFP 2-Sub1-Sub2-T5 RBS 1018 4(28Agt)-T5 RBS 1018 3(2 Sept)- T5 RBS 1018 4(2 Sept)- T5 RBS 1018 4(28 August)

Testified by

Supervisor

(.....)

(.....)

6 September 2014

- a. TolC HlyBD
- b. MalS HlyA
- c. T5 rbs 1018

Nanodrop Isolation Result

Sample	Concentration	A260/280	A260/230
Kontrol (-)	10,8	0,84	0,55
S100 RBS Nuc HlyA 1	67,3	2,04	2,05
S100 RBS Nuc HlyA 3	52,9	2,08	1,8
S100 RBS Nuc HlyA 4	57,6	2,03	1,64
S100 RBS Nuc HlyA 5	58,5	2,09	1,77

Digestion (double cek)

	½ DNA	½ H ₂ O
S100 RBS Nuc HlyA 1	2	6
S100 RBS Nuc HlyA 3	2,4	5,6
S100 RBS Nuc HlyA 4	2,2	5,8
S100 RBS Nuc HlyA 5	2,2	5,8

NEB 3 1,25
BSA 10X 0,25
EcoR1 0,25
Pst1 0,25

Running

M-T5 RBS 3-T5 RBS 1018 3-T5 RBS 1018 4-RFP Amphi-S100 Rose GFP-M

M-S100 RBS Nuc HlyA 1-S100 RBS Nuc HlyA 3- S100 RBS Nuc HlyA 4- S100 RBS Nuc HlyA 5-Nuc HlyA 2-M

Testified by

Supervisor

(.....)

(.....)

7 September 2014

Nanodrop 1/3 EB Bacillus (dilution 1/50)

	Concentration	260/280	260/230
Blank	0,3	0,94	0,38
1/3EB Bac 1	8,2	2,13	2,08
1/3 EB Bac 2			

Nanodrop 1/3 EB Bacillus 1/10 dilution

	Concentration	260/280	260/230
1/3 EB Bac 1 (1/10)	44,5	1,72	1,7
1/3 EB Bac 2 (1/10)	27,9	1,79	1,64
Blank	0,4	1,5	0,33

Testified by

(.....)

Supervisor

(.....)

8 September 2014

Digestion (double) → for check list

	½ DNA	½ H ₂ O
S100-RBS-Nuc-HlyA 1	2	6
S100-RBS-Nuc-HlyA 3	2,4	5,6
S100-RBS-Nuc-HlyA 4	2,2	5,8
S100-RBS-Nuc-HlyA 5	2,2	5,8

Mix

NEB 5

BSA 10X 1 @2

Eco 1

Pst 1

Testified by

(.....)

Supervisor

(.....)

9 September 2014

DNA isolation result 8/9/14

	Concentration	260/280	260/230
S100-RBS-Nuc-HlyA 1	39,4	2,06	1,69
S100-RBS-Nuc-HlyA 3	57,3	2,01	1,59
S100-RBS-Nuc-HlyA 4	70,9	1,99	1,91
S100-RBS-Nuc-HlyA 5	41,6	1,97	1,72
Blank	0,3	0,57	0,23

Digestion (double) DNA isolation 8/9/14

	½ DNA	½ H ₂ O
S100-RBS-Nuc-HlyA 1	3,2	4,8
S100-RBS-Nuc-HlyA 3	2,2	5,8
S100-RBS-Nuc-HlyA 4	1,8	6,2
S100-RBS-Nuc-HlyA 5	3	5

Testified by

Supervisor

(.....)

(.....)

10 September 2014

$V1.M1 = V2.M2$
 $50.M1 = 500 \cdot 0,5$
 $M1 = 5 \text{ mL}$

ElectrophoresisRunning

Marker – SNRH / PCR Colony 1 2 3 4 5 6 7 8 – Peptide 1018 Hifi – Digesti RFP Amphi 1018 Eco Pst –
Marker

Extraordinary occurrence:

Running progress was not finished yet, running was only in the chamber.

Agarose was taken out from the chamber, so we worry that will soften the agarose if it was soaked overnight in the TAE chamber

Testified by

(.....)

Supervisor

(.....)

11 September 2014

Running electrophoresis

SNRH 1 2 3 4 5 6 7 8 – Marker – Subtilisin Hifi (51 is genom, 51 HS, 54 HS) – Peptide 1018 Hifi 1 + 5

- ✓ Make LB Broth 400 mL
- ✓ Make LB agar 300 mL → 15 petri dish [4 antibiotic (-), 10 (7 chloramphenicol, 3 ampicillin)]
- ✓ Making stock: ampicillin 6, chloramphenicol 3, without antibiotic 2
- ✓ Electrophoresis → making new agar
- ✓ Inoculate SNRH
- ✓ PCR PFX Hifi 10 → Subtilis 51 °C, 5 → 1018 53 °C
- ✓ Inoculate 6 pathogen bacteria → 15 mL

Prepare

6 falcon (1 falcon control)

Nanodrop Subtilis Hifi

Concentration 10,3 ng/ul

260/280	2,28	→ heat	11,1
260/230	0,92		2,01
			11,04

Testified by

Supervisor

(.....)

(.....)

12 September 2014

Ligase

		H ₂ O	Mix Ligase
1018 2	Amphi 2	4,5	5 ul buffer 3
TolC 2-HlyBD 2	Amphi 2	2,5	2,5 ul T4 Ligase 1,5
S100 RBS 2 Nuc HlyA 2	Amphi 2	2,5	1,5 ul @ 1,5
MalS 2- HlyA 2	Amphi 2	2,5	@

Testified by

(.....)

Supervisor

(.....)

13 September 2014 → Yesterday transform result is too little.

- Repeat the transform (left 50 ul, right 200ul)
1018-AMPHI
S100-RBS-N H
MalS-HlyA
- Grow the biofilm
- Autoclave SOC

Running Subtilisin

1. 45,3 °C
2. 45,9 °C
3. 46,7 °C
4. 47,8 °C
5. 49,3 °C
6. 551 °C
7. 52,4 °C
8. 53,5 °C
9. 54,3 °C

Testified by

(.....)

Supervisor

(.....)

15 September 2014

iGEM lab schedule

1. PCR HotStar Subtilisin
2. Take out the replica
3. Running PCR Colony Running Hifi Result
4. Inoculate PCR Good Colony, there's starter culture
5. If the running isn't good → ligase again with TolC-HlyBD and the other unwell parts
6. Ligation of TolC HlyBD

Ligation :

TolC 1,6ul

HlyBD 4ul

RFP 2ul

Buffer 1ul

Ligase 0,5ul

Air 0,9ul

PFX –

PCR =

PCR HS → ?

1. Running PCR Colony
Marker – (Amphi-P1018 1 2 3 4 5) – (MalS-HlyA 1 2 3 4) – (S100-RBS-Nuc-HlyA 1 2 3)
2. Pcr HotStar Subtilisin
3. Ligation
 - TolC : HlyBD : amphiRFP = 1,6 : 4 : 2 ∑ H₂O 0,9
 - S100-RBS: Nuc HlyA : amphi RFP = 0,2 : 0,8 : 2 ∑ H₂O 5,5
4. Inoculate
 - 1018 Amphi 1 & 5
 - MalS-HlyA 2 & 3
 - S100 RBS Nuc HlyA 2 & 3 (from replica 3/9/14)
 - S100 RBS Nuc HlyA 1 & 3 (from replica 14/9/14)
5. Change the Broth (bacteria's biofilm)
6. Running HiFI Subtilis Purif & 1018 Hifi

Running electrophoresis 15 Sept 2014

1018 Amphi 336

S100 RBS Nuc HlyA

35 12 561 186 + 300 + 1094

MalS HlyA

2067 186 + 300 = 2553

Testified by

Supervisor

(.....)

(.....)

16 September 2014

Lab iGEM

1. DNA isolation of Running HS
2. Digest, Running → SNRH, MalS HlyA, P1018 Amp, SNRH Rose
3. LMA → PCR PFX P1018
PCR HS Subtilisin
4. Staining biofilm
5. Transform → spread (ampicilin)

Plate 2 Well 3F Chloram S100+RBS Kit

To do list target :

1. LMA P1018 PCRPFX Hifi : delayed → 2%
LMA Subtilisin PCR Hotstar : delayed → 1,5%
2. Isolasi :
1018-Amphi 1 & 5
MalS-HlyA 2 & 3
S100-RBS-Nuc-HlyA 2 & 3 (from replica 3/9/14) → there's band
S100-RBS-Nuc-HlyA 1 & 3 (from replica 14/9/14) → there's no band
3. Nanodrop
P1018 PCR PFX Hifi unpurified
Subtilisin PCR PFX Hifi purified
Subtilisin Hotstar unpurified

Running
Marker – PCR HotStar Subtilisin

Nanodrop

Sample	Concentration	A260	A280	260/280	260/230
PCR HS Subtilisin	3,6	6,673	5,621	1,19	1,43
	131,5	2,627	1,413	1,86	0,34
	135,3	2,709	1,453	1,86	0,24
PCR PFX Hifi 1018	179,5	3,469	1,889	1,84	2,05
	169,8	3,395	1,848	1,84	2,2

Testified by

Supervisor

(.....)

(.....)

17 September 2014

Nanodrop

		Σ H ₂ O	Conc	DNA	260/280	260/230
	K (-)		0,7			
1	SNRH 1 (14)	4,9	94,6	1,5	1,98	1,96
2	SNRH 2 (3)	4,61	83,6	1,79	1,94	1,99
3	SNRH 3 (14)	4,8	93,5	1,6	2,12	1,65
4	SNRH 3 (3)	2,26	36,2	4,14	2,3	1,64
5	1018 Amp 1	4,75	90,8	1,65	2,09	1,9
6	1018 Amp 5	4	73,7	2	2,08	1,99
7	MalS HlyA 2	3,38	49,6	3,02	1,96	1,86
8	MalS HlyA 3	1,05	28	5,35	2,06	1,46

18 September 2014

List of Digestion

29 August 2014	RBS CheZ II GFP 1	1,39	57,7	2,5	2,09	1,19
-------------------	----------------------	------	------	-----	------	------

Testified by

(.....)

Supervisor

(.....)

18 September 2014

Digest the isolation of RBS CheZ GFP ✓ → Running ✓

Take out the transformant ✓ → Inoculate ✓

PCR Hotstar subtilis / evening ✓

Inoculate the transformant S100 RBS Device / evening ✓

Make the bacteria stock -80°C + streak → pathogen's biofilm / evening ✓

Take out the inoculation from incubator shaker ✓

Evening

Running PCR Colony

Marker – SNRH 1 2 3 4 5 6 – TolC HlyBD 1 2 – Hotstar 1 2 – Marker

M – SNRH 1(14) - SNRH 3(14) - SNRH 2(3) - SNRH 3(3) – 1018 Amphi 1 5 – MalS HlyA 2 3

Repeat the inoculate the TolC-HlyBD

Nanodrop

	Concentration	A260/280	260/230
SNRH 3 (from replica 16/9/14)	65,4	2,02	1,96
Blank	0,2		

Inoculate S100-RBS

Testified by

Supervisor

(.....)

(.....)

19 September 2014

1. Running PCR Substilis
2. Isolation
 - TolC HlyBD 1A
 - TolC HlyBD 1B
 - S100 RBS 1
 - S100 RBS 2
 - S100 RBS 3
3. Take out the pathogen bacteria
4. Digestion (repeat) + (new 5)

		6,4 ul		DNA 150 ng MIX ½ Reaction (V=8ul) NEB Buffer 1,5 ul (@1ul) BSA 10X 3 ul (@0,2 ul) EcoR1 3 ul (@0,241) Spe1 3ul (@0,241) @1,6 ul
		DNA	Σ H ₂ O	
1	SNRH 1 (14/9)	2	4,4	
2	SNRH 2 (3/9)	2	4,4	
3	SNRH 3 (14/9)	2	4,4	
4	SNRH 3 (3/9)	2,3	4,1	
5	1018 Amp 1	2	4,4	
6	1018 Amp 5	2	4,4	
7	MalS HlyA 2	3	3,4	
8	MalS HlyA 3	3,4	1	
9	RBS CheZ GFP	2,5	3,9	
10	TolC HlyBD 1A	1,8	4,6	
11	TolC HlyBD 1B	0,7	5,7	
12	S100 RBS 1	2	4,4	
13	S100 RBS 2	1,8	4,6	
14	S100 RBS 3	3	3,4	
15	SNRH 3 (16/9)	2,3	4,1	

5. PCR Substilis
6. Nanodrop

	Concentration (ng/ul)	A260/280	260/230
Blank	-0,1		
TolC HlyBD 1A	102,4	1,95	2,05
TolC HlyBD 1B	288	1,9	2,11
S100 RBS 1	92,5	1,95	2,02
S100 RBS 2	105,2	1,94	2,13
S100 RBS 3	59,5	1,98	2

Testified by

Supervisor

(.....)

(.....)

20 September 2014

Digestion (assembly) 150 ng

				Σ H ₂ O	DNA
	[S100	K	7,96	2,04
		MalS HlyA 2	Chlor	7	3
A B C		SNRH 1 -14	K	8,42	1,58
	[T5 RBS	K	8,4	1,6
		1018 Amphi 5	K	8	2
		RFP chloram		8,2	4,8

TH 1 A, TH 1B, PA, PAS, MN 2, MN 3, SR 2, SR3
SN2-3, SN3,3 , SN3-6, SN3 14, SN1-14, SR, RC

Testified by

Supervisor

(.....)

(.....)

22 September 2014

1. Single Digest

TolC HlyBD 1A 125 ng

$\sum \text{H}_2\text{O}$: 6,8

$\frac{1}{2}$ DNA: 12

NEB Buff: 1,25

BSA 10X: 0,25

Xba I : 0,25

2. Ligation

1. S100 - MalS HlyA – RFP chloramphenicol

0,03 2,18 2

0,5 2,18 2 $\sum \text{H}_2\text{O}$ 3,8

2. T5 RBS – 1018 – RFP chloramphenicol

○ T5-RBS

$(138/2070)*25=1,6$

$(1,6/25)*2\text{ul}=0,13 \text{ ul}$

○ 1018

$(36/2070)*25 = 0,03 \text{ ul}$

T5 RBS – 1018 – RFP chloramphenicol

0,13 0,03 2

0,5 0,5 2 $\sum \text{H}_2\text{O}$ 5,5, Buffer T4 3uL, T4 Ligase 1,5 uL @ 1,5 uL

3. SNRH – RFP chloramphenicol

$(796/2070)*25=9,5$

Buff T4 3ul

T4 Ligase 1,5 ul

@1,5 ul

LIGATION

$(\text{Insert length} / \text{vector length}) * \text{vector mass}(\text{ng}) = \text{insert mass}(\text{ng})$

Molar ration 1:1

Vector 25 ng/2 ul

Example:

TolC –HlyBD –RFP

1500 4000 2000

1,5 4 2

TolC $\rightarrow (1500/2070)*25 = \dots\dots \text{ ng}$

$(\dots\dots\text{ng}/25)*2\text{ul} = \dots\dots \text{ ul}$

Testified by

Supervisor

(.....)

(.....)

HlyBD →

S100 MalS HlyA

S100 : $(35/2070) \times 25 = 10,42$

$(A/25) \times 2 \text{ul} = 0,03 \text{ ul}$

MalS HlyA : $(2253/2070) \times 25 = 27,21 \text{ ng}$

$(27,21/25) \times 2 \text{ul} = 2,18 \text{ul}$

S100 MalS HlyA RFP chloramphenicol

0,03 2,18 2

1,5

PCR Subtilis HotStar (51, 53)

Template Bacillus 1/10

Running Hasil Double Digest kemarin(20/9/14) and Single digest TolC HlyBD (22/9/14)

TolC HlyBD non digest – TH 1A single digest – TH 1B double digest – TolC Hifi – TolC double digest
12/9 – HlyBD double digest 12/9 (Nhe Pst)

S100-T5 RBS – MalS HlyA – 1018 Amp – SNRH –RFP Chloramphenicol

TolC = $(1539/2070) \times 2 = 1,5 \text{ul}$

HlyBD = $(4000/2070) \times 2 = 3,4 \text{ ul}$

RFP Amp = 2ul

$\sum \text{H}_2\text{O} = 1,1$

Testified by

(.....)

Supervisor

(.....)

24 September 2014

1. Take out the transforman → in fridge for 1 hour 11.00 → incubator 12.00
2. Make LB Broth + LB agar
New Tips (5 boxes)
3. Inoculate replica on 12.00 o'clock
4. Revitalize bakteri 10.15 a.m
5. Inoculate bacteria transforman 12.00 a.m

Testified by

(.....)

Supervisor

(.....)

25 September 2014

1. Running PCR Colony
MH1-MH2-SN1-SN2-SN3-SN4-SN5-T181-T182-T183-T184-T185
2. Isolation
T518 1, 2, 4
SNRH 1, 2, 3

Testified by

(.....)

Supervisor

(.....)

26 September 2014

1. Nanodrop

		Concentration	260/280	260/230
Isolation 25/9/14	SNRH 1	66,9	1,91	1,96
	SNRH 2	44,6	1,90	2,07
	SNRH 3	79,9	1,91	2,07
	T518 1	58,7	1,89	2,04
	T518 2	53,7	1,95	2,06
	T518 4	62,4	1,87	1,9
	Blank	0,7	0,67	0,52

2. Ligation

- TolC HlyBD
- MalS HlyA

3. Digestion

	6,4 ul		MIX NEB Buffer 3 6ul BSA10X 1,2ul Eco 0,6ul Pst 0,6ul
	Σ H ₂ O	DNA	
SNRH 1	2,4	3	
SNRH 2	2	4,4	
SNRH 3	3,9	2,5	

Nanodrop Result DNA Isolation September 26th 2014

	Concentration	260/280	250/230
Blank	2		
1018 3	34,1	1,96	2,23
RFP Amphi	26,5	2,16	1,78
SNRH 3	43,3	2	1,74
MalS HlyA	28,7	1,96	1,69
Tr 1018 1	18	1,94	1,43
Nuc HlyA	21,8	2,09	1,36

Running Electrophoresis List

MalS HlyA Insert
Backbone RFP Amphi
MalS HlyA
R 1018 1
RFP Amphi
1018 3
Nuc HlyA
SRNH

Testified by

Supervisor

(.....)

(.....)

29 September 2014

	Digestion	DNA 150	Σ H ₂ O
1	1018-3	4,4	2,2
2	TR 18	6,6	0
3	NH	6,6	0
4	SRNH-3	3,5	3,1
5	MH Amphi	5,2	1,4
6	RFP Amphi	5,5	1,1
		DNA 200	Σ H ₂ O
7	MalS HlyA 2	4	2,6
8	RFP Chloramphenicol	2,4	4,2

Mix

NEB 3	@ 1 uL	x 6 =	6 uL
BSA 10 X	@ 0,2 uL	x 6 =	1,2 uL
EcoR I	@ 0,1 uL	x 6 =	0,6 uL
Pst I	@ 0,1 uL	x 6 =	0,6 uL
Total	@ 1,4 uL		7,4 uL

Testified by

(.....)

Supervisor

(.....)

PCR Touchdown Subtilis

1. 50
2. 50,3
3. 50,9
4. 51,7
5. 52,8
6. 54,3
7. 36
8. 57,9
9. 58,5
10. 39,3

Digestion Receipt 30/9/2014

Mix

NEB 3	@ 1 uL	x 11 =	11 uL
BSA 10 X	@ 0,2 uL	x 11 =	2,2 uL
EcoR I	@ 0,1 uL	x 11 =	1,1 uL
Pst I	@ 0,1 uL	x 11 =	1,1 uL
Total	@ 1,4 uL		15,4 uL

Rumus DNA 150- $\sum H_2O$

Volume template = $150/\text{concentration} = x$

Volume $\sum H_2O = 6,0 - x = y$

$x + y = 6,6 \text{ uL}$

Total per tube = mix 1,4 uL + 6,6 uL = 8 uL

Testified by

(.....)

Supervisor

(.....)

Marker – 1018 3 x x x x – SNRH – Marker – MalS HlyA – RFP Chloram – RFP Amphi – MH Amphi

Nanodrop DNA Isolation 30/9/2014

No	Sample	Concentration	260/280	260/230	x	y ($\sum H_2O$)
1	Blank	0,1	1,02	0,09		
2	T5 RBS 1018	44,8	1,92	2,11	3,35	3,25
3	1018 3	76,5	1,94	2,28	1,96	4,64
4	SNRH 3	104,7	1,89	2,2	1,43	5,17
5	S100 MalS HlyA 4	161,0	1,88	2,25	0,93	5,67
6	S100 MalS HlyA 5	142,1	1,88	2,25	1,05	5,55
7	TolC HlyBD 2	173,1	1,88	2,16	0,86	5,74
8	TolC HlyBD 4	96,7	1,89	2,08	1,55	5,05
9	Nuc HlyA 2	55,4	1,94	1,88	2,7	3,9
10	RFP Amphi	27,8	2,01	1,2	5,39	1,21
11	TolC HlyBD 3	335,1	1,89	2,21	0,45	6,15
12	MalS HlyA 2	72,7	1,9	1,37	2,06	4,54

Testified by

(.....)

Supervisor

(.....)

1 October 2014

iGEM Lab Schedule

1. PCR Substilis 10 reaction scale up

Mix	Primer mix = 5,25
DW 60,38	Sub F = 2,625
SQsol 21,0	Sub R = 2,625
10 x taq 10,5	Taq 0,53
MgCl 2,1	Template 5,25
dNTP 2,1	

Annealing Temperature 51°C

2. PCR Colony Replica 11 Reaction

PCR Tube Code	PCR Coloby Replica
A	Sp100-Rose 1
B	Sp100-Rose 2
C	Sp100-Rose 3
D	Sp100-Rose-1018 1
E	Sp100-Rose-1018 3
F	Sp100-Rose-1018 5
G	Sp100-Rose-1018 6
H	Sp100-Rose-1018 8
I	Sp100-Rose-1018 10
J	RNA rose
K	RNA rose

PCR Colony transforman MalS HlyA 6 colonies
Make the replica

Testified by

Supervisor

(.....)

(.....)

Repeat the *MalS HlyA* 2 digestion

EcoRI = 0,1 uL

PstI = 0,1 uL

NEB 3 = 1 uL

BSA 10x = 0,2 uL

DW = 6,6 – 2,06 = 4,54

Template DNA = 150/72,7 = 2,06

Bikin LB Agar 20 Petri Dishes

6 non antibiotics

4 ampicillin

10 chloramphenicol

Electrophoresis Running

Testified by

(.....)

Supervisor

(.....)

3 October 2014

	Digestion	250 ng DNA	Σ H ₂ O	NEB @ 2 Enzyme @ 0,2 Enzyme @ 0,2
1	T5-RBS 3	2,7	5,3	
2	TolC HlyBD	0,75	7,25	
3	RFP Chloramphenicol	1,5	6,5	

0,238 gram/mL → 1M

Testified by

(.....)

Supervisor

(.....)

5 October 2014

Ligation

T5 RBS	1 uL	T4 Buffer	1 uL
TolC HlyBD	4,8 uL	T4 Ligase	0,5 uL
RFP chloram	2 uL		
Total	7,8 uL	Total	1,5 uL
TOTAL 7,8 + 1,5 = 9,3 uL + 0,7 uL H ₂ O = 10 uL			

1. Ligation T5-RBS II TolC-HlyBD → transform ✓
2. Isolation, nanodrop, running ✓

		Concentration	A260/280	260/230
S Rose 2	K	54,1	1,99	2,12
S Rose 1018 1	K	41,6	1,91	1,85
MalS HlyA 2	K	121,5	1,9	2,18
MalS HlyA 3	K	112,6	1,91	2,09
MalS HlyA 5	K	142,9	1,89	2,16

3. PCR Subtilisin, Running ✓
Running: Marker – PCR Subtilisin 1 2 3 4 5 6 7 8 9 10 – Marker
4. Make LB Broth and SOB Media 500 mL (without MgCl₂) ✓
5. Stock reuse microtube 1,7 mL & new one for induction ✓
Stock cold tips for making competent cell
6. Digestion of Isolation ✓

		200 DNA	Σ H ₂ O	NEB Buffer @ 2 x 5 = 10 EcoRI @ 0,2 x 5 = 1 Pst I @ 0,2 x 5 = 1
1	SR 2	3,7	3,9	→ Running electrophoresis ✓ Marker-SR181-SR2-MH2-MH3-MalS-Marker
2	SH 8	4,8	2,8	
3	MH 2	1,7	5,9	
4	MH 3	1,8	5,8	
5	MH 5	1,5	6,1	

7. Inoculate
 - 14/09 SNRH 3 (K) 15 mL
 - T5 RBS 1018 1 (K) 15 mL
 - S100 MalS HlyA 4 (K) 15 mL
 - Nuc HlyA 2 (K) 15 mL
 - MalS HlyA 2 (A) 15 mL
 - Pirare (K) 10 mL
 - BL21(-AB) 10 mL
 - Top10 (-AB) 10 mL
 - DH5 alpha (-AB) 10 mL

Testified by

Supervisor

(.....)

(.....)

6 October 2014

- Prepare the sample for submit
- Make competent cell
BL21
DH5 alpha
Pirare
Top10

No Tube	Nanodrop	Concentration	A260	A280	260/280	260/230
21	Reverse Cqss	116,9	3,339	1,613	2,07	2,5
22	Reverse Cqss	114,3	2,286	1,080	2,12	2,08
23	Reversed ybcl	125,4	2,509	1,180	2,13	2,08
24	Reversed ybcl	132,9	2,657	1,247	2,13	2,28
25	Reversed yain	79,3	1,586	0,725	2,19	2,62
26	Reversed yain	81,0	1,620	0,733	2,21	2,69

Innoculate biofilm pathogen bacteria 2 well plate

1. Well Plate Non IPTG Gradient Column

	K-	E. coli	V. cholerae	B. subtilis	S. aureus	P. aeruginosa	K. pneumonas
1	O	o	O	o	o	o	o
2	O	o	O	o	o	o	o
3	O	o	O	o	o	o	o
4	O	o	O	o	o	o	o

2. Well plate IPTG Gradient Column for Peptide 1018

K-
E. coli
V. cholerae
B. subtilis
S. aureus
P. aeruginosa
K. pneumonas

3. Revitalize bacteria for innoculating to M9 tomorrow 7 October 2014
4. Make M9 Medium

Testified by

Supervisor

(.....)

(.....)

IPTG

$$V1.M1=V2.M2$$

$$V1.500=0,05.200$$

$$M2=0$$

		M2	V2	V1	M1=0,5M=500mM → 50 mM → 10 mM
1	0,2 uL	0,05	200 uL	0,02	
2	0,4 uL	0,1	200 uL	0,04	
3	1 uL	0,25	200 uL	0,1	
4	2 uL	0,5	200 uL	0,2	
5	3 uL	0,75	200 uL	0,3	
6	4 uL	1	200 uL	0,4	
7	5 uL	1,25	200 uL	0,5	
8	6 uL	1,5	200 uL	0,6	
9	7 uL	1,75	200 uL	0,7	
10	8 uL	2	200 uL	0,8	

$$V1.M1=V2.M2$$

$$V1.500mM=1000uL.50mM$$

$$V1.5=500$$

$$V1=100 \text{ uL}$$

$$1000-100= 900 \text{ uL H}_2\text{O}$$

$$V1.M1=V2.M2$$

$$V1.50=200uL.10$$

$$V1=40 \rightarrow 160 \text{ uL H}_2\text{O}$$

For Tomorrow → 8 October

IPTG Concentration	M9 Shake Added Volume	
0,05 Mm	199 uL	IPTG Stock has used 10 mM
0,1 mM	198 uL	
0,25 mM	197 uL	
0,5 mM	198 uL	IPTG Stock 50 mM
0,75 mM	197 uL	
1 mM	196 uL	
1,25 mM	195 uL	
1,5 mM	194 uL	
1,75 mM	193 uL	
2 mM	192 uL	
0 mM	200 uL	

Extraordinary occurrences

Characterizations activity has encountered obstacle not only biofilm but also enzyme expression.

Details :

- P1018 Characterization
 - ➔ M9 can't be used
 - ➔ Undergo Millard effect
 - ➔ Failure : M9 + Glucose can't be autoclaved "caramelize" → brownish M9
- Wrong receipt for M9

Evaluation :

Testified by

Supervisor

(.....)

(.....)

1. Composition and presentage must be written → write the clear information
 2. Glucose has not mixed in the beginning + has not autoclaved → filtered
 3. There's no sterile bottle → prepare the bottle!
- Can't do spectrofotometri test → sample was not ready yet

Sequencing

T5 RBS Peptide 1018] @ 10ng/uL 8uL + primer forward, primer reverse @ 5 mM 1 uL 2 reaction
Sp100 RBS Nuc HlyA	
Sp100 MalS HlyA	

Testified by

(.....)

Supervisor

(.....)

1. Clean up the biofilm
 2. New inoculation LB Broth and M9
 3. Take out the replica
 4. Staining
-
1. 96 well plate was filled with LB Broth 198 uL
 2. Add pathogen bacteria 2 uL and incubate overnight in static incubator
 3. Well was filled crystal violet 0,5% 200 uL → incubate 10 minute
 4. Crystal violet was took out → put in to the bottle = reused crystal violet
 5. Wash the well, soak it in the water and wah for 4 times
 6. Well is contain crystal violet that pinned in biofilm
 7. Well was filled 96% ethanol (v/v) and acetate acid 2%
96% absolute ethanol + 2% CH₃COOH (v/v)
96 uL absolute ethanol + 2 uL acetate acid + 2 uL H₂O
 8. Wait for 2 minutes → bring it into pass through box
 9. Elisa reader

Testified by

(.....)

Supervisor

(.....)

7 October 2014

Make Stock Agar

1. Prepare LB agar from autoclave. LB starch → LB agar and aquades has heated in the oven
300 mL LB sterile agar + 200 mL sterile aquades
Starch powder 1,5 gr in aluminium foil put in the pass through box
2. 200 mL aquades → pour in 500 mL bottle glass + 1,5 gr starch powder, heat in the oven until dissolve
3. Discard 120 mL. The excess 80 mL starch. Add 120 mL LB agar, then aliquot it to 10 petri dishes (@ 20 mL).
4. Add chloramphenicol 200 ng/uL

Innoculation the bacteria to M9

1. Preapare M9 medium and overnight bacteria culture in LB Broth
2. Take 400 uL bacteria from LB broth, drop into 20 mL M9
1 control (-) + 5 bacteria :
T5 RBS P1018 → M9 without glucose
SRNH
SMH
RNH
MH
3. Put in to the shaker incubator 37°C overnight

Testified by

Supervisor

(.....)

(.....)

8 October 2014

1. Sequencing
Sample template 8 uL
 - SRNH: S100-RBS-Nuc-HlyA
 - SMH
 - T5 RBS PeptidePrimer 5 uM @ H₂O → VF2 & VR @ 1uL per reaction
2. LB revitalization
3. 1018 activation

Testified by

(.....)

Supervisor

(.....)

10 October 2014

V1.M1=V2.M2

96.1000=70.X

X=1371

Plate 1 : 4 x 7

LB(-) E. coli – V. Cholerae – B. Subs – S. Aureus – P. Aeuroginosa- K. Pneumonas

Plate 2

LB(-) E. coli – V. Cholerae – B. Subs – S. Aureus – P. Aeuroginosa- K. Pneumonas

Transformation (NEB/T5)

1. Competent cell + gene of interest (10% comp → ligase 0,5-2 uL)
2. Incubate in ice for 30 minutes
3. Heatshock 55 s
4. Incubate ice 4⁰C for 2 minutes.
5. Add (2x competent cell) SOC
6. Shaker incubation 1 hour
7. Plate M9

Lower cell

Polyacrilamide 30% 1,4 mL = 1400 uL 2800

Tris HCl 1,5 M pH 8,8 0,875 mL = 875 uL 1750

SDS 10% 35 uL

APS 17,9 uL

TEMED 3,5 uL

H₂O

3500 – 2331,4 uL

H₂O = 1168,6 uL

Tris HCl 6,8 375 uL

SDS 10% 15 uL

Polyacrilamide 200 uL

d H₂O 901 uL

APS 10% 7,5 uL

TEMED 1,5 uL

1,5 total

1,5 1500 uL → 901 uL

Testified by

(.....)

Supervisor

(.....)

11 October 2014

Characterization

A. Biofilm Degradation

1. Biofilm growth from 2 days ago (9/10/14) is treated
 - MH SNRH NH
 - SNRH MHX (+)
 - Then, staining
2. Running SDS Page → overnight
 - From Device at M9

B. Antimicrobial Assay

- Dilution 1:20 from inoculation M9 (V=5 mL) was turbid (hope OD still high)
 - a. 1018
 - b. T5 RBS 1018 → in falcon → well plate Tripto → Elisa reader OD 595 nm
- Give IPTG induction

C. Starch Degradation Assay → MaIS

- Make starch agar 200 mL → 4 gr starch 100 mL, 3,....gr LB agar 100 mL
- Centrifuge 1,5 mL inoculation result device (E. coli WT, SMH, MH) in LB
- Take supernatan, drop in to LB and starch agar

E. coli WT	E. coli SMH
E. coli (-)	E. coli MH

Time variable : 16 hours, 24 hours, 32 hours, 48 hours (?)

- Take the pellet to LB and pellet from sonication, drop and streak it on the plate
 - Pellet → live cell that still have metabolic activity
 - primer metabolism → grow, biner
 - secondary metabolism → enzyme secretion
 - Pellet → is worried that will add bias variable
 - because the cell is still live and can produce extracellular (amylase)
 - But, if it is sonicated → pellet → lysis = cell will death
 - Sonication = cell will death → there's no primer metabolic activity, secondary.
 - Cell was lysis → can check enzyme in cell → the lysis cell is intracellular without activity variable

Testified by

Supervisor

(.....)

(.....)

12 October 2014

1. Flooding Iodine
16 Hours [LB agar & LB starch = Supernatan] [LB starch=pelet]
24 Hours [LB agar & LB starch = Supernatan]

2. SDS Page

Supernatan				Marker	Pellet
SMH	MH	SRNH	NH	SMH	SRNH
M9	M9	M9	M9	M9 1	M9

SMH LB	MH 9	MH LB pelet	SMH LB pelet	Marker
--------	------	----------------	-----------------	--------

3. IPTG Induction P1018
Use LB
4. Repeat inoculation 96 well plate (2) → for biofilm
6 pathogen bacteria
5. Revitalize bacteria
E. Coli mutant
SMH MH SRNH NH
6 pathogen bacteria

Testified by

Supervisor

(.....)

(.....)

13 October 2014

1. Flooding Iodine 36 hours
2. Inoculate E. coli mutant
 1. MH2
 2. SRNH 4
 3. NH2
 4. SRNH3
 5. P1018 3
 6. T5 RBS P1018 1
 Pathogen bacteria → bacteria revitalization
 1. E. coli
 2. V. Cholerae
 3. Bacillus subtilis
 4. Pseudomonas aeruginosa
 5. Staphylococcus aureus
 6. Kleoseilla pneumoniae

3. Repeat Biofilm Protocol

- Take out the LB + pathogen bacteria
- Fill new LB 200 uL + 2 uL pathogen based on protocol

	MH2	SMH4	NH	Plate 1
	1 2 3 4	5 6 7 8	9 10 11 12	
E. coli				
V. cholera				
B. Subtilis				
P. aeruginosa				
S. aureus				
K. pneumoniae				
Blank				
LB (-)				

	SRNH	MIX SRNH 3 SMH 4	No E.coli Mutant	Plate 2
	1 2 3 4	5 6 7 8	9 10 11 12	
E. coli				
V. cholera				
B. Subtilis				
P. aeruginosa				
S. aureus				
K. pneumoniae				
Blank				
LB (-)				

○

4. Isolation, digestion

Blank	0,4	1,52	0,67
SMH	31,3	1,99	1,13
SRNP	33,8	2,09	1,47
NH	116,1	1,91	0,85
T5	4	1,94	1,18

Testified by

Supervisor

(.....)

(.....)

MH	130,8	1,99	0,8
1018	127,4	1,97	1

M9

100 mL M9 Salts	150mL
1mL of 1 M MgSO ₄ steril	0,75 mL
10 mL → 20% glucose	7,5 mL
50 uL 1M CaCl ₂	37,5 mL
+ dH ₂ O → sampai 500 mL	375 mL

Na ₂ HPO ₄ ·7H ₂ O	16 gr
KH ₂ PO ₄	3,75 gr
NaCl	0,625 gr
NH ₄ Cl	1,25 gr
MgSO ₄	1,23 gr
CaCl ₂	0,73 gr

Testified by

(.....)

Supervisor

(.....)

Testified by

(.....)

Supervisor

(.....)